Bedienungsanleitung  
 User manual  
 Mode d'emploi  
 Istruzioni d'uso  
 Modo de empleo

**Amann Gírrbach GmbH**  
 Dürrenweg 40  
 75177 Pforzheim | Germany  
 Fon +49 7231 957-100  
 Fax +49 7231 957-159  
 germany@amanngirrbach.com

Herrschaftswiesen 1  
 6842 Koblach | Austria  
 Fon +43 5523 62333-0  
 Fax +43 5523 55990  
 austria@amanngirrbach.com  
 www.amanngirrbach.com

CE 0124

DIN EN ISO 13485

AMANNGIRRBACH

www.heilblau.com 978140 IN 07/05

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Inhalt: 5 Lotstangen  
 Gewicht: 3,6 g

### PRODUKTBESCHREIBUNG

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Contents: 5 soldering rods  
 Weight: 3.6 g

### PRODUCT DESCRIPTION

Girosolder is a dental solder for cobalt-based alloys. It complies with DIN EN 29333.

Melting range:	1112°C - 1138°C	
Working temperature:	1160°C	
Composition in %	cobalt (Co)	70.3
by mass:	chromium (Cr)	18
	molybdenum (Mo)	5
	silicon (Si)	2.5
	manganese (Mn)	2.5
	boron (B)	1.7

No nickel (Ni), beryllium (Be) or cadmium (Cd) content in accordance with DIN EN 29333

### Recommended fluxes

Fluxes suitable for soldering non-precious cobalt-based alloys (Fluxes 25 g 781560).

### Soldering instructions

#### PREPARING THE SURFACES FOR SOLDERING

Remove any obstructions which may hinder the flow of the solder, e.g. grooves caused by grinding. This involves creating very smooth metal surfaces at the solder joint to ensure it is optimally wetted with solder. Smooth the surfaces to be soldered with crosscut tungsten carbide cutters and, if necessary, rubber polishers and polish to a high lustre. The surface area of the solder joints should be as large as possible and parallel. The solder gap should not be greater than 0.2 mm to ensure a good capillary effect for the solder. When soldering cavities or margin extensions, place platinum foil under the defect being soldered (adapt to the model).

#### FABRICATING A SOLDERING MODEL

Ensure that the solder joint is kept as free as possible so that only minimum heat is lost. The soldering investment model should be as small as possible. Bench dry the soldering model (do not dry in the furnace).

#### CHEMICAL CONDITIONING OF THE SOLDER JOINT

##### Formation of oxides:

oxidation, which begins at low temperatures, can hinder efficient diffusion.

##### Prevention

Flux is required to protect the surface against oxidation and prevent the formation of oxide as far as possible. Apply flux to the solder joint and solder: both the solder and solder joint should be generously coated with flux to ensure successful metal bonding.

##### Adequate supply of heat

A hot, large soldering flame is also essential to ensure that the whole soldering area is uniformly heated to the correct flow temperature of the solder.

##### Heating the solder model

After bench drying the solder model, brush heat uniformly and thoroughly with a sufficiently large flame, then concentrate the heat on the solder joint without removing the flame. Apply the solder already precoated with flux to the solder joint. Do not remove the flame as this may cause oxygen to become entrapped. It is essential that the solder flows immediately. If it forms a ball, the metal is not hot enough. Add more flux and repeat the procedure.