
Bedienungsanleitung
User manual
Mode d'emploi
Istruzioni d'uso
Modo de empleo

BEDIENUNGSANLEITUNG CERAMILL ZI

PRODUKTBESCHREIBUNG

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USER MANUAL CERAMILL ZI

PRODUCT DESCRIPTION

Ceramill ZI Blanks are pre-sintered zirconium oxide blanks (Y-TZP ZrO2).

They serve to produce tight-fitting and removable prosthetic supplies (e.g. crowns and bridges, conus/telescopic crowns, supra-constructions etc.) with the aid of manual copy-milling machines (e.g. Ceramill Base) and CAD/CAM systems.

Ceramill ZI complies with the demands of the standard DIN EN ISO 6872:1999 after the above-mentioned final sintering and also with ISO 13356:1997 with the further processing described below.

HANDLING

Material-specific framework parameters

The following material-specific framework parameters have to be complied with when making the Ceramill ZI zirconium oxide framework with the aid of carbide milling cutters.

Table with 5 columns: Framework thickness min. [mm], Connector cross-section [mm²] (Anterior, Posterior), No. of connected bridge spans max. (Anterior, Posterior). Row: Ceramill ZI, 0,5, ≥7, ≥9, 3, 2.

Detailed min. Wall thicknesses and connector cross-sections

Table with 7 columns: Indication, Diagram, No. of units in total, No. of connected bridge spans, Wall thickness [mm] (Incisal/occlusal, Circular), Connector cross-section [mm²]. Rows include Primary parts/double crowns, Individual crowns, Anterior bridge, Posterior bridge, and Cantilever bridge.

Legend: O = tooth abutment, X = Bridge span

Area of indication and procedure when making frameworks

- Conventional crowns and bridges anchored by crowns in the posterior and anterior area.
- Bridging frameworks with a max. of 3 connecting pontics in the anterior and a max. of 2 connecting pontics in the posterior and a maximum length of 50 mm.
- Cantilever bridges with a max. of 2 pontics (cantilever element max. up to premolar)

Further details of framework design and production can be gathered from the training documents on the Ceramill System and the instruction manual relating to the modelling material Ceramill Gel/Pontic (760514/760522).

- The shrinkage factor of Ceramill ZI blanks as ascertained by Amann Girschbach GmbH is given on the blanks coded according to the batch. This has to be regulated on the Ceramill Base milling system or input into the CAD software of the CAD/CAM system.
- The blanks have to be processed in the way and according to the milling strategy as outlined in the training documents. When used in other systems or CAD/CAM systems, the milling parameters for zirconium oxide and the systems' relevant instructions for use have to be observed.
- According to requirements, the frameworks can be dyed with Ceramill Colour dyeing solutions (LL1 - LL5, Set 760411), either partially or completely (cf. processing instructions of the dyeing solutions).

Final sintering

The final sintering of the framework takes place according to the following program:

- Room temperature -> final temperature 1450 °C (heating-up rate 5 - 10 K/min)
- Hold time at final temperature: 2 hours
- Cooling-down phase up to room temperature (at least < 200 °C): approx. 5 K/min (roughly 5 hours)

Put the Ceramill ZI frameworks for sintering into the sintering bowl (I78310) filled with sintering balls (I78311). Lay the frameworks onto the sintering balls with slight pressure to ensure good support. Take care that no balls get jammed in the interdental spaces or find their way into the crown cavities!

For the sintering, the high-temperature Ceramill Therm (I78300) or another similarly high-quality and suitable furnace is recommended (see also the instructions for use of the sintering furnace in question).

Post-processing

After final sintering, the zirconium oxide frameworks can, if necessary, be post-processed, unpreserved, with a water-cooled laboratory turbine and suitable diamond points.

Preparing frameworks for veneering

After checking fitting accuracy, contact points and occlusion, clean the interior surfaces of the crown by corundum jets (aluminium oxide, 110 µm, pressure ≈ 2 bar) and subsequent evaporation.

Regeneration fire

After processing with diamond points or corundum, we recommend a regeneration fire at 1000 °C (starting temperature 500 °C, heating-up rate 100 K/min, hold time at final temperature 15 min)

Veneering ceramic

Ceramill ZI is suitable for veneering with conventional zirconium-oxide veneering ceramics (e.g. Creation ZI). The thermal expansion coefficient of Ceramill ZI blanks can be gathered from the technical data.

Anchoring

Conventional cementing: Due to the high durability and stability of the Ceramill ZI frameworks, in most cases, conventional anchoring with Zinc-oxide phosphate cement or glass ionomer luting cement (e.g. GC Fuji Plus) is possible.

Adhesive anchoring:

For adhesive anchoring, we recommend the adhesive composite Panavia® 21 or Panavia® F 2.0 due to the excellent adhesive bond it forms with the zirconium oxide framework material.

STORAGE

Store in a dry place in its original packaging.

WARNING AND SAFETY INSTRUCTIONS

When processing, please use personal protective clothing (dust mask, protective goggles, ...).

GUARANTEE/DISCLAIMER CLAUSE

Recommendations regarding the application technique, whether they are given orally, in writing or in the course of practical instruction, should be regarded as guidelines. Our products are subject to ongoing development. We therefore reserve the right to make changes in handling and composition.

TECHNICAL DATA

Table with 3 columns: Property, Value, Unit. Rows include Bending strength (4-point), E-module, Grading, Density, Open porosity, Thermal expansion coefficient (WAK) (25-500 °C), Chem. solubility, Radioactivity.

CHEMICAL COMPOSITION

Table with 2 columns: Component, Value. Rows include in mass percent, ZrO2 + HfO2 + Y2O3, Y2O3, HfO2, Al2O3, Other oxides.

PACKAGE SIZES

Table with 4 columns: Code, Ceramill ZI, Zirconium oxide blanks, Pack of. Rows: 760116, 760143, 760181.

SAFETY DATA SHEETS/DECLARATION OF CONFORMITY

Safety data sheets and declaration of conformity can be obtained on demand from www.amanngirschbach.com

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