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MDCG 2019

Summary of safety and clinical performance

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ISO 6872:2015 + Amd.1:2018

DIN EN ISO 9693:2020

Dentistry - Compatibility testing for metal-ceramic and ceramic-ceramic systems

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| | |
|---|--|
| Reg. EU 2017/745 | Regulation (EU) 2017/745 of the European Parliament and of the council on medical devices |
| MEDDEV 2.7/1 Revision 4 | Clinical Evaluation – A Guidance for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC |
| Medical devices act | Austrian Federal Law on Medical Devices Version 2010 German Medical Devices Act Version 2010 |
| DIN EN 1641:2010 EN 1641:2009 | Dentistry - Medical devices for dentistry - Materials |
| DIN EN ISO 13485:2016 ISO 13485:2016 | Medical devices quality management systems |
| DIN EN ISO 14971:2020 EN ISO 14971:2019 ISO 14971:2019 | Medical devices – Application of risk management to medical devices |
| DIN EN ISO 15223-1:2017 EN ISO 15223-1:2016 ISO 15223-1:2016 | Symbols for use in the labeling of medical devices |
| DIN EN 1041:2013 EN 1041:2008+A1:2013 | Information supplied by the manufacturer of medical devices |
| DIN EN 62366:2017 EN 62366-1:2015 + AC:2015 IEC 62366-1:2015 + COR1:2016 | Medical devices – Part 1: Application of usability engineering to medical devices |
| DIN EN ISO 10993-1:2021 EN ISO 10993-1:2020 ISO 10993-1:2018 | Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process |
| DIN EN ISO 10993-5:2009 EN ISO 10993-5:2009 ISO 10993-5:2009 | Biological evaluation of medical devices; Part 5: Tests for in vitro cytotoxicity |
| DIN EN ISO 10993-12:2012 EN ISO 10993-12:2012 ISO 10993-12:2021 | Biological evaluation of medical devices – Part 12: Sample preparation and reference materials |

Öääë. 9

| | |
|--------------------------|---|
| DIN EN ISO 10993-15:2009 | Biological evaluation of medical devices - Part 15: Identification and quantification of degradation products from metals and alloys |
| EN ISO 10993-15:2009 | |
| ISO 10993-15:2019 | |
| DIN EN ISO 10993-18:2009 | Biological evaluation of medical devices Part 18: Chemical characterization of materials |
| DIN EN ISO 10993-18:2021 | Biological evaluation of medical devices – Part 18: Chemical characterization of medical device materials within a risk management process |
| EN ISO 10993-18:2020 | |
| ISO 10993-18:2020 | |
| DIN EN ISO 7405:2019 | Dentistry – Evaluation of biocompatibility of medical devices used in dentistry |
| EN ISO 7405:2018 | |
| ISO 7405:2018 | |
| DIN EN ISO 10993-10:2014 | Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization |
| EN ISO 10993-10:2013 | |
| ISO 10993-10:2010 | |
| DIN EN ISO 10993-11:2018 | Biological evaluation of medical devices – Part 11: Tests for systemic toxicity |
| EN ISO 10993-11:2018 | |
| ISO 10993-11:2017 | |
| DIN EN ISO 14705:2021-04 | |
| EN ISO 14705:2021 | Fine ceramics (advanced ceramics, advanced technical ceramics) - Test method for hardness of monolithic ceramics at room temperature |
| ISO 14705:2016 | |
| DIN EN ISO 843-1:2008 | Advanced technical ceramics – Mechanical properties of monolithic ceramics at room temperature – Part 1: Determination of flexural strength |
| EN 843-1:2006 | |
| ASTM E1876 - 15 | Standard Test Method for Dynamic Young's Modulus, Shear Modulus, and Poisson's Ratio by Impulse Excitation of Vibration |

Öäää. 9

References

- [1] Dogan S, Raigrodski AJ, Zhang H, Mancl LA. Prospective cohort clinical study assessing the 5-year survival and success of anterior maxillary zirconia-based crowns with customized zirconia copings. *J Prosthet Dent* 2017;117:226–32. <https://doi.org/10.1016/j.prosdent.2016.07.019>.
- [2] Worni A, Katsoulis J, Kolgeci L, Worni M, Mericske-Stern R. Monolithic zirconia reconstructions supported by teeth and implants: 1- to 3-year results of a case series. *Quintessence Int* 2017;48:459–67. <https://doi.org/10.3290/j.qi.a38138>.
- [3] Nejatidanesh F, Moradpoor H, Savabi O. Clinical outcomes of zirconia-based implant- and tooth-supported single crowns. *Clin Oral Investig* 2016;20:169–78. <https://doi.org/10.1007/s00784-015-1479-3>.
- [4] Rinke S, Schafer S, Lange K, Gersdorff N, Roediger M. Practice-based clinical evaluation of metal-ceramic and zirconia molar crowns: 3-year results. *J Oral Rehabil* 2013;40:228–37. <https://doi.org/10.1111/joor.12018>.
- [5] Bomicke W, Rammelsberg P, Stober T, Schmitter M. Short-Term Prospective Clinical Evaluation of Monolithic and Partially Veneered Zirconia Single Crowns. *J Esthet Restor Dent* 2017;29:22–30. <https://doi.org/10.1111/jerd.12270>.
- [6] Tartaglia GM, Sidoti E, Sforza C. Seven-year prospective clinical study on zirconia-based single crowns and fixed dental prostheses. *Clin Oral Investig* 2015;19:1137–45. <https://doi.org/10.1007/s00784-014-1330-2>.
- [7] Ortop A, Kihl ML, Carlsson GE. A 5-year retrospective study of survival of zirconia single crowns fitted in a private clinical setting. *J Dent* 2012;40:527–30. <https://doi.org/10.1016/j.jdent.2012.02.011>.
- [8] Pott P-C, Eisenburger M, Stiesch M. Survival rate of modern all-ceramic FPDs during an observation period from 2011 to 2016. *J Adv Prosthodont* 2018;10:18–24. <https://doi.org/10.4047/jap.2018.10.1.18>.
- [9] Beuer F, Sachs C, Groesser J, Gueth J-F, Stimmelmayer M. Tooth-implant-supported posterior fixed dental prostheses with zirconia frameworks: 3-year clinical result. *Clin Oral Investig* 2016;20:1079–86. <https://doi.org/10.1007/s00784-015-1594-1>.
- [10] Guncu MB, Cakan U, Aktas G, Guncu GN, Canay S. Comparison of implant versus tooth-supported zirconia-based single crowns in a split-mouth design: a 4-year clinical follow-up study. *Clin Oral Investig* 2016;20:2467–73. <https://doi.org/10.1007/s00784-016-1763-x>.
- [11] Rinke S, Lange K, Roediger M, Gersdorff N. Risk factors for technical and biological complications with zirconia single crowns. *Clin Oral Investig* 2015;19:1999–2006. <https://doi.org/10.1007/s00784-015-1410-y>.
- [12] Koenig V, Vanheusden AJ, Le Goff SO, Mainjot AK. Clinical risk factors related to failures with zirconia-based restorations: an up to 9-year retrospective study. *J Dent* 2013;41:1164–74. <https://doi.org/10.1016/j.jdent.2013.10.009>.
- [13] Vigolo P, Mutinelli S. Evaluation of zirconium-oxide-based ceramic single-unit posterior fixed dental prostheses (FDPs) generated with two CAD/CAM systems compared to porcelain-fused-to-metal single-unit posterior FDPs: a 5-year clinical prospective study. *J Prosthodont* 2012;21:265–9. <https://doi.org/10.1111/j.1532-849X.2011.00825.x>.

- [14] Monaco C, Llukacej A, Baldissara P, Arena A, Scotti R. Zirconia-based versus metal-based single crowns veneered with overpressing ceramic for restoration of posterior endodontically treated teeth: 5-year results of a randomized controlled clinical study. *J Dent* 2017;65:56–63. <https://doi.org/10.1016/j.jdent.2017.07.004>.
- [15] Burke FJT, Crisp RJ, Cowan AJ, Lamb J, Thompson O, Tulloch N. Five-year clinical evaluation of zirconia-based bridges in patients in UK general dental practices. *J Dent* 2013;41:992–9. <https://doi.org/10.1016/j.jdent.2013.08.007>.
- [16] Ohlmann B, Eiffler C, Rammelsberg P. Clinical performance of all-ceramic cantilever fixed dental prostheses: results of a 2-year randomized pilot study. *Quintessence Int* 2012;43:643–8.
- [17] Pelaez J, Cogolludo PG, Serrano B, Serrano JFL, Suarez MJ. A four-year prospective clinical evaluation of zirconia and metal-ceramic posterior fixed dental prostheses. *Int J Prosthodont* 2012;25:451–8.
- [18] Passia N, Chaar MS, Kern M. Outcome of posterior fixed dental prostheses made from veneered zirconia over an observation period of up to 13 years. *J Dent* 2019;86:126–9. <https://doi.org/10.1016/j.jdent.2019.05.032>.
- [19] Chaar MS, Kern M. Five-year clinical outcome of posterior zirconia ceramic inlay-retained FDPs with a modified design. *J Dent* 2015;43:1411–5. <https://doi.org/10.1016/j.jdent.2015.11.001>.
- [20] Zenthofer A, Ohlmann B, Rammelsberg P, Bomicke W. Performance of zirconia ceramic cantilever fixed dental prostheses: 3-year results from a prospective, randomized, controlled pilot study. *J Prosthet Dent* 2015;114:34–9. <https://doi.org/10.1016/j.prosdent.2015.02.006>.
- [21] Sailer I, Balmer M, Husler J, Hammerle CHF, Kanel S, Thoma DS. 10-year randomized trial (RCT) of zirconia-ceramic and metal-ceramic fixed dental prostheses. *J Dent* 2018;76:32–9. <https://doi.org/10.1016/j.jdent.2018.05.015>.
- [22] Suarez MJ, Perez C, Pelaez J, Lopez-Suarez C, Gonzalo E. A Randomized Clinical Trial Comparing Zirconia and Metal-Ceramic Three-Unit Posterior Fixed Partial Dentures: A 5-Year Follow-Up. *J Prosthodont* 2018. <https://doi.org/10.1111/jopr.12952>.
- [23] Ioannidis A, Bindl A. Clinical prospective evaluation of zirconia-based three-unit posterior fixed dental prostheses: Up-to ten-year results. *J Dent* 2016;47:80–5. <https://doi.org/10.1016/j.jdent.2016.01.014>.
- [24] Larsson C, Vult von Steyern P. Ten-Year Follow-Up of Implant-Supported All-Ceramic Fixed Dental Prostheses: A Randomized, Prospective Clinical Trial. *Int J Prosthodont* 2016;29:31–4. <https://doi.org/10.11607/ijp.4328>.
- [25] Naenni N, Bindl A, Sax C, Hammerle C, Sailer I. A randomized controlled clinical trial of 3-unit posterior zirconia-ceramic fixed dental prostheses (FDP) with layered or pressed veneering ceramics: 3-year results. *J Dent* 2015;43:1365–70. <https://doi.org/10.1016/j.jdent.2015.07.013>.
- [26] Spies BC, Balmer M, Jung RE, Sailer I, Vach K, Kohal R-J. All-ceramic single crowns supported by zirconia implants: 5-year results of a prospective multicenter study. *Clin Oral Implants Res* 2019;30:466–75. <https://doi.org/10.1111/clr.13433>.
- [27] Rathmann F, Bomicke W, Rammelsberg P, Ohlmann B. Veneered zirconia inlay-retained fixed dental prostheses: 10-Year results from a prospective clinical study. *J Dent* 2017;64:68–72. <https://doi.org/10.1016/j.jdent.2017.06.008>.

- [28] Edelhoff D, Werder U. Komplexe Rehabilitation aus hochtransluzentem Zirkonoxid. Restaurationen von der Testphase bis zur definitiven Versorgung; 2020 Dec.
- [29] Haralur SB, Alamrey AA, Alshehri SA, Alzahrani DS, Alfarsi M. Effect of different preparation designs and all ceramic materials on fracture strength of molar endocrowns. *J Appl Biomater Funct Mater* 2020;18:2280800020947329. <https://doi.org/10.1177/2280800020947329>.
- [30] Juntavee N, Kornrum S. Effect of Marginal Designs on Fracture Strength of High Translucency Monolithic Zirconia Crowns. *Int J Dent* 2020;2020:8875609. <https://doi.org/10.1155/2020/8875609>.
- [31] Jansen JU, Lümckemann N, Sener B, Stawarczyk B. Comparison of fracture toughness measurements for zirconia materials using two test methods. *Dental materials journal* 2019;38:806–12. <https://doi.org/10.4012/dmj.2018-361>.
- [32] Almansour HM, Alqahtani F. The Effect of in vitro Aging and Fatigue on the Flexural Strength of Monolithic High-translucency Zirconia Restorations. *J Contemp Dent Pract* 2018;19:867–73.
- [33] Stawarczyk B, Ozcan M, Hallmann L, Ender A, Mehl A, Hämmerlet CHF. The effect of zirconia sintering temperature on flexural strength, grain size, and contrast ratio. *Clin Oral Investig* 2013;17:269–74. <https://doi.org/10.1007/s00784-012-0692-6>.
- [34] Amann Girrbach R&D. Test report 3P bending strength Ceramill ZI; 2015.
- [35] Amann Girrbach R&D. Test report 4P bending strength Ceramill ZI; 2015.
- [36] LMU Munich. Test report 4P bending strength aged Ceramill ZI; 2012.
- [37] Amann Girrbach R&D. Test report 3P bending strength Ceramill Zolid HT+; 2017.
- [38] Amann Girrbach R&D. Test report 4P bending strength Ceramill Zolid HT+; 2017.
- [39] Amann Girrbach R&D. Test report 4P bending strength Ceramill Zolid HT+ aged; 2017.
- [40] Amann Girrbach R&D. Test report 3P bending strength Zolid HT+ preshade A4; 2017.
- [41] Amann Girrbach R&D. Test report 4P bending strength Zolid HT+ preshade A4; 2017.
- [42] Amann Girrbach R&D. Test report 3P bending strength Zolid gen-x A3.5; 2019.
- [43] Amann Girrbach R&D. Test report 4P bending strength Zolid gen-x A3.5; 2019.
- [44] Amann Girrbach R&D. Test report bending strength (3-point) Zolid DRS B4.
- [45] Amann Girrbach R&D. Test report bending strength (4-point) Zolid DRS B4.
- [46] Amann Girrbach R&D. Test report 4P bending strength Zolid HT+ PS A4 aged; 2019.
- [47] Bioserv - Analytics and medical products GmbH. Test report animal studies toxicological and biological safety Ceramill ZI; 2010.
- [48] Bioserv - Analytics and medical products GmbH. Test report cytotoxicity Ceramill ZI; 2007.
- [49] Research Centre Jülich. Test report radioactivity Ceramill ZI; 2007.
- [50] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill ZI; 2017.
- [51] Medical Device Services. Test report cytotoxicity Ceramill Zolid HT+; 2017.
- [52] Medical Device Services. Test report chemical analysis (organic extractables) Ceramill Zolid HT+; 2017.
- [53] Medical Device Services. Test report chemical analysis (inorganic leachables) Ceramill Zolid HT+; 2017.

- [54] Research Centre Jülich. Test report radioactivity Ceramill Zolid HT+; 2016.
- [55] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill Zolid HT+; 2017.
- [56] Medical Device Services. Test report cytotoxicity Zpex 4 Yellow; 2017.
- [57] Medical Device Services. Test report cytotoxicity Zpex Gray; 2017.
- [58] Medical Device Services. Test report cytotoxicity Zpex Yellow and Zpex Pink; 2013.
- [59] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Gray; 2017.
- [60] Medical Device Services. Test report chemical analysis (organic extractables) Zpex4 Yellow; 2017.
- [61] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Pink; 2013.
- [62] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex 4 Yellow; 2017.
- [63] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Gray; 2017.
- [64] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Pink; 2013.
- [65] Research Centre Jülich. Test report radioactivity Zpex 4 Yellow; 2016.
- [66] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex Pink; 2017.
- [67] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex 4 Yellow; 2017.
- [68] Medical Device Services. Biological safety Ceramill Zolid HT+; 2017.
- [69] Medical Device Services. Biological safety Zpex 4 Yellow; 2017.
- [70] Medical Device Services. Biological safety Zpex Gray; 2017.
- [71] Medical Device Services. Certificate biological safety Ceramill Zolid HT+; 2017.
- [72] Medical Device Services. Certificate biological safety Zpex 4 Yellow; 2017.
- [73] Medical Device Services. Certificate biological safety Zpex Gray; 2017.
- [74] Ferrari M, Cagidiaco MC, Garcia-Godoy F, Goracci C, Cairo F. Effect of different prosthetic abutments // Effect of different prosthetic abutments on peri-implant soft tissue. A randomized controlled clinical trial. *Am J Dent* 2015;28:85–9.
- [75] Stober T, Bermejo JL, Schwindling FS, Schmitter M. Clinical assessment of enamel wear caused by monolithic zirconia crowns. *J Oral Rehabil* 2016;43:621–9. <https://doi.org/10.1111/joor.12409>.
- [76] Al Hamad KQ, Abu Al-Addous AM, Al-Wahadni AM, Baba NZ, Goodacre BJ. Surface Roughness of Monolithic and Layered Zirconia Restorations at Different Stages of Finishing and Polishing: An In Vitro Study. *J Prosthodont* 2019;28:818–25. <https://doi.org/10.1111/jopr.13071>.
- [77] Saeed F, Muhammad N, Khan AS, Sharif F, Rahim A, Ahmad P et al. Prosthodontics dental materials: From conventional to unconventional. *Mater Sci Eng C Mater Biol Appl* 2020;106:110167. <https://doi.org/10.1016/j.msec.2019.110167>.
- [78] Farawati FA, Nakaparksin P. What is the Optimal Material for Implant Prosthesis? *Dent Clin North Am* 2019;63:515–30. <https://doi.org/10.1016/j.cden.2019.02.002>.
- [79] Ferrari M, Vichi A, Zarone F. Zirconia abutments and restorations: from laboratory to clinical investigations. *Dent Mater* 2015;31:e63-76. <https://doi.org/10.1016/j.dental.2014.11.015>.

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG

Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
germany@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

P002933878 27.05.2021 Version 2.0

AMANNGIRRBACH