 **ceramill**[®] a-temp
 **ceramill**[®] a-temp multilayer

HR Sažetak o sigurnosnoj i kliničkoj učinkovitosti proizvoda (SSCP)

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MDCG 2019 Summary of safety and clinical performance

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Reg. EU 2017/745 Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices

MEDDEV 2.7/1 Revision 4 Clinical Evaluation – A Guidance for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC

Medical devices act Austrian Federal Law on Medical Devices Version 2010 German Medical Devices Act Version 2010

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DIN EN 1641:2010 Dentistry - Medical devices for dentistry - Materials

EN 1641:2009

DIN EN ISO 13485:2016 Medical devices quality management systems

ISO 13485:2016

DIN EN ISO 14971:2020 Medical devices – Application of risk management to medical devices

EN ISO 14971:2019

ISO 14971:2019

DIN EN ISO 15223-1:2017 Symbols for use in the labeling of medical devices

EN ISO 15223-1:2016

ISO 15223-1:2016

DIN EN 1041:2013 Information supplied by the manufacturer of medical devices

EN 1041:2008+A1:2013

DIN EN 62366:2017 Medical devices – Part 1: Application of usability engineering to medical devices

EN 62366-1:2015 + AC:2015

IEC 62366-1:2015 +

COR1:2016

tab. 8

DIN EN ISO 20795-2:2013	Dentistry – Base polymers; Part 2: Orthodontic base polymers
EN ISO 20795-2:2013	
ISO 20795-2:2013	
DIN EN ISO 10993-1:2021	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
EN ISO 10993-1:2020	
ISO 10993-1:2018	
DIN EN ISO 10993-5:2009	Biological evaluation of medical devices; Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-5:2009	
ISO 10993-5:2009	
DIN EN ISO 10993-12:2012	Biological evaluation of medical devices – Part 12: Sample preparation and reference materials
EN ISO 10993-12:2012	
ISO 10993-12:2021	
DIN EN ISO 10993-15:2009	Biological evaluation of medical devices - Part 15: Identification and quantification of degradation products from metals and alloys
EN ISO 10993-15:2009	
ISO 10993-15:2019	
DIN EN ISO 10993-18:2021	Biological evaluation of medical devices – Part 18: Chemical characterization of medical device materials within a risk management process
EN ISO 10993-18:2020	
ISO 10993-18:2020	
DIN EN ISO 7405:2019	Dentistry – Evaluation of biocompatibility of medical devices used in dentistry
EN ISO 7405:2018	
ISO 7405:2018	
DIN EN ISO 10477:2021	Dentistry - Polymer-based crown and veneering materials
EN ISO 10477:2020	
ISO 10477:2020	

tab. 8

References

- [1] Syed M, Chopra R, Sachdev V. Allergic Reactions to Dental Materials-A Systematic Review. *J Clin Diagn Res* 2015;9:ZE04-9. <https://doi.org/10.7860/JCDR/2015/15640.6589>.
- [2] Souza IR, Pansani TN, Basso FG, Hebling J, Souza Costa CA de. Cytotoxicity of acrylic resin-based materials used to fabricate interim crowns. *J Prosthet Dent* 2020. <https://doi.org/10.1016/j.prosdent.2020.01.030>.
- [3] Carvalho FR, Lentini-Oliveira DA, Prado LBF, Prado GF, Carvalho LBC. Oral appliances and functional orthopaedic appliances for obstructive sleep apnoea in children. *Cochrane Database of Systematic Reviews* 2016.
- [4] Lee J-H, Jun S-K, Moon H-J, Lee H-H. Cytotoxicity and proinflammatory cytokine expression induced by interim resin materials in primary cultured human dental pulp cells. *J Prosthet Dent* 2017;118:524–34. <https://doi.org/10.1016/j.prosdent.2016.11.014>.
- [5] Huettig F, Prutscher A, Goldammer C, Kreutzer CA, Weber H. First clinical experiences with CAD/CAM-fabricated PMMA-based fixed dental prostheses as long-term temporaries. *Clin Oral Investig* 2016;20:161–8. <https://doi.org/10.1007/s00784-015-1475-7>.
- [6] Sari T, Usumez A, Strasser T, Şahinbas A, Rosentritt M. Temporary materials: comparison of in vivo and in vitro performance. *Clin Oral Investig* 2020. <https://doi.org/10.1007/s00784-020-03278-5>.
- [7] Edelhoff D, Beuer F, Schweiger J, Brix O, Stimmelmayer M, Guth J-F. CAD/CAM-generated high-density polymer restorations for the pretreatment of complex cases: a case report. *Quintessence Int* 2012;43:457–67.
- [8] Rosentritt M, Raab P, Hahnel S, Stockle M, Preis V. In-vitro performance of CAD/CAM-fabricated implant-supported temporary crowns. *Clin Oral Investig* 2017;21:2581–7. <https://doi.org/10.1007/s00784-017-2057-7>.
- [9] Digholkar S, Madhav VNV, Palaskar J. Evaluation of the flexural strength and microhardness of provisional crown and bridge materials fabricated by different methods. *J Indian Prosthodont Soc* 2016;16:328–34. <https://doi.org/10.4103/0972-4052.191288>.
- [10] Rayyan MM, Aboushelib M, Sayed NM, Ibrahim A, Jimbo R. Comparison of interim restorations fabricated by CAD/CAM with those fabricated manually. *J Prosthet Dent* 2015;114:414–9. <https://doi.org/10.1016/j.prosdent.2015.03.007>.
- [11] Penate L, Basilio J, Roig M, Mercade M. Comparative study of interim materials for direct fixed dental prostheses and their fabrication with CAD/CAM technique. *J Prosthet Dent* 2015;114:248–53. <https://doi.org/10.1016/j.prosdent.2014.12.023>.
- [12] Yao J, Li J, Wang Y, Huang H. Comparison of the flexural strength and marginal accuracy of traditional and CAD/CAM interim materials before and after thermal cycling. *J Prosthet Dent* 2014;112:649–57. <https://doi.org/10.1016/j.prosdent.2014.01.012>.
- [13] Alt V, Hannig M, Wostmann B, Balkenhol M. Fracture strength of temporary fixed partial dentures: CAD/CAM versus directly fabricated restorations. *Dent Mater* 2011;27:339–47. <https://doi.org/10.1016/j.dental.2010.11.012>.
- [14] Elagra MI, Rayyan MR, Alhomaidhi MM, Alanazy AA, Ainefaie MO. Color stability and marginal integrity of interim crowns: An in vitro study. *Eur J Dent* 2017;11:330–4. https://doi.org/10.4103/ejd.ejd_66_17.
- [15] Amann Girrbaach R&D. 20180731 Test report water absorption and solubility Ceramill A-Temp A3,5; 2018.

- [16] Amann Girrbach R&D. 20180724 Test report bending strength Ceramill A-Temp A3,5; 2018.
- [17] Amann Girrbach R&D. 20190213 Test report wearing period Ceramill A-Temp A-Splint; 2019.
- [18] Nordic Institute of Dental Materials (NIOM). Test report residual monomer Ceramill A-Splint; 2018.
- [19] Medical Device Services. Test report cytotoxicity Ceramill A-Temp; 2018.
- [20] Medical Device Services. Test report chemical analysis inorganic leachables Ceramill A-Temp; 2018.
- [21] Medical Device Services. Test report chemical analysis organic leachables extractables Ceramill A-Temp; 2018.
- [22] Medical Device Services. Certificate biological safety Ceramill A-Temp; 2018.
- [23] Medical Device Services. Biological safety Ceramill A-Temp; 2018.
- [24] Arslan M, Murat S, Alp G, Zaimoglu A. Evaluation of flexural strength and surface properties of prepolymerized CAD/CAM PMMA-based polymers used for digital 3D complete dentures. *Int J Comput Dent* 2018;21:31–40.
- [25] Prpić V, Schauerl Z, Čatić A, Dulčić N, Čimić S. Comparison of Mechanical Properties of 3D-Printed, CAD/CAM, and Conventional Denture Base Materials. *J Prosthodont* 2020. <https://doi.org/10.1111/jopr.13175>.
- [26] Alp G, Murat S, Yilmaz B. Comparison of Flexural Strength of Different CAD/CAM PMMA-Based Polymers. *J Prosthodont* 2018. <https://doi.org/10.1111/jopr.12755>.
- [27] Ayman A-D. The residual monomer content and mechanical properties of CAD/CAM resins used in the fabrication of complete dentures as compared to heat cured resins. *Electron Physician* 2017;9:4766–72. <https://doi.org/10.19082/4766>.
- [28] Stawarczyk B, Ender A, Trottmann A, Özcan M, Fischer J, Hämmerle CHF. Load-bearing capacity of CAD/CAM milled polymeric three-unit fixed dental prostheses. Effect of aging regimens. *Clin Oral Investig* 2012;16:1669–77. <https://doi.org/10.1007/s00784-011-0670-4>.
- [29] Steinmassl P-A, Wiedemair V, Huck C, Klauzner F, Steinmassl O, Grunert I et al. Do CAD/CAM dentures really release less monomer than conventional dentures? *Clin Oral Investig* 2017;21:1697–705. <https://doi.org/10.1007/s00784-016-1961-6>.

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG

Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
germany@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

P002920781 28.05.2021 Version 1.0

AMANNGIRRBACH