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ne ustreza ni na voljo popolna uporaba delna uporaba

MDCG 2019	Summary of safety and clinical performance
-----------	--

Tab. 7

Uporabljene harmonizirane norme (če ustreza):

ne ustreza ni na voljo popolna uporaba delna uporaba

DIN EN ISO 6872:2019	Dentistry – Ceramic materials
----------------------	-------------------------------

EN ISO 6872:2015 + A1:2018

ISO 6872:2015 + Amd.1:2018

DIN EN ISO 9693:2020	Dentistry - Compatibility testing for metal-ceramic and ceramic-ceramic systems
----------------------	---

EN ISO 9693:2019

ISO 9693:2019

Tab. 8

(Izpolnjujemo normo DIN EN ISO 6872 za navedene točke za keramike tipa II, razreda 5. S tem se ne izvajajo kontrole enakomernosti in temperature steklastega prehoda, saj za materiale iz cirkonijevega oksida te niso potrebne).

ne ustreza ni na voljo popolna uporaba delna uporaba

Reg. EU 2017/745	Regulation (EU) 2017/745 of the European Parliament and of the council on medical devices
MEDDEV 2.7/1 Revision 4	Clinical Evaluation – A Guidance for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC
Medical devices act	Austrian Federal Law on Medical Devices Version 2010 German Medical Devices Act Version 2010
DIN EN 1641:2010 EN 1641:2009	Dentistry - Medical devices for dentistry - Materials
DIN EN ISO 13485:2016 ISO 13485:2016	Medical devices quality management systems
DIN EN ISO 14971:2020 EN ISO 14971:2019 ISO 14971:2019	Medical devices – Application of risk management to medical devices
DIN EN ISO 15223-1:2017 EN ISO 15223-1:2016 ISO 15223-1:2016	Symbols for use in the labeling of medical devices
DIN EN 1041:2013 EN 1041:2008+A1:2013	Information supplied by the manufacturer of medical devices
DIN EN 62366:2017 EN 62366-1:2015 + AC:2015 IEC 62366-1:2015 + COR1:2016	Medical devices – Part 1: Application of usability engineering to medical devices
DIN EN ISO 10993-1:2021 EN ISO 10993-1:2020 ISO 10993-1:2018	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
DIN EN ISO 10993-5:2009 EN ISO 10993-5:2009 ISO 10993-5:2009	Biological evaluation of medical devices; Part 5: Tests for in vitro cytotoxicity
DIN EN ISO 10993-12:2012 EN ISO 10993-12:2012 ISO 10993-12:2021	Biological evaluation of medical devices – Part 12: Sample preparation and reference materials

Tab. 9

DIN EN ISO 10993-15:2009	Biological evaluation of medical devices - Part 15: Identification and quantification of degradation products from metals and alloys
EN ISO 10993-15:2009	
ISO 10993-15:2019	
DIN EN ISO 10993-18:2009	Biological evaluation of medical devices Part 18: Chemical characterization of materials
DIN EN ISO 10993-18:2021	Biological evaluation of medical devices – Part 18: Chemical characterization of medical device materials within a risk management process
EN ISO 10993-18:2020	
ISO 10993-18:2020	
DIN EN ISO 7405:2019	Dentistry – Evaluation of biocompatibility of medical devices used in dentistry
EN ISO 7405:2018	
ISO 7405:2018	
DIN EN ISO 10993-10:2014	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization
EN ISO 10993-10:2013	
ISO 10993-10:2010	
DIN EN ISO 10993-11:2018	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
EN ISO 10993-11:2018	
ISO 10993-11:2017	
DIN EN ISO 14705:2021-04	
EN ISO 14705:2021	Fine ceramics (advanced ceramics, advanced technical ceramics) - Test method for hardness of monolithic ceramics at room temperature
ISO 14705:2016	
DIN EN ISO 843-1:2008	Advanced technical ceramics – Mechanical properties of monolithic ceramics at room temperature – Part 1: Determination of flexural strength
EN 843-1:2006	
ASTM E1876 - 15	Standard Test Method for Dynamic Young's Modulus, Shear Modulus, and Poisson's Ratio by Impulse Excitation of Vibration

Tab. 9

References

- [1] Dogan S, Raigrodski AJ, Zhang H, Mancl LA. Prospective cohort clinical study assessing the 5-year survival and success of anterior maxillary zirconia-based crowns with customized zirconia copings. *J Prosthet Dent* 2017;117:226–32. <https://doi.org/10.1016/j.prosdent.2016.07.019>.
- [2] Worni A, Katsoulis J, Kolgeci L, Worni M, Mericske-Stern R. Monolithic zirconia reconstructions supported by teeth and implants: 1- to 3-year results of a case series. *Quintessence Int* 2017;48:459–67. <https://doi.org/10.3290/j.qi.a38138>.
- [3] Nejatidanesh F, Moradpoor H, Savabi O. Clinical outcomes of zirconia-based implant- and tooth-supported single crowns. *Clin Oral Investig* 2016;20:169–78. <https://doi.org/10.1007/s00784-015-1479-3>.
- [4] Rinke S, Schafer S, Lange K, Gersdorff N, Roediger M. Practice-based clinical evaluation of metal-ceramic and zirconia molar crowns: 3-year results. *J Oral Rehabil* 2013;40:228–37. <https://doi.org/10.1111/joor.12018>.
- [5] Bomicke W, Rammelsberg P, Stober T, Schmitter M. Short-Term Prospective Clinical Evaluation of Monolithic and Partially Veneered Zirconia Single Crowns. *J Esthet Restor Dent* 2017;29:22–30. <https://doi.org/10.1111/jerd.12270>.
- [6] Tartaglia GM, Sidoti E, Sforza C. Seven-year prospective clinical study on zirconia-based single crowns and fixed dental prostheses. *Clin Oral Investig* 2015;19:1137–45. <https://doi.org/10.1007/s00784-014-1330-2>.
- [7] Ortop A, Kihl ML, Carlsson GE. A 5-year retrospective study of survival of zirconia single crowns fitted in a private clinical setting. *J Dent* 2012;40:527–30. <https://doi.org/10.1016/j.jdent.2012.02.011>.
- [8] Pott P-C, Eisenburger M, Stiesch M. Survival rate of modern all-ceramic FPDs during an observation period from 2011 to 2016. *J Adv Prosthodont* 2018;10:18–24. <https://doi.org/10.4047/jap.2018.10.1.18>.
- [9] Beuer F, Sachs C, Groesser J, Gueth J-F, Stimmelmayer M. Tooth-implant-supported posterior fixed dental prostheses with zirconia frameworks: 3-year clinical result. *Clin Oral Investig* 2016;20:1079–86. <https://doi.org/10.1007/s00784-015-1594-1>.
- [10] Guncu MB, Cakan U, Aktas G, Guncu GN, Canay S. Comparison of implant versus tooth-supported zirconia-based single crowns in a split-mouth design: a 4-year clinical follow-up study. *Clin Oral Investig* 2016;20:2467–73. <https://doi.org/10.1007/s00784-016-1763-x>.
- [11] Rinke S, Lange K, Roediger M, Gersdorff N. Risk factors for technical and biological complications with zirconia single crowns. *Clin Oral Investig* 2015;19:1999–2006. <https://doi.org/10.1007/s00784-015-1410-y>.
- [12] Koenig V, Vanheusden AJ, Le Goff SO, Mainjot AK. Clinical risk factors related to failures with zirconia-based restorations: an up to 9-year retrospective study. *J Dent* 2013;41:1164–74. <https://doi.org/10.1016/j.jdent.2013.10.009>.
- [13] Vigolo P, Mutinelli S. Evaluation of zirconium-oxide-based ceramic single-unit posterior fixed dental prostheses (FDPs) generated with two CAD/CAM systems compared to porcelain-fused-to-metal single-unit posterior FDPs: a 5-year clinical prospective study. *J Prosthodont* 2012;21:265–9. <https://doi.org/10.1111/j.1532-849X.2011.00825.x>.

- [14] Monaco C, Llukacej A, Baldissara P, Arena A, Scotti R. Zirconia-based versus metal-based single crowns veneered with overpressing ceramic for restoration of posterior endodontically treated teeth: 5-year results of a randomized controlled clinical study. *J Dent* 2017;65:56–63. <https://doi.org/10.1016/j.jdent.2017.07.004>.
- [15] Burke FJT, Crisp RJ, Cowan AJ, Lamb J, Thompson O, Tulloch N. Five-year clinical evaluation of zirconia-based bridges in patients in UK general dental practices. *J Dent* 2013;41:992–9. <https://doi.org/10.1016/j.jdent.2013.08.007>.
- [16] Ohlmann B, Eiffler C, Rammelsberg P. Clinical performance of all-ceramic cantilever fixed dental prostheses: results of a 2-year randomized pilot study. *Quintessence Int* 2012;43:643–8.
- [17] Pelaez J, Cogolludo PG, Serrano B, Serrano JFL, Suarez MJ. A four-year prospective clinical evaluation of zirconia and metal-ceramic posterior fixed dental prostheses. *Int J Prosthodont* 2012;25:451–8.
- [18] Passia N, Chaar MS, Kern M. Outcome of posterior fixed dental prostheses made from veneered zirconia over an observation period of up to 13 years. *J Dent* 2019;86:126–9. <https://doi.org/10.1016/j.jdent.2019.05.032>.
- [19] Chaar MS, Kern M. Five-year clinical outcome of posterior zirconia ceramic inlay-retained FDPs with a modified design. *J Dent* 2015;43:1411–5. <https://doi.org/10.1016/j.jdent.2015.11.001>.
- [20] Zenthofer A, Ohlmann B, Rammelsberg P, Bomicke W. Performance of zirconia ceramic cantilever fixed dental prostheses: 3-year results from a prospective, randomized, controlled pilot study. *J Prosthet Dent* 2015;114:34–9. <https://doi.org/10.1016/j.prosdent.2015.02.006>.
- [21] Sailer I, Balmer M, Husler J, Hammerle CHF, Kanel S, Thoma DS. 10-year randomized trial (RCT) of zirconia-ceramic and metal-ceramic fixed dental prostheses. *J Dent* 2018;76:32–9. <https://doi.org/10.1016/j.jdent.2018.05.015>.
- [22] Suarez MJ, Perez C, Pelaez J, Lopez-Suarez C, Gonzalo E. A Randomized Clinical Trial Comparing Zirconia and Metal-Ceramic Three-Unit Posterior Fixed Partial Dentures: A 5-Year Follow-Up. *J Prosthodont* 2018. <https://doi.org/10.1111/jopr.12952>.
- [23] Ioannidis A, Bindl A. Clinical prospective evaluation of zirconia-based three-unit posterior fixed dental prostheses: Up-to ten-year results. *J Dent* 2016;47:80–5. <https://doi.org/10.1016/j.jdent.2016.01.014>.
- [24] Larsson C, Vult von Steyern P. Ten-Year Follow-Up of Implant-Supported All-Ceramic Fixed Dental Prostheses: A Randomized, Prospective Clinical Trial. *Int J Prosthodont* 2016;29:31–4. <https://doi.org/10.11607/ijp.4328>.
- [25] Naenni N, Bindl A, Sax C, Hammerle C, Sailer I. A randomized controlled clinical trial of 3-unit posterior zirconia-ceramic fixed dental prostheses (FDP) with layered or pressed veneering ceramics: 3-year results. *J Dent* 2015;43:1365–70. <https://doi.org/10.1016/j.jdent.2015.07.013>.
- [26] Spies BC, Balmer M, Jung RE, Sailer I, Vach K, Kohal R-J. All-ceramic single crowns supported by zirconia implants: 5-year results of a prospective multicenter study. *Clin Oral Implants Res* 2019;30:466–75. <https://doi.org/10.1111/clr.13433>.
- [27] Rathmann F, Bomicke W, Rammelsberg P, Ohlmann B. Veneered zirconia inlay-retained fixed dental prostheses: 10-Year results from a prospective clinical study. *J Dent* 2017;64:68–72. <https://doi.org/10.1016/j.jdent.2017.06.008>.

- [28] Edelhoff D, Werder U. Komplexe Rehabilitation aus hochtransluzentem Zirkonoxid. Restaurationen von der Testphase bis zur definitiven Versorgung; 2020 Dec.
- [29] Haralur SB, Alamrey AA, Alshehri SA, Alzahrani DS, Alfarsi M. Effect of different preparation designs and all ceramic materials on fracture strength of molar endocrowns. *J Appl Biomater Funct Mater* 2020;18:2280800020947329. <https://doi.org/10.1177/2280800020947329>.
- [30] Juntavee N, Kornrum S. Effect of Marginal Designs on Fracture Strength of High Translucency Monolithic Zirconia Crowns. *Int J Dent* 2020;2020:8875609. <https://doi.org/10.1155/2020/8875609>.
- [31] Jansen JU, Lümckemann N, Sener B, Stawarczyk B. Comparison of fracture toughness measurements for zirconia materials using two test methods. *Dental materials journal* 2019;38:806–12. <https://doi.org/10.4012/dmj.2018-361>.
- [32] Almansour HM, Alqahtani F. The Effect of in vitro Aging and Fatigue on the Flexural Strength of Monolithic High-translucency Zirconia Restorations. *J Contemp Dent Pract* 2018;19:867–73.
- [33] Stawarczyk B, Ozcan M, Hallmann L, Ender A, Mehl A, Hämmerlet CHF. The effect of zirconia sintering temperature on flexural strength, grain size, and contrast ratio. *Clin Oral Investig* 2013;17:269–74. <https://doi.org/10.1007/s00784-012-0692-6>.
- [34] Amann Girrbach R&D. Test report 3P bending strength Ceramill ZI; 2015.
- [35] Amann Girrbach R&D. Test report 4P bending strength Ceramill ZI; 2015.
- [36] LMU Munich. Test report 4P bending strength aged Ceramill ZI; 2012.
- [37] Amann Girrbach R&D. Test report 3P bending strength Ceramill Zolid HT+; 2017.
- [38] Amann Girrbach R&D. Test report 4P bending strength Ceramill Zolid HT+; 2017.
- [39] Amann Girrbach R&D. Test report 4P bending strength Ceramill Zolid HT+ aged; 2017.
- [40] Amann Girrbach R&D. Test report 3P bending strength Zolid HT+ preshade A4; 2017.
- [41] Amann Girrbach R&D. Test report 4P bending strength Zolid HT+ preshade A4; 2017.
- [42] Amann Girrbach R&D. Test report 3P bending strength Zolid gen-x A3.5; 2019.
- [43] Amann Girrbach R&D. Test report 4P bending strength Zolid gen-x A3.5; 2019.
- [44] Amann Girrbach R&D. Test report bending strength (3-point) Zolid DRS B4.
- [45] Amann Girrbach R&D. Test report bending strength (4-point) Zolid DRS B4.
- [46] Amann Girrbach R&D. Test report 4P bending strength Zolid HT+ PS A4 aged; 2019.
- [47] Bioserv - Analytics and medical products GmbH. Test report animal studies toxicological and biological safety Ceramill ZI; 2010.
- [48] Bioserv - Analytics and medical products GmbH. Test report cytotoxicity Ceramill ZI; 2007.
- [49] Research Centre Jülich. Test report radioactivity Ceramill ZI; 2007.
- [50] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill ZI; 2017.
- [51] Medical Device Services. Test report cytotoxicity Ceramill Zolid HT+; 2017.
- [52] Medical Device Services. Test report chemical analysis (organic extractables) Ceramill Zolid HT+; 2017.
- [53] Medical Device Services. Test report chemical analysis (inorganic leachables) Ceramill Zolid HT+; 2017.

- [54] Research Centre Jülich. Test report radioactivity Ceramill Zolid HT+; 2016.
- [55] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill Zolid HT+; 2017.
- [56] Medical Device Services. Test report cytotoxicity Zpex 4 Yellow; 2017.
- [57] Medical Device Services. Test report cytotoxicity Zpex Gray; 2017.
- [58] Medical Device Services. Test report cytotoxicity Zpex Yellow and Zpex Pink; 2013.
- [59] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Gray; 2017.
- [60] Medical Device Services. Test report chemical analysis (organic extractables) Zpex4 Yellow; 2017.
- [61] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Pink; 2013.
- [62] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex 4 Yellow; 2017.
- [63] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Gray; 2017.
- [64] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Pink; 2013.
- [65] Research Centre Jülich. Test report radioactivity Zpex 4 Yellow; 2016.
- [66] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex Pink; 2017.
- [67] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex 4 Yellow; 2017.
- [68] Medical Device Services. Biological safety Ceramill Zolid HT+; 2017.
- [69] Medical Device Services. Biological safety Zpex 4 Yellow; 2017.
- [70] Medical Device Services. Biological safety Zpex Gray; 2017.
- [71] Medical Device Services. Certificate biological safety Ceramill Zolid HT+; 2017.
- [72] Medical Device Services. Certificate biological safety Zpex 4 Yellow; 2017.
- [73] Medical Device Services. Certificate biological safety Zpex Gray; 2017.
- [74] Ferrari M, Cagidiaco MC, Garcia-Godoy F, Goracci C, Cairo F. Effect of different prosthetic abutments // Effect of different prosthetic abutments on peri-implant soft tissue. A randomized controlled clinical trial. *Am J Dent* 2015;28:85–9.
- [75] Stober T, Bermejo JL, Schwindling FS, Schmitter M. Clinical assessment of enamel wear caused by monolithic zirconia crowns. *J Oral Rehabil* 2016;43:621–9. <https://doi.org/10.1111/joor.12409>.
- [76] Al Hamad KQ, Abu Al-Addous AM, Al-Wahadni AM, Baba NZ, Goodacre BJ. Surface Roughness of Monolithic and Layered Zirconia Restorations at Different Stages of Finishing and Polishing: An In Vitro Study. *J Prosthodont* 2019;28:818–25. <https://doi.org/10.1111/jopr.13071>.
- [77] Saeed F, Muhammad N, Khan AS, Sharif F, Rahim A, Ahmad P et al. Prosthodontics dental materials: From conventional to unconventional. *Mater Sci Eng C Mater Biol Appl* 2020;106:110167. <https://doi.org/10.1016/j.msec.2019.110167>.
- [78] Farawati FA, Nakaparksin P. What is the Optimal Material for Implant Prosthesis? *Dent Clin North Am* 2019;63:515–30. <https://doi.org/10.1016/j.cden.2019.02.002>.
- [79] Ferrari M, Vichi A, Zarone F. Zirconia abutments and restorations: from laboratory to clinical investigations. *Dent Mater* 2015;31:e63-76. <https://doi.org/10.1016/j.dental.2014.11.015>.

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG

Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
germany@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

P002933878 27.05.2021 Version 2.0

AMANNGIRRBACH