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MDCG 2019	Summary of safety and clinical performance
-----------	--

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EN ISO 6872:2015 + A1:2018

ISO 6872:2015 + Amd.1:2018

DIN EN ISO 9693:2020	Dentistry - Compatibility testing for metal-ceramic and ceramic-ceramic systems
----------------------	---

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Reg. EU 2017/745	Regulation (EU) 2017/745 of the European Parliament and of the council on medical devices
MEDDEV 2.7/1 Revision 4	Clinical Evaluation - A Guidance for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC
Medical devices act	Austrian Federal Law on Medical Devices Version 2010 German Medical Devices Act Version 2010
DIN EN 1641:2010	Dentistry - Medical devices for dentistry - Materials
EN 1641:2009	
DIN EN ISO 13485:2016	Medical devices quality management systems
ISO 13485:2016	
DIN EN ISO 14971:2020	Medical devices - Application of risk management to medical devices
EN ISO 14971:2019	
ISO 14971:2019	
DIN EN ISO 15223-1:2017	Symbols for use in the labeling of medical devices
EN ISO 15223-1:2016	
ISO 15223-1:2016	
DIN EN 1041:2013	Information supplied by the manufacturer of medical devices
EN 1041:2008+A1:2013	
DIN EN 62366:2017	Medical devices - Part 1: Application of usability engineering to medical devices
EN 62366-1:2015 + AC:2015	
IEC 62366-1:2015 + COR1:2016	
DIN EN ISO 10993-1:2021	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-1:2020	
ISO 10993-1:2018	
DIN EN ISO 10993-5:2009	Biological evaluation of medical devices; Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-5:2009	
ISO 10993-5:2009	
DIN EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
EN ISO 10993-12:2012	
ISO 10993-12:2021	

Tab. 9

DIN EN ISO 10993-15:2009	Biological evaluation of medical devices - Part 15: Identification and quantification of degradation products from metals and alloys
EN ISO 10993-15:2009	
ISO 10993-15:2019	
DIN EN ISO 10993-18:2009	Biological evaluation of medical devices Part 18: Chemical characterization of materials
DIN EN ISO 10993-18:2021	Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process
EN ISO 10993-18:2020	
ISO 10993-18:2020	
DIN EN ISO 7405:2019	Dentistry - Evaluation of biocompatibility of medical devices used in dentistry
EN ISO 7405:2018	
ISO 7405:2018	
DIN EN ISO 10993-10:2014	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
EN ISO 10993-10:2013	
ISO 10993-10:2010	
DIN EN ISO 10993-11:2018	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
EN ISO 10993-11:2018	
ISO 10993-11:2017	
DIN EN ISO 14705:2021-04	
EN ISO 14705:2021	Fine ceramics (advanced ceramics, advanced technical ceramics) - Test method for hardness of monolithic ceramics at room temperature
ISO 14705:2016	
DIN EN ISO 843-1:2008	Advanced technical ceramics - Mechanical properties of monolithic ceramics at room temperature - Part 1: Determination of flexural strength
EN 843-1:2006	
ASTM E1876 - 15	Standard Test Method for Dynamic Young's Modulus, Shear Modulus, and Poisson's Ratio by Impulse Excitation of Vibration

Tab. 9

## References

- [1] Dogan S, Raigrodski AJ, Zhang H, Mancl LA. Prospective cohort clinical study assessing the 5-year survival and success of anterior maxillary zirconia-based crowns with customized zirconia copings. *J Prosthet Dent* 2017;117:226-32. <https://doi.org/10.1016/j.prosdent.2016.07.019>.
- [2] Worni A, Katsoulis J, Kolgeci L, Worni M, Mericske-Stern R. Monolithic zirconia reconstructions supported by teeth and implants: 1- to 3-year results of a case series. *Quintessence Int* 2017;48:459-67. <https://doi.org/10.3290/j.qi.a38138>.
- [3] Nejatidanesh F, Moradpoor H, Savabi O. Clinical outcomes of zirconia-based implant- and tooth-supported single crowns. *Clin Oral Investig* 2016;20:169-78. <https://doi.org/10.1007/s00784-015-1479-3>.
- [4] Rinke S, Schafer S, Lange K, Gersdorff N, Roediger M. Practice-based clinical evaluation of metal-ceramic and zirconia molar crowns: 3-year results. *J Oral Rehabil* 2013;40:228-37. <https://doi.org/10.1111/joor.12018>.
- [5] Bomicke W, Rammelsberg P, Stober T, Schmitter M. Short-Term Prospective Clinical Evaluation of Monolithic and Partially Veneered Zirconia Single Crowns. *J Esthet Restor Dent* 2017;29:22-30. <https://doi.org/10.1111/jerd.12270>.
- [6] Tartaglia GM, Sidoti E, Sforza C. Seven-year prospective clinical study on zirconia-based single crowns and fixed dental prostheses. *Clin Oral Investig* 2015;19:1137-45. <https://doi.org/10.1007/s00784-014-1330-2>.
- [7] Ortorp A, Kihl ML, Carlsson GE. A 5-year retrospective study of survival of zirconia single crowns fitted in a private clinical setting. *J Dent* 2012;40:527-30. <https://doi.org/10.1016/j.jdent.2012.02.011>.
- [8] Pott P-C, Eisenburger M, Stiesch M. Survival rate of modern all-ceramic FPDs during an observation period from 2011 to 2016. *J Adv Prosthodont* 2018;10:18-24. <https://doi.org/10.4047/jap.2018.10.1.18>.
- [9] Beuer F, Sachs C, Groesser J, Gueth J-F, Stimmelmayer M. Tooth-implant-supported posterior fixed dental prostheses with zirconia frameworks: 3-year clinical result. *Clin Oral Investig* 2016;20:1079-86. <https://doi.org/10.1007/s00784-015-1594-1>.
- [10] Guncu MB, Cakan U, Aktas G, Guncu GN, Canay S. Comparison of implant versus tooth-supported zirconia-based single crowns in a split-mouth design: a 4-year clinical follow-up study. *Clin Oral Investig* 2016;20:2467-73. <https://doi.org/10.1007/s00784-016-1763-x>.
- [11] Rinke S, Lange K, Roediger M, Gersdorff N. Risk factors for technical and biological complications with zirconia single crowns. *Clin Oral Investig* 2015;19:1999-2006. <https://doi.org/10.1007/s00784-015-1410-y>.
- [12] Koenig V, Vanheusden AJ, Le Goff SO, Mainjot AK. Clinical risk factors related to failures with zirconia-based restorations: an up to 9-year retrospective study. *J Dent* 2013;41:1164-74. <https://doi.org/10.1016/j.jdent.2013.10.009>.
- [13] Vigolo P, Mutinelli S. Evaluation of zirconium-oxide-based ceramic single-unit posterior fixed dental prostheses (FDPs) generated with two CAD/CAM systems compared to porcelain-fused-to-metal single-unit posterior FDPs: a 5-year clinical prospective study. *J Prosthodont* 2012;21:265-9. <https://doi.org/10.1111/j.1532-849X.2011.00825.x>.

- [14] Monaco C, Llukacej A, Baldissara P, Arena A, Scotti R. Zirconia-based versus metal-based single crowns veneered with overpressing ceramic for restoration of posterior endodontically treated teeth: 5-year results of a randomized controlled clinical study. *J Dent* 2017;65:56-63. <https://doi.org/10.1016/j.jdent.2017.07.004>.
- [15] Burke FJT, Crisp RJ, Cowan AJ, Lamb J, Thompson O, Tulloch N. Five-year clinical evaluation of zirconia-based bridges in patients in UK general dental practices. *J Dent* 2013;41:992-9. <https://doi.org/10.1016/j.jdent.2013.08.007>.
- [16] Ohlmann B, Eiffler C, Rammelsberg P. Clinical performance of all-ceramic cantilever fixed dental prostheses: results of a 2-year randomized pilot study. *Quintessence Int* 2012;43:643-8.
- [17] Pelaez J, Cogolludo PG, Serrano B, Serrano JFL, Suarez MJ. A four-year prospective clinical evaluation of zirconia and metal-ceramic posterior fixed dental prostheses. *Int J Prosthodont* 2012;25:451-8.
- [18] Passia N, Chaar MS, Kern M. Outcome of posterior fixed dental prostheses made from veneered zirconia over an observation period of up to 13 years. *J Dent* 2019;86:126-9. <https://doi.org/10.1016/j.jdent.2019.05.032>.
- [19] Chaar MS, Kern M. Five-year clinical outcome of posterior zirconia ceramic inlay-retained FDPs with a modified design. *J Dent* 2015;43:1411-5. <https://doi.org/10.1016/j.jdent.2015.11.001>.
- [20] Zenthofer A, Ohlmann B, Rammelsberg P, Bomicke W. Performance of zirconia ceramic cantilever fixed dental prostheses: 3-year results from a prospective, randomized, controlled pilot study. *J Prosthet Dent* 2015;114:34-9. <https://doi.org/10.1016/j.prosdent.2015.02.006>.
- [21] Sailer I, Balmer M, Husler J, Hammerle CHF, Kanel S, Thoma DS. 10-year randomized trial (RCT) of zirconia-ceramic and metal-ceramic fixed dental prostheses. *J Dent* 2018;76:32-9. <https://doi.org/10.1016/j.jdent.2018.05.015>.
- [22] Suarez MJ, Perez C, Pelaez J, Lopez-Suarez C, Gonzalo E. A Randomized Clinical Trial Comparing Zirconia and Metal-Ceramic Three-Unit Posterior Fixed Partial Dentures: A 5-Year Follow-Up. *J Prosthodont* 2018. <https://doi.org/10.1111/jopr.12952>.
- [23] Ioannidis A, Bindl A. Clinical prospective evaluation of zirconia-based three-unit posterior fixed dental prostheses: Up-to ten-year results. *J Dent* 2016;47:80-5. <https://doi.org/10.1016/j.jdent.2016.01.014>.
- [24] Larsson C, Vult von Steyern P. Ten-Year Follow-Up of Implant-Supported All-Ceramic Fixed Dental Prostheses: A Randomized, Prospective Clinical Trial. *Int J Prosthodont* 2016;29:31-4. <https://doi.org/10.11607/ijp.4328>.
- [25] Naenni N, Bindl A, Sax C, Hammerle C, Sailer I. A randomized controlled clinical trial of 3-unit posterior zirconia-ceramic fixed dental prostheses (FDP) with layered or pressed veneering ceramics: 3-year results. *J Dent* 2015;43:1365-70. <https://doi.org/10.1016/j.jdent.2015.07.013>.
- [26] Spies BC, Balmer M, Jung RE, Sailer I, Vach K, Kohal R-J. All-ceramic single crowns supported by zirconia implants: 5-year results of a prospective multicenter study. *Clin Oral Implants Res* 2019;30:466-75. <https://doi.org/10.1111/clr.13433>.
- [27] Rathmann F, Bomicke W, Rammelsberg P, Ohlmann B. Veneered zirconia inlay-retained fixed dental prostheses: 10-Year results from a prospective clinical study. *J Dent* 2017;64:68-72. <https://doi.org/10.1016/j.jdent.2017.06.008>.

- [28] Edelhoff D, Werder U. Komplexe Rehabilitation aus hochtransluzentem Zirkonoxid. Restaurationen von der Testphase bis zur definitiven Versorgung; 2020 Dec.
- [29] Haralur SB, Alamrey AA, Alshehri SA, Alzahrani DS, Alfarsi M. Effect of different preparation designs and all ceramic materials on fracture strength of molar endocrowns. *J Appl Biomater Funct Mater* 2020;18:2280800020947329. <https://doi.org/10.1177/2280800020947329>.
- [30] Juntavee N, Kornrum S. Effect of Marginal Designs on Fracture Strength of High Translucency Monolithic Zirconia Crowns. *Int J Dent* 2020;2020:8875609. <https://doi.org/10.1155/2020/8875609>.
- [31] Jansen JU, Lümkemann N, Sener B, Stawarczyk B. Comparison of fracture toughness measurements for zirconia materials using two test methods. *Dental materials journal* 2019;38:806-12. <https://doi.org/10.4012/dmj.2018-361>.
- [32] Almansour HM, Alqahtani F. The Effect of in vitro Aging and Fatigue on the Flexural Strength of Monolithic High-translucency Zirconia Restorations. *J Contemp Dent Pract* 2018;19:867-73.
- [33] Stawarczyk B, Ozcan M, Hallmann L, Ender A, Mehl A, Hämmerlet CHF. The effect of zirconia sintering temperature on flexural strength, grain size, and contrast ratio. *Clin Oral Investig* 2013;17:269-74. <https://doi.org/10.1007/s00784-012-0692-6>.
- [34] Amann Girrback R&D. Test report 3P bending strength Ceramill ZI; 2015.
- [35] Amann Girrback R&D. Test report 4P bending strength Ceramill ZI; 2015.
- [36] LMU Munich. Test report 4P bending strength aged Ceramill ZI; 2012.
- [37] Amann Girrback R&D. Test report 3P bending strength Ceramill Zolid HT+; 2017.
- [38] Amann Girrback R&D. Test report 4P bending strength Ceramill Zolid HT+; 2017.
- [39] Amann Girrback R&D. Test report 4P bending strength Ceramill Zolid HT+ aged; 2017.
- [40] Amann Girrback R&D. Test report 3P bending strength Zolid HT+ preshade A4; 2017.
- [41] Amann Girrback R&D. Test report 4P bending strength Zolid HT+ preshade A4; 2017.
- [42] Amann Girrback R&D. Test report 3P bending strength Zolid gen-x A3.5; 2019.
- [43] Amann Girrback R&D. Test report 4P bending strength Zolid gen-x A3.5; 2019.
- [44] Amann Girrback R&D. Test report bending strength (3-point) Zolid DRS B4.
- [45] Amann Girrback R&D. Test report bending strength (4-point) Zolid DRS B4.
- [46] Amann Girrback R&D. Test report 4P bending strength Zolid HT+ PS A4 aged; 2019.
- [47] Bioserv - Analytics and medical products GmbH. Test report animal studies toxicological and biological safety Ceramill ZI; 2010.
- [48] Bioserv - Analytics and medical products GmbH. Test report cytotoxicity Ceramill ZI; 2007.
- [49] Research Centre Jülich. Test report radioactivity Ceramill ZI; 2007.
- [50] Amann Girrback R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill ZI; 2017.
- [51] Medical Device Services. Test report cytotoxicity Ceramill Zolid HT+; 2017.
- [52] Medical Device Services. Test report chemical analysis (organic extractables) Ceramill Zolid HT+; 2017.
- [53] Medical Device Services. Test report chemical analysis (inorganic leachables) Ceramill Zolid HT+; 2017.

- [54] Research Centre Jülich. Test report radioactivity Ceramill Zolid HT+; 2016.
- [55] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Ceramill Zolid HT+; 2017.
- [56] Medical Device Services. Test report cytotoxicity Zpex 4 Yellow; 2017.
- [57] Medical Device Services. Test report cytotoxicity Zpex Gray; 2017.
- [58] Medical Device Services. Test report cytotoxicity Zpex Yellow and Zpex Pink; 2013.
- [59] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Gray; 2017.
- [60] Medical Device Services. Test report chemical analysis (organic extractables) Zpex4 Yellow; 2017.
- [61] Medical Device Services. Test report chemical analysis (organic extractables) Zpex Pink; 2013.
- [62] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex 4 Yellow; 2017.
- [63] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Gray; 2017.
- [64] Medical Device Services. Test report chemical analysis (inorganic leachables) Zpex Pink; 2013.
- [65] Research Centre Jülich. Test report radioactivity Zpex 4 Yellow; 2016.
- [66] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex Pink; 2017.
- [67] Amann Girrbach R&D. Test report chemical solubility DIN EN ISO 6872 Zpex 4 Yellow; 2017.
- [68] Medical Device Services. Biological safety Ceramill Zolid HT+; 2017.
- [69] Medical Device Services. Biological safety Zpex 4 Yellow; 2017.
- [70] Medical Device Services. Biological safety Zpex Gray; 2017.
- [71] Medical Device Services. Certificate biological safety Ceramill Zolid HT+; 2017.
- [72] Medical Device Services. Certificate biological safety Zpex 4 Yellow; 2017.
- [73] Medical Device Services. Certificate biological safety Zpex Gray; 2017.
- [74] Ferrari M, Cagidiaco MC, Garcia-Godoy F, Goracci C, Cairo F. Effect of different prosthetic abutments // Effect of different prosthetic abutments on peri-implant soft tissue. A randomized controlled clinical trial. *Am J Dent* 2015;28:85-9.
- [75] Stober T, Bermejo JL, Schwindling FS, Schmitter M. Clinical assessment of enamel wear caused by monolithic zirconia crowns. *J Oral Rehabil* 2016;43:621-9. <https://doi.org/10.1111/joor.12409>.
- [76] Al Hamad KQ, Abu Al-Addous AM, Al-Wahadni AM, Baba NZ, Goodacre BJ. Surface Roughness of Monolithic and Layered Zirconia Restorations at Different Stages of Finishing and Polishing: An In Vitro Study. *J Prosthodont* 2019;28:818-25. <https://doi.org/10.1111/jopr.13071>.
- [77] Saeed F, Muhammad N, Khan AS, Sharif F, Rahim A, Ahmad P et al. Prosthodontics dental materials: From conventional to unconventional. *Mater Sci Eng C Mater Biol Appl* 2020;106:110167. <https://doi.org/10.1016/j.msec.2019.110167>.
- [78] Farawati FA, Nakaparksin P. What is the Optimal Material for Implant Prosthesis? *Dent Clin North Am* 2019;63:515-30. <https://doi.org/10.1016/j.cden.2019.02.002>.
- [79] Ferrari M, Vichi A, Zarone F. Zirconia abutments and restorations: from laboratory to clinical investigations. *Dent Mater* 2015;31:e63-76. <https://doi.org/10.1016/j.dental.2014.11.015>.

Manufacturer | Hersteller  
Distribution | Vertrieb

**Amann Girrbach AG**

Herrschaftswiesen 1  
6842 Koblach | Austria  
Fon +43 5523 62333-105  
Fax +43 5523 62333-5119  
austria@amanngirrbach.com  
germany@amanngirrbach.com  
[www.amanngirrbach.com](http://www.amanngirrbach.com)

Made in the European Union

P002933878 27.05.2021 Version 2.0

**AMANNGIRRBACH**