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Original-Gebrauchsanweisung

PRODUKTBESCHREIBUNG

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INSTRUCTIONS FOR USE

Translation of the original instructions for use

PRODUCT DESCRIPTION

Milling wax for use in the CNC Ceramill Motion milling machine.

TECHNICAL DATA

Colour:	white to cream colours	
Dripping point:	100 - 130 °C	ASTM D 3954
Viscosity at 140 °C:	~ 650 mPa	DIN 53018
Solubility in water:	insoluble	
Acidity:	0 mg KOH/g	ISO 2114
Saponification number:	0 mg KOH/g	ISO 3681
Density at 23 °C:	0.92 - 0.96 g/cm ³	DIN 53479
Odour:	odourless	

Milling wax can be burned out without residue.

CHEMICAL COMPOSITION (PERCENTAGE BY WEIGHT)

100% Polyethylene wax

PACKAGE SIZES

760313	Ceramill WAX white 71 L (1 piece)	ht = 20 mm	Dental arch shape
760314	Ceramill WAX white 71 L (8 pieces)	ht = 20 mm	Dental arch shape
760315	Ceramill WAX white 71 XS (1 piece)	ht = 13 mm	Dental arch shape
760316	Ceramill WAX white 71 XS (8 pieces)	ht = 13 mm	Dental arch shape

STORAGE

Store in a dry place in its original packaging.
 (< 25 °C). Keep away from direct sunlight.

WARNING AND SAFETY INSTRUCTIONS

When processing, please use personal protective clothing (dust mask, protective goggles, ...).

INDICATIONS

_ Can be cast, pressed or over-pressed.
 _ No clinical dental product, not suitable for use in the oral cavity.

DESIGN OF THE SUBSTRUCTURE

Follow the manufacturer´s instructions for materials in which the substructure is cast or pressed.

GUARANTEE

Application recommendations, no matter whether verbally, in writing, or by means of practical instructions are based on our experience and tests and are thus only a recommendation. Our products are subject to continuous further development. We reserve the right to alter the products in terms of design or composition.

SAFETY DATA SHEETS AND DECLARATION OF CONFORMITY

Safety data sheets and declaration of conformity are available on request and under www.amanngirrbach.com

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