Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones de empleo

Manufacturer | Hersteller
 Distribution | Vertrieb

Amann Girschbach AG
 Herrschaftswiesen 1
 6842 Koblach | Austria
 Fon +43 5523 62333-105
 Fax +43 5523 62333-5119
 austria@amanngirschbach.com
 www.amanngirschbach.com

DIN EN ISO 13485
 ISO 9001

AMANNGIRSBACH

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Original-Gebrauchsanweisung

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INSTRUCTIONS FOR USE

Translation of the original instructions for use

Dental bake-on alloy on NiCrMo basis, type 3 according to DIN ISO 22674 for fixed dental prostheses.

For crowns and bridges.

Girobond CBS

Technical and physical data

(ISO 22674 and 9693)	
Proof stress (Rp 0.2)	400 MPa
Tensile strength (Rm)	650 MPa
Modulus of elasticity (E)	180
Elongation at rupture (A)	45 %
Vickers hardness VH 10	185
Melting range	1270 °C (Solidus) - 1356 °C (Liquidus)
Casting temperature	1410 °C
Density	8.4g/cm ³
WAK (25-500 °C)	13,8 x 10 ⁻⁶ K ⁻¹
(25-600 °C)	14,0 x 10 ⁻⁶ K ⁻¹
Colour	silver
Oxide colour	grey
Biocompatibility	corrosion resistance complies with ISO/FDIS 10271 cytotoxicity complies with ISO 10993-5

Packs		
Laboratory pack	1000 g	Order no. 781690

Chemical composition in % by mass

Nickel	Ni	63,5
Chromium	Cr	24
Molybdenum	Mo	10
Silicon	Si	1,5
Niobium	Nb	< 1
Manganese	Mn	< 1
Does not contain beryllium or carbon.		

Shelf life

Up to 5 years from the date of manufacture.

FABRICATING A CROWN OR BRIDGE TO THE PORCELAIN FACING STAGE

- Waxing up**
 - Dip the coping to a wall thickness of 0.3 - 0.5 mm.
 - A maximum uniform porcelain thickness of 2 mm above the wax pattern is ideal.
- Attaching sprues**
 - Use a 5 mm thick horizontal or ring sprue.
 - Attach a 3 mm Ø x 2 mm long feeder sprue to the highest point. Shape the junctions to ensure a smooth flow of the metal.
- Investing**
 - Line the ring. Use two layers with larger rings (sizes 6 and 9).
 - Place the pattern in the centre of the ring and the horizontal sprue in the heat centre at a uniform distance from the edge of the ring.
 - Invest according to the manufacturer's instructions. Recommended investment Giroinvest Speed and Giroinvest Super.
- Preheating**
 - Burn out the wax.
 - Set the heat rate according to the instructions for the investment.
 - Select the holding time at casting temperature according to the size of the ring. Casting temperature according to manufacturer's instructions (950 °C - 1050 °C).
- Casting (casting temperature approx. 1420 °C, melting range 1270-1356 °C)**
 - Place the Girobond CBS ingots (1 ingot - 6 g) into the ceramic crucible.
 - 1.) Casting with the flame**
 - Gas-oxygen ratio: Propane-oxygen ~ 0.7:1.4 bar
Acetylene-oxygen ~ 0.7:0.7 bar
 - Cast immediately after the last ingot has melted into the molten metal and the surface of the molten metal splits open and becomes shiny.
Note: do not remove the flame from the alloy during melting.
 - 2.) Induction casting machine**
Cast immediately at the moment when the surface of the molten metal splits open and becomes shiny.
 - 3.) Vacuum /pressure casting machines**
Do not melt the alloy in vacuum. Set to the chrome cobalt programme, resp. switch off the vacuum or reduce it as far as possible. Cast immediately as soon as the last ingot has melted into the molten metal and the surface of the molten metal splits open, then trigger off the casting machine.
 - Allow the ring to air cool slowly after casting.
Note: the melting procedure can be shortened if the crucibles are preheated in the furnace. Do not use the same crucibles for other alloys. Do not overheat the alloy. Cast only new metal so that each batch can be clearly identified.
- Devesting and preparing**
 - 6.1 Sandblast with 110 - 250 µm aluminium oxide.
 - 6.2 Cut off the sprues with a cutting disc.
 - 6.3 Prepare the framework with a tungsten carbide cutter (cross-cut).
- Metal treatment before facing with porcelain**
 - 7.1 Fine sandblast the framework with 110 - 250 µm pure aluminium oxide with a pencil blaster at 3-4 bar.
 - 7.2 Grip the framework with haemostatic forceps. Do not handle again.
 - 7.3 Steam clean the framework with a steam cleaner (Steamer X3 116910).
 - 7.4 Oxidise, if required, at 1040 °C with a 1 - 2 minute holding stage to ensure a clean surface. Then sandblast again (refer to 7.1 - 7.3). **Note:** do NOT dry with compressed air due to possible contamination with oil or dirt. Do not immerse the framework in hydrofluoric acid.
- Porcelain facing**
 Follow the porcelain manufacturer's instructions. When using Creation porcelain (788000), carry out the first opaque firing at 1000 °C with a 1 min. holding stage. After building up the porcelain, allow it to dry for 6 min. before firing. **Note:** a longer cooling phase may be required with larger span bridges.
- Finishing**
 Polish the metal edges to a high lustre with silicone polishers avoiding any heat build-up.
- Laser welding**
 - 10.1 Sandblast the areas to be welded with 110 - 250 µm aluminium oxide and then clean them.
 - 10.2 Prepare the connector areas in the form of a modified cross-stitch.
 - 10.3 After welding the Girobond CB framework, it can be faced with porcelain.
- Soldering**
 Presolder using NiCrMo solder (781630 Giro solder) and the appropriate flux. Post solder in the furnace at 860 °C with white gold solder and the appropriate flux.
- Storage**
 Store in a dry place at room temperature.
- Side effects**
 Allergies to the constituents of the alloy or discomfort caused by electrochemical reaction are possible in rare cases. Girobond CBS should not be used in patients with known allergy to nickel.
- Guarantee**
 Recommendations for technical use, whether given orally, in writing or in the course of practical instructions, are based on our own findings and tests and are therefore only to be used as guidelines. Our products are subject to continuous development. We therefore reserve the right to make any changes to the structure and composition which improve the product.
- Contraindications**
 Known intolerance to the constituents.
- Warning/safety information**
 Damage to health caused by dust produced when trimming alloys.
- Declaration of conformity and safety data sheet** are available on request at www.amanngirschbach.com

Dieses Produkt enthält Nickel!
 This product contains nickel!
 Ce produit contient du nickel !
 Questo prodotto contiene nichel!
 Este producto contiene níquel!

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