



<b>DE</b> Bedienungsanleitung	4-5
<b>EN</b> User manual	6-7
<b>FR</b> Mode d'emploi	8-9
<b>IT</b> Istruzioni d'uso	10-11
<b>ES</b> Modo de empleo	12-13

**AMANNGIRRBACH**

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**OPEN CENTRIC**

Guide pivoting lever **(2)** up to stop, moving the centering axes (2 semi-axes) in, thus releasing the centric.

Centric opened **(3)**. Fixing of condyle balls released. The top may be freely moved.

**CLOSE CENTRIC**

Move pivoting **(4)** lever down to stop. The two semi-axes are guided out until softly stopping against the condylar centre, fixing the Artex top section in centered position. Centric closed. Visible fixation on condylar centre by semi-axis. 2 Support in inclined

Centric closed **(5)**. Visible fixation on condylar centre by semi-axis.

**ADJUSTMENT OF BENNETT ANGLE**

Open centric **(6)**. Release locking screw of Bennett angle with allen key. Adjust condylar unit from -5° to +30° by rotary button.

Fig.: Bennett angle = 0°.

Rotate Bennett angle **(7)** to required scale setting by means of the rotary knob, then lock by allen key.

Fig.: Bennett angle = 30°.

**ADJUSTMENT OF CONDYLE PATH INCLINATION**

Release condyle **(8)** path inclination (HCI) set screw with allen key. Adjust required angle by rotating the condylar box. then fix HCI set screw with allen key.

Adjustment **(9)** of HCI from -20° to +60°.

Read angle of inclination from front edge of rear section (congruent with the scale).

Fig.: HCI = +60°.

**SEQUENCE OF MOVEMENT**

Condyle balls in centered position and during movement with different Bennett angles (view into condyle)

**ATTENTION**

Sequences of movement in Arcon articulators must be performed in permanent contact of condyle balls (bottom part of articulator) with the condylar box (top part of articulator). Only thus they can be reproduced.

**(10)** Bennett angle adjusted to 0°. Condylar spheres in centered position = initial position for movements.

**(11)** Bennett angle adjusted to 0°. Condyle ball guided in all three dimensions in permanent contact with the condylar unit. The condyle ball may never lose contact with the fossa, neither in cranial nor mesial direction.

**(12)** Bennett angle adjusted to 30°. Condyle ball in centered position = initial position for movements.

**(13)** Bennett angle adjusted to 30°. Condyle ball guided in all three dimensions in permanent contact with the condylar box.

**ARBITRARY PIN / CENTRE OF AXIS**

Maxillary cranium/axis-based transfer of the maxilla can be effected directly on the articulator by adapting the facebow. Both the arbitrary pin **(14)** (which reproduces the position of the ear opening/bore of the porus buttons) and the axis mean point can be used as a reference.

**DISTRACTION**

Loose screw from frontal. **(15)** Install distraction to the desired value (0-3 mm) and fix it with the screw **(16)**.

**RETRUSION**

**(17)** Open centric, release locking screw of retrusion insert using allen key. Guide retrusion insert **(18)** to required individual value, then relock using allen key.

Fig.: Retrusion = 2 mm.

**! WARNINGS AND SAFETY ASPECTS:****\_ Anterior guidance**

Angle of protrusion track =  $0^\circ$

**\_ General information**

Articulators are precision instruments. Therefore, it is important to handle them with care. Please read the following notes:

- \_ Do not place the articulator in a boil-out unit, a pressure pot or an ultrasonic bath.
- \_ For cleaning and care, use only cleaning agents that are suitable for light metals. Regarding the use of the cleaning agent, please follow the instructions of the according manufacturer. Do not use cleaning agents with a pH  $<5$  or  $>8$  since those might damage the anodized surfaces.

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**Amann Girrbach AG**