Registrieren | Artikulieren
Registration | Articulation
Enregistrement de l'occlusion
Montage sur articulateur
Registrazione | Messa in
articolatore
Registro | Articulación

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User Manual

Mode d'emploi


Istruzioni d'uso

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Deutsch	04 - 04
English	05 - 05
Française	06 - 06
Italiano	07 - 07
Español	08 - 08

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ARTEX® CT

OPEN CENTRIC

Swivel locking (1) bar upwards with the thumb; the top part of the Artex can be guided backwards freely. The locking bar limits travel in a posterior direction and prevents the top part of the Artex falling out.

Centric open (2), condyle ball released.

CLOSE CENTRIC

Press locking (3) bar down; the plastic tongue centres the condyle ball in the lowest position of the condyle path. Only rotation is possible with the top part.

Centric closed (4): the condyle ball is fixed in position by the plastic tongue.

SETTING BENNETT ANGLE

Open the Bennett setting screw (5), turn condyle to desired position and fix position with the setting screw. Bennett angle can optionally be set from 0° to 20°.

Fig.: Bennett angle = 0°.

SETTING CONDYLE PATH

Open the condyle setting screw (6). Set condyle to desired inclination and fix position with setting screw. HCI can be set from -15° to +60°.

Fig.: HCI = +15°.

ARBITRARY PIN

If the facebow is to be adapted "arbitrarily" directly on the articulator, i.e. via the bore on the ear balls, the condylar box has to be set at 60° and the Bennett angle at 0°. In this position, the arbitrary pin reproduces the (average-value) position of the ear opening at which the facebow was held up to the patient. (Adaptation at the axis mean point is always preferable).

! WARNINGS AND SAFETY ASPECTS:

- Anatomical anterior guidance

Angle of protrusion track = 40°;
on lateral movements the angle varies between 30 and 40°, according to the predetermined Bennett angle

- How to check the centric adjustment of Artex CN and Artex CT

If the upper part of the articulator has lateral play, this deviation to the side can be heard and felt by a "click". In this case, it is necessary to readjust the centric.

Then, please proceed as follows:

- Set the Bennett angle to 10°.
- Close the centric lock.
- Unscrew the thread pins at the bottom of the rear part using an allen wrench.
- In the anatomically shaped anterior guidance plate, the anterior guidance pin is centrally adjusted at zero.
- Hold the upper articulator part in this position. At the same time, pull the left semi-axis outwards until reaching the stop, and use the allen wrench to tighten the thread pin. Now proceed in the same way with the right semi-axis.

- General information

Articulators are precision instruments. Therefore, it is important to handle them with care.

Please read the following notes:

- Periodically check if the axis-joint system and the condylar paths are clean in order to maintain the smooth mobility of the condyle in the axis and thus maintain the functionality of the articulator.
- Remove dirt in this area with a steamcleaner, then blow dry.
- Do not place the articulator in a boil-out unit, a pressure pot or an ultrasonic bath.
- For cleaning and care, use only cleaning agents that are suitable for light metals. Regarding the use of the cleaning agent, please follow the instructions of the according manufacturer. Do not use cleaning agents with a pH <5 or >8 since those might damage the anodized surfaces.

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Amann Girrbach AG
Amann Girrbach GmbH