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
DE Aufbereitungsanweisung	2 - 9
EN Instructions for processing	10 - 17
FR Instructions de traitement	18 - 25
IT Istruzioni per la preparazione	26 - 33
ES Instrucciones de procesamiento	34 - 41

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
PROCESSING (CLEANING, DISINFECTION, AND STERILIZATION) OF COMPONENTS (SEMI-CRITICAL CONTACT)

Fundamental points

All components with semi-critical contact to the patient are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile components (cleaning and disinfection after removal of the protective packaging, sterilization after packaging). An effective cleaning and disinfection is an indispensable requirement for an effective sterilization of the components.

You are responsible for the sterility of the components with semi-critical contact to the patient. Therefore, please ensure that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilization, that the used devices (WD, sterilizer) will be maintained and checked regularly, as well as that the validated parameters will be applied for each cycle.

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA).

For all further components (without contact to the patient) standard surface disinfection may be applied (not content of this reprocessing instruction).

CLEANING AND DISINFECTION

Basics

If possible, an automated procedure (WD (Washer-Disinfector)) should be used for cleaning and disinfection of the components. A manual procedure - even in case of application of an ultrasonic bath - should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered.

The pre-treatment step is to be performed in both cases.

Pre-treatment

Please remove coarse impurities of the components directly after application (within a maximum of 2 h).

Procedure:

1. Disassemble the components as possible (see chapter Specific aspects).
2. Rinse the components at least 1 min under running water¹ (temperature < 35 °C/95 °F). If applicable (see chapter Specific aspects): Rinse all lumens of the components five times by application of a single-use syringe (minimum volume 5 ml).
3. Remove manually all visible impurities by use of a clean and soft brush (or a clean, soft, and lint-free cloth) only to be for this, in no case metal brushes or steel wool. If applicable (see chapter Specific aspects)²: Brush holes by use of a conical interdental brush.
4. Rinse again at least 1 min under running water. If applicable (see chapter Specific aspects): Rinse all lumens of the components five times by application of a single-use syringe (minimum volume 5 ml).

¹ In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety) please consider, that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the components (see chapter „material resistance“).

Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety, but cannot replace the disinfection step later to be performed after cleaning.

² Still remaining remnants of registration material easily can be removed after cooling-down in a refrigerator.

Automated cleaning/disinfection (WD (Washer-Disinfector))

Pay attention to following points during selection of the WD:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- possibility for an approved program for thermal disinfection (A0 value > 3000 or - in case of older devices - at least 5 min at 90 °C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the components)
- fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional application - in case of non-application of a thermal disinfection - of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible to the used cleaning detergent
- compatibility of the used detergents with the components (see chapter „material resistance,,)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

Procedure:

1. Transfer the disassembled components in the WD (pay attention that the components have no contact). If applicable (see chapter Specific aspects): Use a small pieces basket for this.
2. Start the program.
3. Remove the components of the WD after end of the program.
4. Check and pack the components immediately after the removal (see chapters „check,, „maintenance,, and “packaging,, if necessary after additional post-drying at a clean place).

The fundamental suitability of the components for an effective automated cleaning and disinfection was demonstrated by an independent accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the cleaning detergent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

Manual cleaning and disinfection

Pay attention to following points during selection of the cleaning and disinfection detergents:

- fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- in case of application of an ultrasonic bath: suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the components (see chapter „material resistance,,)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

- Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

Procedure:

Cleaning

1. Soak the disassembled components for the given soaking time in the cleaning solution so that the components are sufficiently covered. Pay attention that there is no contact between the components and no remaining air in the cavities. Assist cleaning by careful brushing with a soft brush and ultrasonic treatment (duration at least once the required soaking time in addition to brushing).

If applicable (see chapter Specific aspects):

Rinse all lumens of the components at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5 ml). Brush holes by use of a conical interdental brush.

2. Then, remove the components of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.

If applicable (see chapter Specific aspects): Rinse all lumens of the components at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5 ml).

3. Check the components (see chapters „check,, and „maintenance,,).

Disinfection

4. Soak the disassembled components for the given soaking time in the disinfectant solution so that the components are sufficiently covered. Pay attention that there is no contact between the components and no remaining air in the cavities.

If applicable (see chapter Specific aspects):

Rinse all lumens of the components at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5 ml).

- Then, remove the components of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water.

If applicable (see chapter Specific aspects): Rinse all lumens of the components at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5 ml).

- Dry and pack the components immediately after the removal (see chapter „packaging,,“ if necessary after additional post-drying at a clean place).

The fundamental suitability of the components for an effective cleaning and disinfection was demonstrated by an independent accredited and recognized (§ 15 (5) MPG) test laboratory by application of the cleaning detergent Cidezyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt) considering the specified procedure.

CHECK

Check all components after cleaning or cleaning/disinfection, respectively, on corrosion, damaged surfaces, and impurities. Do not further use damaged components (for limitation of the numbers of re-use cycles see chapter „reusability,,“). Still dirty components are to be cleaned and disinfected again.

MAINTENANCE

Instrument oils or grease must not be used.

Do not mount prior to packaging and sterilization (no sterilization of the other components without contact to the patient admitted)

PACKAGING

Please insert the cleaned and disinfected components single-use sterilization packagings (single packaging), which fulfill the following requirements (material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 142 °C (288 °F), sufficient steam permeability)
- sufficient protection of the components as well as of the sterilization packagings to mechanical damage

STERILIZATION

Please use for sterilization only the listed sterilization procedures; other sterilization procedures must not be applied.

Do not sterilize the other components without contact to the patient.

Steam sterilization

- fractionated vacuum/dynamic air removal procedure³ or gravity displacement procedure⁴ (with sufficient product drying⁵)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min ⁵	at least 30 min at 121 °C (250 °F), drying time at least 20 min ⁵
other countries	at least 3 min ⁶ at 132 °C (270 °F)/134 °C (273 °F)	at least 20 min at 121 °C (250 °F)

STORAGE

Please store the components after sterilization in the sterilization packagings at a dry and dust-free place.

MATERIAL RESISTANCE

Please take care that the listed substances are not ingredients of the cleaning or disinfection detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 6.5)
- lyes (maximum admitted pH-value 8.5, neutral/enzymatic cleaner recommended)
- organic solvents (for example: acetone, ether, alcohol, benzene)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

Please do not clean any components by use of metal brushes or steel wool.

Please do not applied acidic neutralization agents or any rinse aids.

Please do not expose any components to temperatures higher than 142 °C (288 °F)!

REUSABILITY

The components can be reused - in case of adequate care and if they are undamaged and clean - 50times. The user is responsible for each further use as well as for the use of damaged and dirty components (no liability in case of disregard).

The fundamental suitability of the components for an effective steam sterilization was demonstrated by an independent accredited and recognized (§ 15 (5) MPG) test laboratory by application of the steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and the fractionated vacuum/dynamic air removal procedure as well the gravity displacement procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered.

The flash/immediate use sterilization procedure must not be used.

Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

3 at least three vacuum steps

4 The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure.

5 The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

6 respectively 18 min (inactivation of prions, not relevant for USA)

SPECIFIC ASPECTS


Art. no	article specification	rinsing volume	brush	specific/additional procedure in case of			packaging		sterilization	maximum admitted cycle number	recommended classification according to RKI/BAfM/RINKE guideline (only Germany, with respect to intended use)
				pretreatment	manual cleaning/ disinfection	automated cleaning/ disinfection	Maintenance				
218607 218609	Porus buttons	5 ml (single-use syringe)	interdental brush conical (3-6.5 mm), tooth brush	after removing the support bracket: Brush off outside and inside	five-times rinsing inside, brushing outside and inside	small parts basket	lubrication not admitted	standard procedure	50	semi-critical B	
217695	Nasion adapter	5 ml (single-use syringe with needle)	interdental brush conical (3-6.5 mm), tooth brush	after removing the nasion bar: Brush off outside and inside	Rinse five times; brush off outside and inside	small parts basket	lubrication not admitted	standard procedure	50	semi-critical B	
217928	Tray handle	n.a.	tooth brush	after removing the universal joint: Brush off outside and inside	Rinse five times; brush off outside and inside	small parts basket	lubrication not admitted	standard procedure	50	semi-critical B	
217611	Bite fork partial	n.a.	tooth brush	after removing the universal joint: Brush off outside and inside	Rinse five times; brush off outside and inside	standard basket	lubrication not admitted	standard procedure	50	semi-critical B	
217661	Bite fork full	n.a.	tooth brush	after removing the universal joint: Brush off outside and inside	Rinse five times; brush off outside and inside	standard basket	lubrication not admitted	standard procedure	50	semi-critical B	
218606	Set screws	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	not admitted	n.a.	without contact to the patient	
218610	Axis plane indicator optional	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	not admitted	n.a.	without contact to the patient	
218635	Support bracket right/left (pair)	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	not admitted	n.a.	without contact to the patient	
218608	Removable axis pins (pair)	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	not admitted	n.a.	without contact to the patient	

218616	Support bracket screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218692	Fixing screw for nasion	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218693	Set screw for nasion	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218613	Nasion bar	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218695	Nasion support	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218632	Adapter to bow	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
217696	Adapter screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
217633	Fork mounting	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
217634	Collar screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218696	Toggle nut, complete	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
218697	Toggle screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)	not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient

For drawings and further information about articles see IFU of Artrex facebook on www.amannjrbsch.com/instruction-manuals.

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Amann Girrbach AG