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Registrieren | Artikulieren
Registration | Articulation
Enregistrement de l'occlusion
Montage sur articulateur
Registrazione | Messa in
articolatore
Registro | Articulación

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
User Manual

Mode d'emploi

Istruzioni d'uso

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AMANNGIRRBACH

_Deutsch	04
_English	05
_Française	06
_Italiano	07
_Español	08

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ARTEX® TRANSFER STAND

PRODUCT DESCRIPTION

For safe transport of the facebow registration to the lab. The universal joint with bite fork is removed from the facebow and inserted into the reception of the transfer stand. There, the registration is fixed in plaster. The detachable transfer table safely and unshakeably takes this important patient information to the lab.

OPERATION/HANDLING

The transfer stand before mounting the 3-D universal joint. The transfer table has been moved down into its lowest position. Hexagon socket screw is open. Allen key is inserted in the rubber bushing so that it is always available. Guide the 3-D universal joint into the socket of the transfer stand and fix in position using adapter screw **(Fig. 1)**. Support toggle pin **(Fig. 2)** with the little finger. The gap between the transfer platform and the bite fork can be seen.

Move transfer table platform to bite fork. Retention screws are at a distance from the bite fork. Fix table platform in position by tightening up hexagon socket screw. **(Fig. 3)**

Remove 3-D universal joint from transfer stand. To do so open adapter screw completely. Tip bite fork upwards. Release toggle pin with little finger. **(Fig. 4)**

Mix articulation plaster (Artifix 711217) to a creamy consistency and apply lengthwise over the retention screws on the transfer table platform in two portions. **(Fig. 5)**

Guide 3-D universal joint at an angle towards the socket of the transfer stand. Release toggle pin with the little finger. **(Fig. 6)**

Fix 3-D universal joint in socket of transfer stand at the stop and sink the bite fork into the soft pattern resin by tightening up the adapter screw. **(Fig. 7)**

When the plaster has set, open the collar screw at the bite fork mounting. Open the adapter screw and remove the 3-D universal joint from the transfer stand. The facebow registration is securely fixed on the transfer platform and ready for transport to the laboratory. **(Fig. 8)**

Note

The illustrations here show the transfer in articulators fitted with the Splitex® calibration system (magnetic profile plates). The procedure for screwed transfer tables is the same.

Individual parts:

Order No.:	
218670C	Transfer Stand with table for carbon version
216240	Splitex Transfer Stand complete

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Amann Girrbach AG
Amann Girrbach GmbH