<b>DE</b> Aufbereitungsanweisung	3 - 6
<b>EN</b> Instructions for processing	7 - 10
<b>FR</b> Instructions de préparation	11 - 14
<b>IT</b> Istruzioni per la preparazione	15 - 18
<b>ES</b> Instrucciones de preparación	19 - 22

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- Translation of the original Instructions for processing -

## 1 Preparing the Components

All parts in semi-critical contact with the patient (see chapter 7) must be cleaned, disinfected and sterilized before each use. This is also necessary for the first use after delivery of the non-sterile components (cleaning and disinfection after removal from the protective packaging, sterilisation after packaging). Effective cleaning and disinfection are indispensable requirements for effective sterilisation of components.

For all other components (without patient contact) common surface disinfection can be applied (not included in these instructions).

## 2 Cleaning and Disinfection

### 2.1 Pretreatment

Coarse contamination of the components must be removed immediately after use (within a maximum of two hours).

If a cleaning and disinfection agent is used instead of water for pretreatment (e.g. to protect personnel), it must have the following properties:

- \_ General suitability for cleaning instruments made of metal or plastic
- \_ Compatibility with the components (see section 5)
- \_ Tested efficiency (e.g. approval/authorisation/registration with/from VAH/DGHM or FDA/EPA or CE marking)
- \_ Aldehyde-free (otherwise fixation of impurities with blood)

A disinfectant used during pre-treatment does not replace disinfection after cleaning.

For pretreatment of components:

- ▷ If possible: Disassemble components (see section 7).
- ▷ Rinse the (disassembled) components under running water for at least one minute (temperature < 35° °C).
- ▷ Remove all visible contamination with a clean and soft brush.  
The brush may only be used for this purpose. Do not clean with a metal brush or steel wool.
- ▷ Again, rinse the components under running water for at least one minute.

## 2.2 Selection of the Cleaning and Disinfection Agents

The cleaning and disinfection agents must have the following properties:

- \_ Cleaning agents
  - General suitability for cleaning instruments made of metal or plastic
  - Compatibility with the components (see section 5)
  - When using an ultrasonic bath: Suitability of the cleaning agent for ultrasonic cleaning (no foaming)
- \_ Disinfection agents
  - Tested efficiency (e.g. approval/authorisation/registration with/from VAH/DGHM or FDA/EPA or CE marking)
  - Compatible to the cleaning agent being used

Combined cleaning and disinfection agents may only be used in the case of exceptionally minor contamination (no visible contamination).

## 2.3 Cleaning/Disinfection

Cleaning and disinfection differ only by the agents used, not by the procedure.

### NOTE:

Insufficient results due to wrong agents or procedures!

- ▷ Observe the instructions of the manufacturers of the cleaning and disinfection agents regarding concentration, temperature, immersion time and rinsing.
- ▷ Only use freshly prepared solutions.
- ▷ Only use sterile or slightly contaminated water (e.g. purified/highly-purified water):
  - max. 10 germs/ml
  - max. 0.25 endotoxin units/ml
- ▷ For drying, use a soft, clean and lint-free cloth and/or filtered air.
- ▷ If possible: Disassemble components.
- ▷ Immerse the (disassembled) components completely in the cleaning/disinfection solution for the duration of the specified immersion time.
  - The components may not touch each other.
  - No air may remain in cavities or hollow spaces.
- ▷ Thoroughly clean the components with a soft brush and/or apply ultrasonic treatment for the duration of the specified immersion time. The brush may only be used for this purpose. Do not clean with a metal brush or steel wool.
- ▷ Remove the components from the cleaning/disinfection solution and rinse intensively with water at least three times (for at least one minute).

### 3 Checking

- ▷ After cleaning/disinfection, check all components for corrosion, damaged surfaces and contamination.
- ▷ Do not continue to use damaged components.
- ▷ Clean and disinfect components again that are still contaminated.
- ▷ When the maximum number of preparation cycles has been reached: Dispose of components.

### 4 Maintenance

#### NOTE:

Contamination of the cleaned/disinfected components!

- ▷ Do not apply oils or greases to the components.

- ▷ Do **not** assemble components before packing and sterilisation.

Components not in contact with patients may not be sterilised!

### 5 Material Resistance

The cleaning/disinfection agents may **not** contain the following substances:

- \_ Organic, mineral and oxidizing acids (lowest permissible pH value 6.5)
- \_ Alkaline solutions (maximum permissible pH value 8.5); neutral/enzymatic cleaner recommended
- \_ Organic solvents (e.g. acetone, ether, alcohol, petrol)
- \_ Oxidizing agents (e.g. peroxides)
- \_ Halogens (chlorine, iodine, bromine)
- \_ Aromatic, halogenated hydrocarbons

#### NOTE:

Possible damage to components due to improper cleaning!

- ▷ Do **not** clean components with metal brushes or steel wool.
- ▷ Do **not** use acid neutralisers or rinsing agents.
- ▷ Do **not** subject components to temperatures above 142 °C.

### 6 Reusability

Undamaged and clean components can be reused 50 times.

The user is responsible for any use beyond this (no liability in the event of non-observance).

## 7 Special Factors

### Specific/additional procedure for

Article No.	Product Designation	Rinsing Specification	Brush	Manual Cleaning	Disinfection	Maintenance	Packing, Sterilisation	Max. Reusability	Recommended classification <sup>1)</sup>
242700	Centrofix								
242710	Tracing tablet, wide	At least 1 minute under running water	Clean, soft brush	Rinse five times on inside; brush off on outside and inside	After removing the support bracket: Brush off outside and inside	Lubrication not permitted	Standard procedure	50	Semi-critical B
242720	Tracing tablet, narrow								
242740	Tracer mounting								
242750	Tracing ball, long								
242751	Tracing ball, short								
242760	Alignment clamp								
242790	Universal key								
242780	Alignment cross	Not specified	Not specified	Not specified	Not specified	Not specified	Not permitted	Not specified	No contact with patient
242810	Copper tracing tablet								
242820	Spirit level								
242830	Space holder (plastic hose)								

Tab. 1

1) according to Directive RKI/BfArM/KRINKO (Germany only, in view of intended use)

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**Amann Girrbach AG**