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
User manual

Mode d'emploi



Istruzioni d'uso

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Deutsch	04 - 05
English	06 - 07
Française	08 - 09
Italiano	10 - 11
Español	12 - 13
EU-Declaration of conformity	14 - 15

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USER MANUAL

- | | |
|---|--------------------------|
| 1 | Main switch |
| 2 | Magnetic disk |
| 3 | Vertical adjusting screw |
| 4 | Fork wrench SW 17 |
| 5 | Fork wrench SW 15 |
| 6 | Allen key |

DESCRIPTION OF APPARATUS

The Laserpin is an ergonomically designed pin drilling unit featuring a combination of precision mechanics and new laser technology. Its 30° angle of design is the key to optimum use, as it enables the operator to remain comfortably seated while working - yet have a clear view of the model. The inclined position of the drilling table prevents plaster dust from affecting the drilling area, as it simply falls below during work on the model. All mechanical parts are designed for maximum precision and long service life. The drills are available in five different diameters and can be easily changed by using the two wrenches supplied. This makes Laserpin to an absolute necessity for the dental technician.

INSTALLATION

After unpacking the unit, please check it immediately for possible damage incurred during the transport. Should the unit have been damaged during the transport, please complain to your supplier immediately. Although this unit is supplied ready for use, always check the voltage setting (115/230V) immediately above the power socket. If you are not certain which voltage is used in your country, please ask your local supplier of electrical appliances.

MOUNTING AND CHANGING THE DRILL

Always switch the unit off first. Use the open-end wrenches (4+5) supplied with the unit to open the chuck. Hold the bottom wrench firmly in place and turn the top wrench clockwise. This will open the chuck so that the drill can be inserted or replaced through the top opening (B). (Measure the diameter of the pins used with your system) Insert the desired drill into the chuck and push it approx. 1-2 mm beneath the white magnetic disk.

Tighten the chuck again by holding the bottom wrench firmly and turning the top wrench anti-clockwise.

ADJUSTING THE DRILLING DEPTH

Loosen the vertical adjusting screw (3) with the Allen key (6). Only tighten the vertical adjusting screw very slightly (it should remain slightly mobile, (D)). Press the platform down slowly to establish the desired drilling depth (depending on the pin system).

Drill depth = retention part of pin + 1mm

Use the Allen key to tighten the vertical adjusting screw again. (A)

The unit is now ready for use. Switch the power switch on - the laser beam is aimed at the hole.

HANDLING

The user can either stand or sit to operate the Laserpin. When operating the unit in sitting, ensure that the model can be seen from above and that the unit can be operated with both hands. Place the model or dental arch on the white magnetic disk.

The laser beam aims exactly at the spot on the bottom of the model where the hole is to be drilled. When the model or white magnetic disk is depressed, the motor switches on and the hole is drilled. (E)

IMPORTANT: To prevent the model or drill slipping off-center or even fracturing, always press the model firmly onto the magnetic disk with both thumbs! (C)

AMBIENT CONDITIONS/OPERATING CONDITIONS

- Ambient temperature 0 C to +40 C
- max. atmospheric humidity 90%

SERVICING

Like most AmannGirrbach equipment, Laserpin is designed to be maintenance-free. Simply keep the unit clean. Always remove plaster dust completely. To remove the plaster dust, simply lift the magnetic disk. The used E-fuse is a type 0,8A slow. When changing the fuse make sure to use one of the same type. If you use a different type E-fuse the safety of the Laserpin cannot be granted any more.

! WARNINGS AND SAFETY ASPECTS:

- never push your hand on the base plate
- switch off the apparatus when not using it
- while cleaning the unit, switch of the mains
- never grasp a drill while it is rotating
- use only the original drill
- when drilling, hold model securely and press on the base plate
- do not look directly into the laser beam
- when drilling performance declines, change the drill do not place finger on drill hole and press down (special care should be taken during start-up)
- when pushing down the magnetic disc the drill appears
- Precaution with long hair and sleeves

Technical Data:

D/W/H:	360 x 140 x 240 mm
Weight:	6 kg
Mains supply:	230V 50Hz
E-fuse:	T 0,8 A
Power:	160 W
Motor speed:	3000 min ⁻¹
Accuracy:	< 0,03 mm
Drill stroke:	continuously adjustable 0-12,5 mm
Collet capacity:	ø 2,5-3,00 mm
Laser class:	2
Laser wavelength:	680 nm

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Amann Girrbach AG
Amann Girrbach GmbH