# Sliding guide

<b>DE</b> Gebrauchsanweisung	3 - 7
<b>EN</b> Instruction Manual	8 - 12
<b>FR</b> Instructions d'utilisation	13 - 17
<b>IT</b> Istruzioni per l'uso	18 - 22
<b>ES</b> Instrucciones de uso	23 - 26

- Original Gebrauchsanweisung -

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## Table of Contents

<b>1</b>	<b>Explanation of Symbols</b> .....	<b>9</b>
<b>2</b>	<b>General Safety Instructions</b> .....	<b>9</b>
<b>3</b>	<b>Suitable Personnel</b> .....	<b>10</b>
<b>4</b>	<b>Product Specifications</b> .....	<b>10</b>
4.1	Intended use .....	10
4.2	Product Description .....	10
4.3	Delivery Scope .....	10
4.4	Safety data sheet/Declaration of conformity .....	10
<b>5</b>	<b>Preparation before using on the patient</b>	<b>10</b>
<b>6</b>	<b>Application</b> .....	<b>11</b>
6.1	Selection of the Sliding Guide after checking the sagittal level .....	11
6.2	Deprogramming the engrams (occlusal contacts) .....	11
6.3	Registration .....	11
6.4	Checking .....	12
<b>7</b>	<b>Warranty/Exclusion of Liability</b> .....	<b>12</b>
<b>8</b>	<b>Environmental Protection</b> .....	<b>12</b>

# 1 Explanation of Symbols

## Warning indications

Warning indications in the text are marked with a triangle and boxed.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

## Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

## Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

## Other symbols on the product

Symbol	Meaning
	Medical device
	Unique Device Identifier
	Article number
	Batch code
	Manufacturer
	Date manufacturing (YYYYMM)
	Do not use again

# 2 General Safety Instructions

When installing, starting-up and operating the product, always observe the following safety instructions:

### CAUTION:

- A defective product can lead to malfunctions!  
 On detection of damage of the product:
- ▷ Label the product as defective.
  - ▷ Prevent further use.
  - ▷ Submit complaint, if required.

### NOTE:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### 3 Suitable Personnel

**NOTE:**

The product may only be used by trained specialised personnel.

### 4 Product Specifications

#### 4.1 Intended use

Sliding Guide is used to determine and encode the jaw relation in dentate patients.

#### 4.2 Product Description

Sliding Guides are so-called jigs made of plastic by injection moulding. The jig is tapered from the handle towards the tip. The jigs are available in different thicknesses so that, depending on the sagittal level of the patient, it is possible to block the rows of upper and lower teeth for determination of the jaw relation in dentate patients.

The jig is placed between the anterior front teeth for minimal blocking of the teeth rows. The patient holds this jig between the teeth for a certain time and is deprogrammed (with regard to the occlusal contacts when assigning upper to lower jaw). This leads to unmanipulated guidance of the temporomandibular joints into their rotation center.

The resulting positional relationship is fixed with bite registration material and later used to mount the lower jaw model against the upper jaw model in the articulator.

The shelf life of the Sliding Guide is 10 years.

### 4.3 Delivery Scope

**Art. No.**

504611	Woelfel-Sliding Guide 4 mm, 10 pce.
504621	Woelfel-Sliding Guide 9 mm, 10 pce.

### 4.4 Safety data sheet/Declaration of conformity

The product's declaration of conformity is available upon request from the manufacturer.

### 5 Preparation before using on the patient

Sliding Guide is delivered non-sterile. Before using:

- ▷ Clean and disinfect Sliding Guide with a disinfectant suitable for polycarbonate (PC).
- ▷ To remove disinfectant residue: rinse with running water or clean with steamer.

Sliding Guide is a disposable article. Repeated use and the required sterilisation is not intended for.

## 6 Application

### 6.1 Selection of the Sliding Guide after checking the sagittal level

The Sliding Guide cannot be used for mandibular prognathism and end-to-end bite.

To select the appropriate Sliding Guide:

- ▷ In case of habitual intercuspitation, determine the size of the sagittal level by eye.  
For teeth that have not been prepared, the 4 mm Sliding Guide is usually sufficient.
- ▷ Slide in the Sliding Guide in the final-bite position in the sagittal level until there is resistance.
- ▷ Read the size of the sagittal level on the Sliding Guide scale of the incisal edge.

### 6.2 Deprogramming the engrams (occlusal contacts)

The patient should sit relaxed in the chair. The head is slightly overstretched to the rear. This position stabilizes the lower jaw.

- ▷ To lock, slide the Sliding Guide between the front teeth (sagittal level + 1 - 1.5 mm).
- ▷ Check posterior region to see if all molars disclude (separate).

The Sliding Guide remains in this position in the mouth for at least 5 minutes.

In case of premature contacts during the deprogramming:

- ▷ To undo the contact: Increase the blocking by repositioning the Sliding Guide.  
From the time of the last tooth contact until the final deprogramming, the Sliding Guide must remain in the patient's mouth for at least 5 minutes.

### 6.3 Registration


#### NOTE:

Loss of the deprogramming!

If the patient closes the mouth, the deprogramming is lost and must be repeated.

- ▷ Do not let patients close their mouth.

- ▷ Apply bite registration silicone directly to the patient's lower jaw teeth.
- ▷ In the first part of the closing movement, guide the patient by the lower jaw until the anterior teeth touch the Sliding Guide.
- ▷ Have the patient to hold the bite registration evenly and relaxed with the lower jaw until the registration material has set.

- ▷ When the registration material has set: Remove the bite registration.
- ▷ Note down the mm value of the anterior tooth lock that was selected for registration on the order form.

## 6.4 Checking

- ▷ Inspect registration for correctness and explicit cusp impressions for model orientation.

It is important to understand that the bite registration aligns to the opposite cusp, but is held without contact by the anterior guidance (no direct contact of antagonistic cusps). This guides the patient from the centric relation into the centric occlusion.

## 7 Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

## 8 Environmental Protection

### Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally friendly and recyclable.

### Disposal of the products

Disposal must take place in accordance with local/regional/national regulations.

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**Amann Girrbach AG**