 **af350**

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OPERATING INSTRUCTIONS

AF350

Translation of the original operating instructions

- ❶ Check light for magnet
- ❷ On/Off switch for magnet
- ❸ Motor revolution regulator
- ❹ Switch for clockwise/ anticlockwise rotation of motor
- ❺ Motor On/Off switch
- ❻ Base
- ❼ Locking screw for jointed arm
- ❽ Clamping screw for motor unit + adapter sleeve
- ❾ Locking screw for vertical adjuster
- ❿ Revolution indicator
- ⓫ Spring for balance
- ⓬ Micrometer
- ⓭ Motor cable
- ⓮ Main switch On/Off
- ⓯ Pivoted lever for drilling
- ⓰ Fuse holder with voltage selector
- ⓱ Magnetic plate
- ⓲ Height fixation
- ⓳ Modelling fixation screw
- ⓴ Centring ring
- ⓵ Fixation screw for inclination
- ⓶ Modelling table MT3
- ⓷ Fixation screw for height adjustment
- ⓸ Adapter cable AF350 adjustable
- ⓹ Spindle column
- ⓺ Clamp
- ⓻ LED light
- ⓼ Socket for electrical connection
- ⓽ Left/right running indicator
- ⓾ Socket for foot-operated switch
- ⓿ Motor unit

PRODUCT DESCRIPTION

The AF350 is a universal milling device developed in cooperation with specialists in milling technology. The result was an apparatus which was distinguished for its easy and ergonomic handling. Since only proven materials were used in finishing, the apparatus is technically at the highest level.

UNPACK AND CHECK

Immediately after unpacking, check the appliance for any transport damage. If any damage has occurred during transport, please submit a complaint to your supplier immediately.

INTENDED USE

The AF350 is suitable for the following areas of application: milling and drilling of various technical dental works.

For example:

- _ Circular milling
- _ Supports for clasps
- _ Interlocking milling
- _ Individually-made RSS milling attachments
- _ Fabrication of tapered and telescopic crowns
- _ Parallel bars
- _ Attachments

Operating conditions:

- _ temperature 0°C to 40°C
- _ max. atmospheric humidity 90%
- _ Noise level with max. speed and working distance of appr. 20 cm - 73dB (A)

INITIAL START-UP

Electric connection

Site the appliance and plug the power cable (supplied) into the socket on the back of the appliance ❸ (check correct voltage setting at the fuse holder ❶ below the power cable connection). Now plug the power plug into the mains socket.

Changing the voltage

The appliance can be switched from 230 Volt to 100 Volt or vice versa. To that end, pull out the fuse holder 16, remove fuses and position the insert (voltage selector) in such a way that the desired voltage is visible in the small window. Afterwards, reinsert the fuse box complete with fuses.

CONTROL ELEMENTS IN THE BASE

The base 6 houses various actuators and switches to control the milling device. At the back of the base is the main On/Off switch 14, the connection for the power supply 28, the fuse holder with the voltage selector 16 and the socket for the plug of the foot-operated switch. On the front are the switch for the motor 5, the switch to activate the magnet to the activation of the magnet to lock the modelling table 22 and the control to regulate the motor's direction 4 and speed 3. The current revs are digitally displayed in the display window 10. Two LED lights are mounted on the base, lighting up automatically when the mains switch is switched on.

INSTALLING THE JOINTED ARM

Assembling the jointed arm

The jointed arm is fitted onto the vertical adjuster 9; the first step is to fit the height adjustment onto the vertical adjuster with the aid of the centring wedge (pointing upwards) and to tighten it with an Allen key. Then attach the jointed arm and tighten in place with the locking screw 7. For the assembly, it is important for the jointed arm to be attached parallel and untilted. The next step is to insert the motor unit 31 into the designated holder and fixed at the side with the Allen screw 8. The power lead is plugged into the top side of the spindle column, and the spiral spring secured.

Height adjustment

The height can be adjusted in several places, primarily however at the spindle column 25. The first step is to open the lever; the desired height is then set with the aid of the large plastic screw. The height can also be changed at the vertical adjuster 9. In addition, this provides the option of moving the jointed arm away to the side and returning it to the definite position with the aid of the height fixation ring. Drilling depth can be roughly set with the graduation at the milling head; it is set or changed by tightening or loosening the locking screw. The jointed arm can be completely locked in place with the locking screw for the jointed arm fixation 7. The spring tension adjuster 11 is used to create different counterpressure for free milling and drilling. To that end, move the plastic part to the side and adjust the black ribbed part to the top or bottom. Height for drilling is set with the micrometer screw 12. This can be adjusted over a total range of 25 mm. One revolution equates to 1 mm. The graduation is given in tenths of a millimetre.

The milling head is changed by opening the spring chuck 25. The choice of milling head depends on the nature of the milling and the construction. It is important that the motor is switched off. A 2.35 mm diameter spring chuck is fitted as standard. To change the spring chuck, open the spring and unscrew it manually using a chuck key. Screw the new spring chuck in to the stop. Insert the milling head or mandrel into the chuck and close. Then set the main switch at the back 14 to ON and check all functions of the AF350.

The drilling lever 15 is only plugged on and can be pulled off to the side when not in use. By switching between clockwise/ anticlockwise rotation of the motor 4, the milling head can be used for smoothing the object by reversing the motor direction of the original cut. Any blockage during the switch of direction causes the appliance to switch off automatically. The motor restarts when using the On/Off switch;

make sure that the reason for the blockage has been ascertained before switching the motor back on.

ASSEMBLY OF THE MODELLING TABLE

The patented MT3 modelling table has been specifically developed for milling. The model is fixed with the finger screw ⑱. Pulling the centring ring ⑳ downwards allows the MT3 to move freely. Activating the magnetic plate ㉑ locks the MT3 in its inclination and position. Switching the magnetic plate ㉒ off releases the MT3 to move freely again. The MT3 can be locked in any position mechanically using the screw on the side of the appliance ㉓, even without magnetic locking. To unlock the table, open the screw slightly and press it in; the MT3 can move freely again. To return the MT3 into a definite position (0 degrees), push the centring ring upwards.

ACCESSORIES (p. 27)

Marking stick and adjustable adapter cable

The marking stick and the adjustable cable are available as accessories for this appliance. To use the marking stick, first remove the motor unit ㉔: with an Allen key, loosen the screw fixing the motor to the milling arm, the motor removed and the power cable at the spindle column unplugged. The adjustable adapter cable consists of two parts:

1. the adapter with temperature control;
2. the adapter cable. The marking stick is inserted into the adapter, both together are inserted into the milling arm and fastened with the Allen screw. The adapter cable is plugged into the spindle column and the spiral spring secured.

Blockout set B2

The set includes a heating element with a connecting cable, five inserts and a key. Inserts are shaped in 0°, 2°, 4° and 6°, with the 0° insert available in two diameters (1.3 mm and 2.0 mm).

To use the blockout set, the heating element is attached to the lower part of the marking stick. The inserts are attached to the heating element by twisting, and the connecting cable of the heating element is plugged into the adjustable adapter cable. To release the insert, insert the key between the heating element and the insert and tilt downwards.

ATTENTION: Check all plugs for proper connections.

ATTENTION

Pull the plug!

Basically, the apparatus is maintenance-free. Dirty marks should be cleaned dry. No steaming and do not clean with water or solvents.

WARNING AND SAFETY INSTRUCTIONS:

A.S.M. products have been developed for use in the dental industry. Start-up and operation must be done according to the Instructions for Use. The manufacturer cannot accept responsibility for damages occurred in other uses or incorrect handling. The user is furthermore obliged on his own responsibility to test the apparatus, prior to use, for its suitability and utilization in objectives intended if these objectives are not stated in the Instructions for Use.

_ The fixation screw for height adjustment must always be in touch with the swivel arm to prevent the arm from falling down.

_ When milling metal, plastic or wax wear protection goggles.

_ Persons with long hair should wear it tied up or with a hairnet.

_ By using maximum heat with the block-out there is a danger of burning the skin.

_ While cleaning or servicing the apparatus, switch off the mains.

_ Do not steam clean and do not clean with water or solvents.

_ Repairs must be done only by authorized persons.

Overload protection:

The motor is switched off automatically when overloaded (e.g. seized by milling tool, open clamp etc.). Red lamp shows overload.

Reset:

_ Rectify the reason

_ Start motor again with on-/off-switch

Technical data:

L/B/H:	340 x 250 x 450 mm
Weight:	15 kg
Electrical connection:	100/115/230V, 50/60 Hz
Foot pedal connection:	With plug connection on/off switch
Block-out connection:	24 V
E-fuse:	T1A
Motor:	1.500 to 40.000 rpm
Performance:	125 W
Revolutions:	7,5 Ncm
Clamping:	2,35 mm

Packing list:

Order No.:	
177605	AF350 incl. vertical adjustment 1 Milling arm, incl. motor unit 1 Clamp Ø 2,35 mm 1 Modelling table MT3 Halogen lamp 1 Dust cover

Accessories/ Individual parts:

Order No.:	
177998	Block out Upgrade Kit AF350 adjustable* 3 parts (for the latest version see image on title page)
177995	Upgrade kit AF350 4 parts
177994	Adjustable adapter cable
177800	Marker
177661	W&H® Chuck set
873001	Milling Starter Kit

*Blockout Set 177990, Adjustable adapter cable AF350
177994, Marking stick 177800

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Amann Girrbach AG
Amann Girrbach GmbH