CAD/CAM Accessories

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3 - 10

EN Instruction Manual

11 - 18

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D-Set Production (181216) / BD Production (181217)

Calibration of blank holder

D-Set Production (181216) / BD Production (181217)

- Original Gebrauchsanweisung -

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
Bild 7 Autokalibrierung starten

- Translation of the original Operating Instructions -

Table of Contents

1	Explanation of Symbols	11
2	Calibration of Blank Holder D-Set	
	Production	12
2.1	Product description	12
2.2	Delivery scope	12
2.3	Required system components	13
2.4	Calibrating the blank holder	13
2.4.1	Preparation of the production machine	13
2.4.2	Calibration of the blank holder	14
3	Calibration of Blank Holder BD	
	Production	15
3.1	Product description	15
3.2	Delivery scope	15
3.3	Required system components	16
3.4	Calibrating the blank holder	16
3.4.1	Preparation of the production machine	16
3.4.2	Calibration of the blank holder	17
4	Auto-calibration	18
5	Article Numbers	18

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
_	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 Calibration of Blank Holder D-Set Production

Before using the blank holder D-Set Production (181216) in the production machine, the blank holder must be calibrated.

NOTE:

The blank holder is individually calibrated to the respective production machine. Afterwards, the holder can no longer be used in other production machines.

2.1 Product description

Calibratable blank holder D-Set Production for use in production machines. The blank holder is used for fixation of prefabricated D-set blanks.

The blank holder consists of two components:

- _ Blank holder
- _ Fixation plate

2.2 Delivery scope


Fig. 1 Blank holder D-Set Production

- 1 Blank holder (uncalibrated)
- 2 Fixation plate
- 3 Allen screws M4x10 (7 pce.)

The blank holder is delivered completely assembled.

2.3 Required system components

- _ CAD/CAM software for removable prosthetics (D-Flow)
- _ Production machine Ceramill Matik
- _ Roto RFID 1.2 ST (751001)

In order to ensure the software functionality:

- ▷ Before using, install the current update of all software components.

2.4 Calibrating the blank holder

The calibration process consists of two steps that must be performed one after the other.

1. Preparation of the production machine
2. Calibration of the blank holder

 **NOTE:**

Insufficient milling results when used in other production machines!

The blank holder is individually calibrated to the respective production machine.

- ▷ Afterwards, do not use holder in other production machines.

2.4.1 Preparation of the production machine

To ensure correct calibration of the blank holder:

- ▷ Perform auto-calibration (please refer to section 4 on page 18).
- ▷ Perform collet cleaning (see maintenance booklet of production machine).

2.4.2 Calibration of the blank holder

NOTE:

Danger of tool breakage!

Calibration with the fixation plate screwed on will lead to cutter breakage and incorrect calibration.

- ▷ Loosen screws and remove fixation plate.


Fig. 2 Blank holder D-Set Production

- ▷ Insert blank holder (without fixation plate and screws) into the material changer.
- ▷ Carry out a blank holder inventory run.
- ▷ Insert tool Roto RFID 1.2 ST (751001) into the tool magazine and close the magazine. The production machine carries out an inventory run of the tool magazine.
- ▷ Start the production machine software on the PC.
- ▷ Switch to calibration folder by clicking on *Calibration* button.

- ▷ Select *Calibration matik D-Set* file and then click on *Transfer* button to transfer data to production machine.


Fig. 3 Transfer "Calibration matik D-Set" file

- ▷ If *Start jobs after transfer* was not selected: Start job manually at the production machine. 16 calibration marks are milled for fixation of the tooth blanks in the correct position. The calibration process takes approx. 35 minutes.
- ▷ When the calibration is completed: Mount fixation plate back to blank holder with the seven screws.

The blank holder is now calibrated for this production machine and can immediately be used for machining D-Set blanks.

3 Calibration of Blank Holder BD Production

Before using the blank holder BD Production (181217) in the production machine, the blank holder must be calibrated.

NOTE:

The blank holder is individually calibrated to the respective production machine. Afterwards, the holder can only be used in another production machine after replacing the sliding blocks and recalibration.

3.1 Product description

Calibratable blank holder BD Production for use in production machines. The blank holder is used to affix prefabricated blanks of the Baltic Denture System.

The blank holder consists of four components:

- _ Blank holder
- _ Sliding blocks 10 mm
- _ Sliding block 8 mm
- _ Fastening ring

3.2 Delivery scope


Fig. 4 Blank holder BD Production

- 1** Sliding block 10 mm
- 2** Allen screws M3x6 (4 pce.)
- 3** Blank holder
- 4** Sliding block 8 mm
- 5** Fastening ring
- 6** Allen screws M3x8 (4 pce.)

The blank holder is delivered completely assembled.

3.3 Required system components

- _ CAD/CAM software for Baltic Denture
- _ Production machine Ceramill Matik
- _ Roto RFID 1.2 ST (751001)

In order to ensure the software functionality:

- ▷ Before using, Install the current update of all software components.

3.4 Calibrating the blank holder

The calibration process consists of two steps that must be performed one after the other.

1. Preparation of the production machine
2. Calibration of the blank holder

 **NOTE:**

Insufficient milling results when used in other production machines!

The blank holder is individually calibrated to the respective production machine.

- ▷ Afterwards, the holder can only be used in another production machine after replacing the sliding blocks and recalibration.

3.4.1 Preparation of the production machine

To ensure correct calibration of the blank holder:

- ▷ Perform auto-calibration (please refer to section 4 on page 18).
- ▷ Perform collet cleaning (see maintenance booklet of production machine).

3.4.2 Calibration of the blank holder

NOTE:

 Danger of tool breakage!
 Calibration with the fixation ring screwed on will lead to cutter breakage and incorrect calibration.

▷ Loosen screws and remove fixation ring.


Fig. 5 Blank holder BD Production

- ▷ Insert blank holder (without fixation ring and screws) into the material changer.
- ▷ Carry out a blank holder inventory run.
- ▷ Insert tool Roto RFID 1.2 ST (751001) into the tool magazine and close the magazine.
 The production machine carries out an inventory run of the tool magazine.
- ▷ Start the production machine software on the PC.
- ▷ Switch to calibration folder by clicking on *Calibration* button.

- ▷ Select *Calibration matik BDS* file and then click on *Transfer* button to transfer data to production machine.


Fig. 6 Transfer "Calibration matik BDS" file

- ▷ If *Start jobs after transfer* was not selected:
 Start job manually at the production machine.
 The sliding blocks are milled for fixation of the Baltic-Denture blanks in the correct position.
 The calibration process takes approx. 20 minutes.
- ▷ When the calibration is completed: Mount fixation ring back to blank holder with the four screws.

The blank holder is now calibrated for this production machine and can immediately be used for machining Baltic-Denture blanks.

4 Auto-calibration

NOTE:

Incorrect calibration!
 A soiled production machine leads to insufficient results.

- ▷ Carry out auto-calibration only on a dust- and grease-free, completely dried production machine.

- ▷ Insert auto-calibration blank holder (181204) into the material changer.
- ▷ Carry out a blank holder inventory run.
- ▷ Insert Roto RFID calibration pin (181205) into tool magazine.
- ▷ Carry out inventory run of tool magazine.
- ▷ Start auto-calibration in the service menu by clicking on *Calibrate machine* button.


Fig. 7 Starting the auto-calibration.

5 Article Numbers

Art. No.	Designation
181,216	Blank holder D-Set Production
181217	Blank holder BD Production
751001	Roto RFID 1.2 ST
181204	Auto-calibration blank holder
181205	Roto RFID calibration pin

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG

Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

ISO 9001

36228-FB_2020-08-26