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
Important information for using the new Argovent components

General

- After initial use, the grey sintering tray has whitish streaks on the surface from burning out the binder. This is normal and does not indicate any impairment.
- The sintering beads adhere to the sintering cover during the first sintering cycles. This behaviour decreases over time.
- Beads adhering to the edge of the sintering tray can be wiped off when removing the sintering cover. Any excess adhering beads can be wiped off manually, e.g. using the reverse side of the shielding gas retort.
- Replace the loosened beads into the sintering tray.
- A maximum of 100 g of sintering beads may be filled in the sintering tray. This supplied amount should never be exceeded as the sintering tray may fracture as a result of pressure from the beads.
- Generally, 70 g of sintering beads is recommended to prevent any possible damage to the sintering tray.
- It is essential to fit the sintering cover as shown in the information sheet supplied.
- The sintering tray combination has a 2.5 mm smaller diameter in the usable space (sintering cover 2S).
- The new components do not require any maintenance.
- The instructions for use of the sintering furnace should always continue to be observed!

Use of the float sintering disc

- The float sintering disc represents the usable surface.
- The beads must be completely removed beforehand if the float sintering disc has to be used for sintering!
- It is sufficient when the float sintering disc is supported on the separating disc and the sintering cover sits on the floor of the tray.
- Otherwise there is the risk of getting wedged too tightly by the beads between the sintering cover and float sintering disc!
- Use of a separate sintering cover is recommended for frequent float sintering, as the lower section of the sintering cover can move slightly inwards due to the permanent pressure of the beads. This results in the float sintering disc not fitting perfectly in the sintering cover. This effect can be prevented by using a separate sintering cover.
- Where possible, it is recommended to work alternately between float sintering and bead sintering to avoid oxidation on the sintered restorations. Residual powder on the separating disc/sintering tray after bead sintering promotes oxide-free sintering.
- Partial discoloration is non-critical, but should not exceed the permitted amount. (See Ceramill Argotherm sintering chart)

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EN

The sinter cover must not simply be pressed into the beads in the tray. It must be grasped on its cross handle and rotated clockwise/anti-clockwise into the beads several times until it touches the tray bottom. This fully loosens the beads all around to prevent them jamming and exerting excessive pressure on the sinter tray. Otherwise there is an acute risk that the parts may break due to thermal expansion!