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**GB IMPORTANT NOTE ON UPGRADING THE CERAMILL BASE TO CERAMILL MULTI-X**

- A Because of production tolerances in the manufacture of Ceramill Base the scale of the milling guide may be positioned differently. In some devices, during upgrade to Multi-x this may cause the scaling to interfere with the socket of the milling guide (see photo) so this is not completely seated solidly. In this case the scaling needs repositioning.
- B Screw off the milling guide and mark the original position of the scaling on the left and on the right with a water-soluble pen (IMPORTANT!!!). Unscrew the two screws of the socket and take the socket off the bearing.
- C Unfasten the two screws at the underside of the socket (see arrows) and shift the scaling so that the milling guide can be positioned upon the socket without leaving a chink - do not tighten the screws.
- D Place the socket, together with the milling guide, at the edge of a table so that the fastening screws of the scaling are accessible. Push the scaling towards the guide - as a spacer, place a strip of paper in between - and position the scaling horizontally according to the markings applied before. Tighten the screws again. Screw the socket onto the bearing. Attach the milling guide to the socket and control the placement of the scaling by shifting the milling guide to the right and to the left.

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