Betriebsanleitung
Operating instructions
Mode de fonctionnement
Istruzioni sul
funzionamento
Instrucciones de servicio

Deutsch	02 - 04
English	05 - 07
Français	08 - 10
Italiano	11 - 13
Español	14 - 16

EG-Declaration of conformity	18
---------------------------------	----

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PRODUKTBESCHREIBUNG

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OPERATING INSTRUCTIONS CERAMILL UV

Translation of the original operating instructions

PRODUCT DESCRIPTION

Ceramill UV is an LED-lamp for the polymerisation of light-hardening dental materials. Its wavelength and controls are designed so that they ensure fast and optimal polymerisation especially for the modelling material Ceramill Gel and Pontic.

UNPACKING AND CHECKING

After unpacking the machine, check immediately for any possible damage that occurred in transit. If any damage did occur in transit, then please inform your supplier at once. The device is delivered in a fully operational state.

Please keep the transport packaging!
The packaging the device came in protects the device from damage in transit.

SETTING UP THE DEVICE/ASSEMBLY

Assemble the orange visor as shown in **Fig. 1**.

COMMISSIONING

Check that the voltage specified on the label on the back of the casing and on the mains adaptor corresponds to the mains voltage available.

Do not attempt to carry out any alterations on the machine or on the accompanying mains adaptor by yourself! The protective earthing of the supply

point must be fully functioning.

CAUTION! The manufacture accepts no liability for any damage arising from improper installation.

OPERATION/HANDLING

Plug in mains adaptor and connect to the power supply - the polymerisation LED will briefly light up to check functionality. Switch on the Ceramill UV using the ON/OFF button ①. The control LED will shine green and the machine is ready for operation.

Ceramill UV is fitted with a movement sensor. If an object is held about 2 cm under the lamp, it will switch on automatically and contact-free. After 7 seconds, an acoustic signal will sound: this signals that the optimum polymerisation time of Ceramill Gel and Pontic modelling material has been reached. If the object remains under the lamp, the next 7-second polymerisation cycle will be automatically triggered.

When the polymerisation cycle is triggered, please ensure that neither the LED nor the sensor is fouled with the synthetic material.

If the polymerisation cycle is not triggered for a period of 5 minutes, the lamp automatically switches itself off. If an object remains under the lamp for long-term polymerisation (at a distance of approx. 2 cm), the lamp will also automatically switch itself off after 5 minutes. After pressing the ON/OFF button again, the machine is ready for operation once more.

Adjustment of sensitivity

The sensitivity of the sensor depends on the lighting conditions in the workplace. It might be necessary, therefore, to adjust the sensitivity. Please carry out any adjustments carefully using an Allen key (2 mm), as shown in **Fig. 3**.

_ More sensitive: turn in a counter-clockwise direction

_ Less sensitive: turn in a clockwise direction

The sensor is optimally set when the polymerisation-LED switches on at a distance of approx. 20-30 mm. The final position (last setting in a clockwise direction) should not be selected since this setting is not adequately sensitive.

MAINTENANCE/CARE

Pull out the plug!

The machines are basically maintenance-free.

When cleaning, keep the device dry.

CAUTION! Do not allow evaporation and do not clean with water or solvents.

Technical data:

D/B/H:	145 x 300 x 30 mm (gooseneck 800 mm)
Weight:	2,6 kg
Main supply:	110 - 230 V 50/60 Hz
Initial value mains adaptor:	7,5 V, 1,2 A
Polymerisation-LED:	600 mW, 450 - 470 nm

Packing list:

Order No.:	
178200	Ceramill UV
	inc. mains adaptor
	visor
	Polymerisation-LED

Accessories/Components:

Order No.:	
178201	Mains adaptor 110 - 230 V
178210	Ersatz LED
178211	Blendschutz

WARNING AND SAFETY INSTRUCTIONS:

Ceramill UV was developed for use in the dental industry. Commissioning and operation must be performed precisely according to the instructions for use. The manufacture accepts no liability for damage caused by inappropriate use or improper handling. Furthermore, the operator is obliged to check the machine before using it. On his own responsibility, he should check that the machine is suitable and fit for use for the designated purposes, insofar as these purposes are not specified in the instructions for use.

All repairs and/or servicing jobs, apart from those described under Maintenance/care, may only be carried out by the manufacturer or authorised service providers.

- _ During the polymerisation process, do not look at the light without a visor!**
- _ Do not foul the sensors and the LED with synthetic material**

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Amann Girrbach AG
Amann Girrbach GmbH