ceramil map200

DE Betriebsanleitung	3- 18
EN User Manual	19- 33
FR Mode d'emploi	34- 48
IT Istruzioni d'uso	49- 63
ES Modo de empleo	64- 78

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Tab. 3

- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	20	9	Environmental Protection	32
2	General Safety Instructions	21	10	Technical data	33
3	Suitable Personnel	21	11	Accessories	33
4	Machine Specifications	22			
	4.1 Delivery scope	22			
	4.2 Intended Use	22			
	4.3 CE Declaration of Conformity	22			
	4.4 Components and interfaces	23			
5	Installation	24			
	5.1 Transport	24			
	5.2 Unpacking	24			
	5.3 Setting up	25			
	5.4 Installation of the software on the PC	25			
	5.5 Installation of the scanner	25			
	5.6 Activating additional modules	25			
	5.7 Calibration of axes	26			
6	Application Tips	27			
	6.1 Universal holder	27			
	6.2 Positioning and affixing single models	28			
	6.3 Positioning and affixing model pairs	29			
7	Cleaning and Maintenance	31			
	7.1 Cleaning	31			
	7.2 Maintenance	31			
8	Malfunctions, Repairs and Warranty ..	32			
	8.1 Malfunctions	32			
	8.2 Repairs	32			
	8.3 Warranty	32			

EN

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Fuse
	USB camera
	Mains switch

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

 **CAUTION:**
A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

 **NOTE:**

- ▷ In case of ± 15 degree temperature deviations, allow the device to acclimate at least four hours at room temperature.
- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.
- ▷ When lifting and transporting the scanner, observe the transport and storage instructions.
- ▷ Please do not place objects on top of the scanner, as there is a risk that they may fall off due to the surface properties and the resulting vibrations during the scanning operation.

The scanner comprises a motion unit of two electrically-driven rotary axes, which are used to position the object during the scan process. To enable the clamping of an object in a defined position and a more swift workflow during the scan process, the motion axes are held electrically in position. An electronic protection mechanism prevents unintentional starting of the motors as well

as continued running of motors when reaching in. Since the theoretical possibility of a failure of the protection mechanism is given, there is a residual risk that you can counteract against by observing the following guidelines:

 **CAUTION:**

Never reach below the area marked by the warning labels, between rocker and housing! This also applies when the rocker is pivoted back, in order to remove fallen in objects from the scanner or to clean the scanner interior.

- ▷ It is imperative to switch the main switch of the device off.

When the device moves in one or more motion axes or otherwise reacts unusual when reaching in:

- ▷ Do not reach into the scanner interior.
- ▷ Switch the device off immediately.
- ▷ Do not continue to use the device; mark the device as defective until the fault has been corrected and ensure the device is not restarted.

3 Suitable Personnel

 **NOTE:**
Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Delivery scope

- _ Desktop scanner
 - USB cable (2 x)
 - Mains cable
 - Operating Instructions
 - CD with installation software and calibration data
 - Calibration model
 - Universal model holder
 - Modelling wax
- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.2 Intended Use

The desktop scanner is a PC-controlled 3D scanner for scanning dental tooth models.

NOTE:

The scan result depends on the materials being scanned. Metal surfaces, for example, can lead to artefacts. Scan spray can improve the scan capability of materials.

For improper use of the system, the manufacturer shall not assume any liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

4.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety. In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking. The product's Declaration of Conformity can be requested from the manufacturer.

EC Guidelines	Applied harmonised standards
2006/42/EC (Machinery Directive)	DIN EN ISO 12100
2014/35/EC (Low Voltage Directive)	DIN EN 61326-1
2014/30/EC (EMV directive)	DIN EN 61010-1

Tab. 1

4.4 Components and interfaces

4.4.1 Desktop scanner


Fig. 1 Machine overview, front

- 1 Scan base plate
- 2 Swivelling axis for scan focus adjustment


Fig. 2 Machine overview, rear

- 1 USB port (control system)
- 2 Power supply connection
- 3 USB port (camera)
- 4 Main switch
- 5 Functional earthing

4.4.2 Model holder


Fig. 3 Model fixation

- 1 Universal holder
- 2 Modelling wax
- 3 Fixation pins incl. fixation spring levers
- 4 Scan base plate

The scan base plate is stationary in the desktop scanner. The universal holder or alternatively models with Splitex counter plate are placed on it, depending on requirement.

4.4.3 PC system requirements

Minimum system requirements:

- _ Resolution of 1680 × 1050 pixels (22" monitor)
- _ Operating system, Windows 8 64 bit
- _ 2 GB of free space for the installation
- _ 2 × USB 2.0
- _ Internet connection

System recommendation:

- _ CPU Intel Xeon E5 2.8 GHz
- _ 16 GB memory
- _ 1 TB free disk space
- _ Workstation graphics card with at least 2 GB memory and OpenGL 4.4
- _ Resolution of 1920 × 1080 pixels

4.4.4 Software system requirements

Ceramill software version above version 3.3.X

5 Installation

5.1 Transport

- ▷ For transport, have two persons lift the equipment by the floor plate.
- ▷ If the packaging is damaged: Promptly notify the delivering forwarding agent as well as the dealer.

5.2 Unpacking

The device is delivered in a sturdy carton. A foam element inside the carton protects the scanner during transportation.

- ▷ Upon delivery, immediately check the packaging for any visible damage.
- ▷ Open upper folding lid of the carton.
- ▷ Firstly, lift foam element out of the carton.
- ▷ Afterwards, remove complete carton by lifting it upward and away.
- ▷ Carefully remove device from packaging and place it down at the intended, prepared workstation (see Section 5.3).

The scanner accessories are at the bottom of the packaging.

- ▷ Check completeness of delivery scope.

5.3 Setting up

- _ The machine is intended exclusively for use within dry, closed rooms.
- _ Minimum space requirements for the desktop scanner (incl. connections):
 - Width: 360 mm
 - Depth: 390 mm (excl. connection cable)
 - Height: 310 mm
- _ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Place device on a sturdy substructure (table, work bench, etc.) that can hold at least twice the scanner weight, as listed in the technical specification. Apart from the carrying capacity of the substructure, its stability and sturdiness is also decisive for safe operation of the scanner.
- ▷ The opening of the scan chamber must face away from windows and strong artificial lighting.

5.4 Installation of the software on the PC

 As each scanner is individually factory-calibrated, the data on the installation CD are only valid for the respective scanner. Prior to installation:

- ▷ Check if the serial number of the scanner on the type plate is identical with the serial number stated on the CD.

- ▷ Carry out the installation routine of the supplied software as administrator, and follow the installer instructions.

5.5 Installation of the scanner

 Prior to installation:

- ▷ Check if the serial number of the scanner on the type plate is identical with the serial number stated on the CD.

- ▷ Make sure that the main switch of the scanner is positioned on **0**.
- ▷ Connect the two enclosed USB cables to the scanner and the rear USB ports of the PC (USB 2.0).
- ▷ Connect the mains cable to the power supply connection of the scanner as well as to a mains socket outlet.
- ▷ Switch device on (set main switch of scanner to **1**).

5.6 Activating additional modules

Additional modules must be activated in the scanner software before they can be used.

- ▷ Click on *Application > Add Activation Code...* to enter the received scanner-specific activation code.


Fig. 4

- ▷ Restart the software to activate the module.

5.7 Calibration of axes

Following the set-up and connection to the PC, the desktop scanner must be calibrated using the supplied calibration model.

 The scanner should be recalibrated at regular intervals.
The device must be calibrated each time after transporting.

▷ Place calibration model on the scan base plate.


Fig. 5 Calibration model in scanner

- ▷ Create a new patient case in the data base and start the scanner software from the data base.
- ▷ Cancel the subsequent command prompt for the 2D scan.
- ▷ Click on *Scanner Options* > *Axis calibration* to start the axis calibration.


Fig. 6

In the subsequent window, the values #1 and #2 must correspond with those of the calibration model. If the values do not correspond, see Section 5.7.1.

▷ Start the axis calibration.


Fig. 7

5.7.1 Entering calibration object values

▷ To enter the calibration object values, click on *Application... > Register Calibration Object...*


Fig. 8

▷ Enter calibration object values and confirm with *OK*.


Fig. 9

▷ After entering the values, a calibration process must be carried out.

6 Application Tips

6.1 Universal holder

The universal holder is used for transferring models without Splitex holder/plate into the scanner. It can be used both to affix individual models as well as model pairs in occlusion.


Fig. 10 Universal holder

- 1 Top side
- 2 Metal pins
- 3 Holes for fixation pins

- _ The top side of the universal holder [1] is level; it is used to apply modelling wax onto which the dental models are affixed.
- _ Two metal pins are located at the rear side of the level surface [2]. They represent the rear limit for placement of the models.
- _ The universal holder comprises 16 holes [3]. These are used for placement of the fixation pins including the fixation spring levers.

6.1.1 Placement in the scanner

- ▷ Position the universal holder directly on the scan base plate in the scanner.

The integrated magnet provides sufficient fixation on the scan base plate without the requirement of further fixation elements.

- ▷ To remove the universal holder from the scanner, grasp it at the rear or at the sides and lift it out.

6.1.2 Fixation pins incl. fixation spring levers

The fixation pins incl. the fixation spring levers are used to affix model pairs in occlusion onto the universal holder.

- ▷ For this, insert the fixation pins into optimally matching holes until the metal ring faces against the universal holder.

The fixation pins are placed loosely into the holes and achieve their stability through proper use with the fixation spring levers.

NOTE:

Only insert the fixation pins into the universal holder when being used for model fixation. Otherwise, they may fall out during the scanning process and damage the scanner.

Function method

The fixation via the fixation spring levers takes place by means of canting/chocking the spring levers with the fixation pins.

- ▷ Position and align the fixation spring lever with respect to the object.
- ▷ Using one finger, centrally press/push down the on the spring lever.

This causes the fixation spring lever to deform. After pressing/pushing down the fixation spring lever, it cants/chocks with the fixation pins and thus firmly clamps the object against the universal holder.


Fig. 11 Use of fixation spring levers

6.2 Positioning and affixing single models

Always position models in such a manner that they are within the scan focus. In the scanner, the scan focus is located on the level of the tilting axis (see figure 12).


Fig. 12 Positioning in the scan focus

6.2.1 Single model fixation using the universal holder

The placement of the model on the universal holder should be performed outside of the scanner.

- ▷ Apply modelling wax uniformly on the universal holder.
- ▷ Place the model centrally on the universal holder and press it firmly against the modelling wax.

To ensure that the model being scanned is sufficiently affixed to the universal holder and does not become loose and fall off during the scan process, the following must be observed:

- ▷ Press the model with both hands firmly against the modelling wax.
- ▷ Pay attention that the complete bottom surface of the model base has contact with the modelling wax.

NOTE:

The model must not project beyond the rear limitation (metal pins) of the universal holder.

- ▷ Afterwards, place the universal holder with the model onto the scan base plate of the scanner and press it against the scan base plate.


Fig. 13 Model fixation on universal holder

NOTE:

When using modelling wax not supplied by the manufacturer, it can occur that the model becomes loose/detaches during the scan process.

Insufficient scan results or damage of mechanical components and/or the model are possible.

- ▷ Make sure that the modelling wax provides sufficient holding force during the scan process.

6.2.2 Single model fixation directly on the scan base plate

Models with Splitex counter plate can be placed directly on the scan base plate in the scanner.

NOTE:

- ▷ Observe scan focus!

6.3 Positioning and affixing model pairs

Always position model pairs in such a manner that they are within the scan focus. In the scanner, the scan focus is located on the level of the tilting axis (see figure 14).


Fig. 14 Positioning in the scan focus

Model pairs are affixed exclusively using the universal holder. Model pairs can be affixed either with a rubber band or with use of the fixation spring levers.

The placement of the model pairs on the universal holder should be performed outside of the scanner.

6.3.1 Placement of model pairs with integrated articulation elements

Model pairs with integrated articulation elements (e.g. Triple-Tray), at which the jaw articulation is already determined, can be directly attached to the universal plate.

- ▷ Apply modelling wax uniformly on the universal holder.

- ▷ Place the model pair centrally on the universal holder and press it firmly against the modelling wax.

6.3.2 Placement of model pairs without integrated articulation elements

For model pairs without integrated articulation elements, it is recommended to firstly fasten only the lower jaw to the universal plate.

- ▷ Apply modelling wax uniformly on the universal holder.
- ▷ Place the lower jaw model centrally on the universal holder and press it firmly against the modelling wax.
- ▷ Afterwards, bring the upper jaw model in occlusion to the affixed lower jaw model.

6.3.3 Fixation with rubber band

For the fixation of model pairs with a rubber band, a suitable and appropriate household rubber band can be used. The contact surface of the universal plate to the scan base plate has an opening that provides sufficient space for a rubber band.

- ▷ Place the model pair onto the universal plate.
- ▷ Fasten/tension suitable rubber band around model pair and universal plate.

NOTE:

The rubber band must be tensioned tight enough, to prevent the models from loosening or moving during the scan process.

The contact surfaces of the Splitex plate may not be enclosed by the rubber band, as this would influence the hold or grip of the universal plate in the scanner.


Fig. 15 Fixation with rubber band

NOTE:

Tearing rubber bands during scanning can lead to damage of the scanner and the models. Especially old rubber bands can become brittle and easily tear.

- ▷ It is recommended to use several strong, new rubber bands.

6.3.4 Fixation with spring levers

After placing the model pairs on the universal plate, mount the fixation spring levers (see Section 6.1.2).

NOTE:

It is recommended to use at least two fixation spring levers to achieve a uniform distribution of the hold-down pressure and to prevent the model pairs from tilting.

The position of the fixation pins and the fixation spring levers can be adjusted.

- ▷ Insert the respective fixation pin into a hole of the universal holder that is suitable for fixation of the model.


Fig. 16 Use of fixation spring levers

7 Cleaning and Maintenance

7.1 Cleaning

The scanner must be cleaned regularly.

- ▷ Switch the scanner off.

NOTE:

Possible damage to the 3D sensor:

- ▷ Do **not** clean the optics of the 3D sensor in the top of the scanner!

- ▷ Carefully vacuum out the interior of the scanner.
- ▷ Wipe off the housing with a moist cloth. Do not use aggressive cleaning agents.

7.2 Maintenance

As the scanner is a sensitive, optical device, the maintenance may be carried out only by trained, specialized personnel.

8 Malfunctions, Repairs and Warranty

8.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the scanner.
- ▷ Restart the PC.
- ▷ Check the cable connections; disconnect and then re-connect, as required.

In case of malfunctions that are not rectified by restarting the scanner and the software:

- ▷ Contact the customer service.

8.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

8.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

9 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the scanner via the public disposal systems in an environmentally-friendly manner.

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Amann Girrbach AG
Amann Girrbach GmbH