Hand Scanner & IntelliCode

DE	Betriebsanleitung Handscanner und Data Matrix Code	3-10
EN	Operating instructions hand scanner, Data Matrix Code	11-18
FR	Mode d'emploi scanner à main, Data Matrix Code	19-26
IT	Istruzioni per l'uso Serie scanner manuale, Data Matrix Code	27-34
ES	Instrucciones de uso Serie escáner de mano, Data Matrix Code	35-42

1. INFORMATION ZUR ANLEITUNG

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
1. INFORMATION ABOUT THE INSTRUCTIONS

These are brief instructions for integrating the hand scanner into the CAD/CAM System from Amann Girrbach. The original instructions of the manufacturer Honeywell are also included with the product.

2. INFORMATION ABOUT THE HAND SCANNER

The hand scanner is used for scanning blank-specific codes (so-called IntelliCodes) and for transferring the data contained in the code to the CAM software.

2.1 Installation

- _ Switch off the PC before connecting the scanner.
- _ Connect the hand scanner to the PC via the USB port.
- _ Restart the PC once the scanner has been properly connected.

Fig. 1: Connect the hand scanner

2.2 Scan the barcode (reading process)

Fig. 2: Reading process of the aiming beam

An aiming beam is emitted via the viewfinder. This should be directed over the middle of the barcode.

Hold the scanner over the barcode, press the key and direct the aiming beam on the middle of the barcode.

Hold the scanner closer to the barcode with smaller barcodes and further away from the barcode with larger barcodes.

With a highly reflective barcode (e.g. due to a coating) it may be necessary to hold the scanner at an inclined angle, so that the barcode will be read.

2.3 Country-specific keypad

Fig. 3: Programme country-specific keypad

Scan the barcode **“Programme country-specific keypad”**. (Fig. 3)

Country code	Scanning
Belgium	1
Denmark	8
Finland	2
France	3
Germany / Austria	4
Great Britain	7
Italy	5
Norway	9
Switzerland	6
USA (Standard)	0

Fig. 4: Country codes

Then scan the **numeric barcode** (according to the country code, see Fig. 4) in the programme chart (Page 13).

Note:

A complete list of the country codes can be found in the operating instructions of the product on the Honeywell website www.honeywell.com/search#q=VG1400

Fig. 5: Save

Then scan the **save** barcode (Fig. 5), to programme the keypad for your country.

As a rule, the following scanner symbols for countries other than the USA are not supported: @|{ } [] = / \ < > ~

Programming chart

0

1

2

3

4

5

6

7

8

9

2.4 Resetting hand scanners to standard product settings

Fig. 6: Factory settings

All settings for standard products are reset to factory settings using the **Factory settings** (Fig. 6) barcode.

3. INTELLICODE

The IntelliCode (Data Matrix Code) is on all Amann Girrbach materials. It contains blank-specific information such as serial number, batch number, material type, height, shrinkage factor, shade and shape.

The IntelliCode is either on a large packaging label (see Fig. 7) or on a small packaging label (see bottom Fig.). The choice of which label is used is based on internal process technological reasons.

Fig. 7: Packaging label

The IntelliCode (1) can be scanned using the hand scanner (recommended). Alternatively the letters/number sequence (2) below or next to the IntelliCode can be entered manually.


Fig. 8: Packaging label

Note:

The IntelliCode on the packaging label must be scanned for easy, quick creation of blanks in the CAM software. This is the only way to transfer blank information to the CAM software.

3.3 Creating new blanks in the CAM software


Fig. 10

1) Creating a round or horseshoe-shaped blank in the CAM software

Open the CAM software and click in the field with the cursor **3** (Fig. 10)


Fig. 11

Scanning the IntelliCode on the blank packaging:

- _ Scan the IntelliCode using the hand scanner
- _ Click on "Generate" and "Use immediately" (Fig. 11)


Fig. 12

(2) Creating a block (e.g. VITABLOCS Mark II) in the CAM software:

- _ Select the required material in the database (Fig. 12)
- _ Click in the field with the cursor in **3**
- _ Scan the IntelliCode on the packaging


Fig. 13

(3) Creating a titanium abutment blank in the CAM software:

- _ Select the CAM software in the “Product code” menu with a right mouse click in the background. (Fig. 13)
- _ Click in the field with the cursor **3** (Fig. 14)
- _ Scan the IntelliCode on the blank packaging


Fig. 14

3.4 Find blanks again that have already been created in the CAM software


Fig. 15

Open the CAM software and click in the field with the cursor ③ (Fig. 15)

Fig. 16

The message that the product code already exists appears in the software when scanning the code on the blank label (Fig. 16). The blank already created can now be loaded using the "Load" button. (Fig. 17)

Fig. 17

Note: The same message would appear if the code on the packaging label had been scanned, if still available.

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Amann Girrbach AG