SCAN Accessories

DE Betriebsanleitung	3- 7
EN User Manual	8- 13
FR Mode d'emploi	14- 19
IT Istruzioni d'uso	20- 25
ES Modo de empleo	26- 30

Universalhalter

Universal holder

Support universel

Supporto universale

Soporte universal

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Tab. 1

Table of Contents

1	Explanation of Symbols	8
2	Suitable Personnel	8
3	Product Description	9
4	Application Tips	9
4.1	Universal holder	9
4.2	Positioning and affixing single models	11
4.3	Positioning and affixing model pairs	11
5	Accessories	13

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

3 Product Description


Fig. 6 Model fixation

- 1 Universal holder (delivery scope)
- 2 Modelling wax
- 3 Fixation pins incl. fixation spring levers
- 4 Scan plate

The scan plate is stationary in the desktop scanner. The universal holder or alternatively models with Splitex counter plate are placed on it, depending on requirement.

4 Application Tips

4.1 Universal holder

The universal holder is used for transferring models without Splitex holder/plate into the scanner. It can be used both to affix individual models as well as model pairs in occlusion.


Fig. 7 Universal holder

- 1 Top side
- 2 Metal pins
- 3 Holes for fixation pins

- _ The top side of the universal holder [1] is level; it is used to apply modelling wax onto which the dental models are affixed.
- _ Two metal pins are located at the rear side of the level surface [2]. They represent the rear limit for placement of the models.
- _ The universal holder comprises 16 holes [3]. These are used for placement of the fixation pins including the fixation spring levers.

4.1.1 Placement in the scanner

- ▷ Position the universal holder directly on the scan plate in the scanner.

The integrated magnet provides sufficient fixation on the scan plate without the requirement of further fixation elements.

- ▷ To remove the universal holder from the scanner, grasp it at the rear or at the sides and lift it out.

4.1.2 Fixation pins incl. fixation spring levers

The fixation pins incl. the fixation spring levers are used to affix model pairs in occlusion onto the universal holder.

- ▷ For this, insert the fixation pins into optimally matching holes until the metal ring faces against the universal holder.

The fixation pins are placed loosely into the holes and achieve their stability through proper use with the fixation spring levers.

NOTE:

Only insert the fixation pins into the universal holder when being used for model fixation. Otherwise, they may fall out during the scanning process and damage the scanner.

Function method

The fixation via the fixation spring levers takes place by means of canting/chocking the spring levers with the fixation pins.

- ▷ Position and align the fixation spring lever with respect to the object.
- ▷ Using one finger, centrally press/push down the on the spring lever.

This causes the fixation spring lever to deform. After pressing/pushing down the fixation spring lever, it cants/chocks with the fixation pins and thus firmly clamps the object against the universal holder.


Fig. 8 Use of fixation spring levers

4.2 Positioning and affixing single models

Always position models in such a manner that they are within the scan focus.

The placement of the model on the universal holder should be performed outside of the scanner.

- ▷ Apply modelling wax uniformly on the universal holder.
- ▷ Place the model centrally on the universal holder and press it firmly against the modelling wax.

To ensure that the model being scanned is sufficiently affixed to the universal holder and does not become loose and fall off during the scan process, the following must be observed:

- ▷ Press the model with both hands firmly against the modelling wax.
- ▷ Pay attention that the complete bottom surface of the model base has contact with the modelling wax.

NOTE:

The model must not project beyond the rear limitation (metal pins) of the universal holder.

- ▷ Afterwards, place the universal holder with the model onto the scan plate of the scanner and press it against the scan plate.

NOTE:

When using modelling wax not supplied by the manufacturer, it can occur that the model becomes loose/detaches during the scan process.

Insufficient scan results or damage of mechanical components and/or the model are possible.

- ▷ Make sure that the modelling wax provides sufficient holding force during the scan process.

4.3 Positioning and affixing model pairs

Always position model pairs in such a manner that they are within the scan focus.

Model pairs can be affixed either with a rubber band or with use of the fixation spring levers.

The placement of the model pairs on the universal holder should be performed outside of the scanner.

4.3.1 Placement of model pairs with integrated articulation elements

Model pairs with integrated articulation elements (e.g. Triple-Tray), at which the jaw articulation is already determined, can be directly attached to the universal plate.

- ▷ Apply modelling wax uniformly on the universal holder.
- ▷ Place the model pair centrally on the universal holder and press it firmly against the modelling wax.

4.3.2 Placement of model pairs without integrated articulation elements

For model pairs without integrated articulation elements, it is recommended to firstly fasten only the lower jaw to the universal plate.

- ▷ Apply modelling wax uniformly on the universal holder.
- ▷ Place the lower jaw model centrally on the universal holder and press it firmly against the modelling wax.
- ▷ Afterwards, bring the upper jaw model in occlusion to the affixed lower jaw model.

4.3.3 Fixation with rubber band

For the fixation of model pairs with a rubber band, a suitable and appropriate household rubber band can be used. The contact surface of the universal plate to the scan plate has an opening that provides sufficient space for a rubber band.

- ▷ Place the model pair onto the universal plate.
- ▷ Fasten/tension suitable rubber band around model pair and universal plate.

NOTE:

The rubber band must be tensioned tight enough, to prevent the models from loosening or moving during the scan process. The contact surfaces of the Splitex plate may not be enclosed by the rubber band, as this would influence the hold or grip of the universal plate in the scanner.


Fig. 9 Fixation with rubber band

NOTE:

Tearing rubber bands during scanning can lead to damage of the scanner and the models. Especially old rubber bands can become brittle and easily tear.

- ▷ It is recommended to use several strong, new rubber bands.

4.3.4 Fixation with spring levers

After placing the model pairs on the universal plate, mount the fixation spring levers (see Section 4.1.2).

NOTE:

It is recommended to use at least two fixation spring levers to achieve a uniform distribution of the hold-down pressure and to prevent the model pairs from tilting.

The position of the fixation pins and the fixation spring levers can be adjusted.

- ▷ Insert the respective fixation pin into a hole of the universal holder that is suitable for fixation of the model.


Fig. 10 Use of fixation spring levers

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Amann Girrbach AG