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- Translation of the original Operating Instructions -

Table of Contents

Explanation of Symbols	20
General Safety Instructions	22
Suitable Personnel	22
Properties	22
Warranty/Exclusion of Liability	26
Testing the Blanks	26
Mounting to the Worktable	26
Application	27
Environmental Protection	34
Accessories	34
Download information	34

EN

Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Date of manufacture
	Used by YYYY-MM-DD or YYYY-MM
Rx only	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.

General Safety Instructions

**CAUTION:**

Possible health impairment from CoCrMo dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).

Suitable Personnel

**NOTE:**

The product may only be used by trained dental technicians.

Properties

Intended use

CoCrMo blanks for manufacturing permanent and removable dental prosthetics.

Product description

Ceramill Sintron are blanks made of a CoCrMo alloy.

Ceramill Sintron blanks are used for manufacturing carrier structures for permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/ telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g., Ceramill Motion).

The material dries in the blank condition, meaning, it is machined without the use of coolant/lubricant (C/L); afterwards, it is final-sintered exclusively in the high-temperature sintering furnace (Ceramill Argotherm), inclusive of a special sintering structure (Ceramill Argovent) at a pre-programmed temperature that is adapted to this material, while under admission of argon.

**NOTE:**

A different sintering furnace than the Ceramill Argotherm may not be used. Should a different sintering furnace be used, it cannot be ensured that the technical properties of the sintered frames are achieved.

Ceramill Sintron is a Class IIa medical product; upon completion of the specified final-sintering, it meets the requirements of a Type 4 dental alloy according to DIN ISO 22674 for permanent and removable dentures.

Safety data sheet/Declaration of conformity

The data sheet on safety and the declaration of conformity are available on request under *www.amanngirrbach.com*.

Technical data

	Unit	Value
Tensile strength (R_m)	MPa	900
0.2 % proof stress ($R_{p0,2}$)	MPa	450
E-module (E)	GPa	180
Ductile yield	%	30
Vickers hardness	HV 10	270
Thermal expansion coefficient (CTE) (25 - 500°C)	1/K	14.5×10^{-6}
Density	g/cm^3	7.9
Open porosity	%	0
Colour	–	Silver
Oxidation colour	–	Grey-green
Body-compatibility:		
Corrosion resistance	–	DIN EN ISO 10271:2011-10
Bio-compatibility	–	ISO 10993-1

Chemical composition

Sintron	Mass percentage
Cobalt (Co)	66
Chromium (Cr)	28
Molybdenum (Mo)	5
Further elements (Mn, Si, Fe)	< 1
Further elements (C)	< 0.1
Organic binder (for blanks in blank condition)	1 - 2

The alloy is free of nickel, beryllium, gallium and cadmium according to DIN EN ISO 22674.

Shelf life of the blanks

The Ceramill Sintron blanks can be used 3 years after the manufacturing date.

Storage

Store Ceramill Sintron blanks in the original packaging in a dry location.

Store already milled blanks in a plastic bag together with a desiccant bag.

Packaging sizes

Art. No.	Designation	Description	Package
761132	Ceramill Sintron 98x10 N	CoCr blank, 98 × 10 mm	Pack of 1
761130	Ceramill Sintron 98x12 N	CoCr blank, 98 × 12 mm	Pack of 1
761128	Ceramill Sintron 98x14 N	CoCr blank, 98 × 14 mm	Pack of 1
761126	Ceramill Sintron 98x16 N	CoCr blank, 98 × 16 mm	Pack of 1
761124	Ceramill Sintron 98x18 N	CoCr blank, 98 × 18 mm	Pack of 1
761122	Ceramill Sintron 98x20 N	CoCr blank, 98 × 20 mm	Pack of 1
761120	Ceramill Sintron 98x25 N	CoCr blank, 98 × 25 mm	Pack of 1

Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the Ceramill Sintron blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

Application

Course information and an instruction video on manufacturing large-span bridges are available under www.amanngirrbach.com.

Indication ranges

- _ Anatomically reduced and fully anatomical crown and bridge frames in the anterior and posterior tooth range
- _ Bridge frames with a maximum of two connected intermediate units in the anterior and posterior region, and a maximum anatomical length of 50 mm
- _ Cantilever bridges with a maximum of one bridge unit (maximum one free-end pontic and no further than the second premolar)

Contraindications

- _ Known incompatibilities with respect to the components

Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating Ceramill Sintron frames in high-density sintered condition:

	Minimum frame thickness in mm	Connector cross-section in mm ²		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill Sintron	0.4	≥ 7	≥ 9	2	2

Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm ²
				incisal/occlusal	Circular	
Primary components/Double crowns	–	1	–	0.5	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	> 7
	OXXO	4	2	0.7	0.5	> 9
Posterior bridge	OXO	3	1	0.5	0.5	> 9
	OXXO	4	2	0.7	0.5	> 12
Cantilever bridge	OOX	3	Free-end pontic	0.7	0.5	> 12
	OOOX	4				
	OXOX	4	1 + free-end pontic	0.7	0.5	> 12

O Abutment crown

X Bridge unit

Procedure when fabricating dental framework

The enlargement factor determined by Amann Girrbach for the Ceramill Sintron blanks is provided on the blanks with a batch-specific coding.

NOTE:

Insufficient milling results!

- ▷ Machine the Ceramill Sintron blanks exclusively with the dry-milling procedure, without the use of coolant/lubricant (C/L).

Preparation for the final sintering

- ▷ Remove the Ceramill Sintron frames from the blank using a cross-toothed carbide cutter.
- ▷ After removal from the blank, brush off adhering milling dust from the frame with a suitable brush.

NOTE:

In the blank condition (unsintered), the frame may not be steamed off or subjected to other moist environments (e.g. ultrasonic bath) under any circumstances.

- ▷ Place the Ceramill frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.
- ▷ Mount sintering bowl and sintering aids together.
- ▷ Place sintering bowl incl. the sintering aids into the furnace.

Final sintering

For the sintering, use only the Ceramill Argotherm high-temperature furnace with protective gas connection (please also refer to Ceramill Argotherm Operating Instructions).

NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

NOTE:

Faulty sintering results!

- ▷ To ensure proper sintering results, pay attention to the correct arrangement of the sintering aids. Please also refer to the Ceramill Argotherm Operating Instructions.
- ▷ Perform the final sintering of the frames in the Ceramill Argotherm.
 - Program duration: approx. 5 hours
 - ▷ Upon completion of the program, remove the sintering aids using pliers.

**NOTE:**

Faulty sintering results!

- ▷ Normally, the sintered bridges have a silver colour. When the frame or certain areas show signs of oxidation, select the sintering result using the document for evaluation of the sintering results (see Ceramill M-Center homepage, protected download area).

Post-processing

Post-processing fully anatomical dentures

After the final sintering has taken place:

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 μm aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework/adapt the frame using a cross-toothed carbide cutter.
- ▷ For fully anatomical restorations, polish using commercially available polishing aids (rubber polishers, brushes, NEM polishing paste).

Post-processing/preparation for the veneering

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 μm aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework the frame using a cross-toothed carbide cutter.
- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 μm aluminium oxide at a pressure of 3 - 4 bar.

- ▷ Grasp the frame with artery forceps; do not touch with fingers.

- ▷ Steam off frame with steam cleaner (e. g. Steamer X3 (116910)).
- ▷ Oxidation firing is not necessary. Should an optical check nevertheless be desired, perform oxidation firing at 980 °C with 1 minute dwell time. Afterwards, thoroughly blast off one more time.

Finish

- ▷ Gloss up metal edges using silicone polish. Prevent heat generation.

Laser welding

- ▷ Prepare connecting location in the form of a modified X seam.
- ▷ Blast off and clean the spot to be laser welded (using 110 - 250 µm pure aluminium oxide).
- ▷ Use CoCr laser welding wire (721130 ff).

After laser welding the Ceramill Sintron, the workpiece can be ceramic veneered.

Soldering

- ▷ Solder using CoCrMo filler (781630 Giosolder) and the corresponding flux.

After soldering the Ceramill Sintron, the workpiece can be ceramic veneered.

Veneer ceramics

As a general rule, all conventional metal veneer ceramics can be used.

- ▷ Observe and adhere the respective manufacturer's information!

- ▷ Match the thermal expansion coefficient (CTE) of the ceramic to the CTE of the alloy.

When using Creation (788000) ceramic material, do not apply bonder:

- ▷ Use Pulveropaker opaque powder.

Fixation

Conventional cementation

Due to the high strength and stability of Ceramill Sintronframes, conventional fixation with commercially available cement is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Environmental Protection

Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Disposal of the blanks

May not be disposed of together with household waste.

Do not dispose of into sewer systems.

Always dispose of according to official regulations.

Accessories

For more information on product-specific accessories, please refer to www.amanngirrbach.com.

Download information

Further instructions/manuals can be downloaded under www.ceramill-m-center.com and/or www.amanngirrbach.com.

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The logo for Ceramill Sintron, featuring a stylized icon of three horizontal curved lines to the left of the text "ceramill sintron" with a registered trademark symbol.
Amann Girrbach AG
Amann Girrbach GmbH