ceramill sintron® multi

**Betriebsanleitung**  
**Instruction Manual**  
**Instructions d'utili-**  
**sation**  
**Istruzioni per l'uso**  
**Instrucciones de uso**

Deutsch	3 - 17
English	18 - 32
Français	33 - 47
Italiano	48 - 62
Español	63 - 77

AMANN GIRRIBACH

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

## **Table of contents**

<b>Explanation of Symbols</b> .....	<b>19</b>
<b>General Safety Instructions</b> .....	<b>20</b>
<b>Suitable Personnel</b> .....	<b>20</b>
<b>Properties</b> .....	<b>21</b>
<b>Warranty/Exclusion of Liability</b> .....	<b>24</b>
<b>Testing the Blanks</b> .....	<b>24</b>
<b>Inserting into Lava™ Milling Machines</b> .....	<b>24</b>
<b>Application</b> .....	<b>25</b>
<b>Download Infos</b> .....	<b>32</b>

## Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

\_ **NOTE** means that property damage can occur.

\_ **CAUTION** means that light to fairly serious personal injury can occur.

\_ **WARNING** means that serious personal injury can occur.

\_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

**Other symbols in the Manual**

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

**General Safety Instructions****CAUTION:**

Possible health impairment from CoCrMo dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).

**Suitable Personnel****NOTE:**

The product may only be used by trained dental technicians.

## Properties

### Product description

Ceramill Sintron are blanks made of a CoCrMo alloy.

Ceramill Sintron blanks are used for manufacturing carrier structures for permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/ telescopic crowns, supra-constructions, etc.) using Lava™ milling machines.

The material dries in the blank condition, meaning, it is machined without the use of coolant/lubricant (C/L); afterwards, it is final-sintered exclusively in the high-temperature sintering furnace (Ceramill Argo-therm), inclusive of a special sintering structure (Ceramill Argovent) at a pre-programmed temperature that is adapted to this material, while under admission of argon.

#### NOTE:

A different sintering furnace than the Ceramill Argotherm may not be used. Should a different sintering furnace be used, it cannot be ensured that the technical properties of the sintered frames are achieved.

Ceramill Sintron is a Class IIa medical product; upon completion of the specified final-sintering, it meets the requirements of a Type 4 dental alloy according to DIN ISO 22674 for permanent and removable dentures.

**Safety data sheet/Declaration of conformity**

The data sheet on safety and the declaration of conformity are available on request under [www.amanngirrbach.com](http://www.amanngirrbach.com).

**Technical Data**

	<b>Unit</b>	<b>Value</b>
Tensile strength ( $R_m$ )	MPa	900
0.2 % proof stress ( $R_{p0,2}$ )	MPa	450
E-module (E)	GPa	180
Ductile yield	%	30
Vickers hardness	HV 10	270
Thermal expansion coefficient (CTE) (25 - 500°C)	1/K	$14.5 \times 10^{-6}$
Density	$g/cm^3$	7.9
Open porosity	%	0
Colour	–	Silver
Oxidation colour	–	Grey-green
Body-compatibility:		
Corrosion resistance	–	DIN EN ISO 10271:2011-10
Bio-compatibility	–	ISO 10993-1

## Chemical composition

Sintron	Mass percentage
Cobalt (Co)	66
Chromium (Cr)	28
Molybdenum (Mo)	5
Further elements (Mn, Si, Fe)	< 1
Further elements (C)	< 0.1
Organic binder (for blanks in blank condition)	1 - 2

The alloy is free of nickel, beryllium, gallium and cadmium according to DIN EN ISO 22674.

## Shelf life of the blanks

The Ceramill Sintron blanks can be used 1.5 years after the manufacturing date.

## Storage

Store Ceramill Sintron blanks in the original packaging in a dry location.

Store already milled blanks in a plastic bag together with a desiccant bag.

## Packaging sizes

Art. No.	Designation	Description	Package
761151	Ceramill Sintron Multi (15.3 mm)	CoCr blank, height = 15.3 mm	Pack of 2
761153	Ceramill Sintron Multi (18 mm)	CoCr blank, height = 18 mm	Pack of 2

## Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

## Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the Ceramill Sintron blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

## Inserting into Lava™ Milling Machines

Before processing Ceramill Sintron, the milling chamber of the milling machine must be cleaned.

The blanks are already bonded in the normal blank holder for Lava™ systems.

- ▷ Remove the label with the scaling factor from the blank.
- ▷ Insert the blank as usual into Lava™ systems.

## Application

### Indication ranges

- \_ Anatomically reduced and fully anatomical crown and bridge frames in the anterior and posterior tooth range.
- \_ Bridge frames with a maximum of two connected intermediate units in the anterior and posterior region, and a maximum anatomical length of 50 mm.
- \_ Cantilever bridges with a maximum of one bridge unit (maximum one free-end pontic and no further than the second premolar).

### Contraindications

- \_ Known incompatibilities with respect to the components

### Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating Ceramill Sintron frames in high-density sintered condition:

	Minimum frame thickness in mm	Connector cross-section in mm <sup>2</sup>		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill Sintron	0.4	≥ 7	≥ 9	2	2

## Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm <sup>2</sup>
				incisal/occlusal	circular	
Primary components/double crowns	–	1	–	0.5	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	>7
	OXXO	4	2	0.7	0.5	>9
Posterior bridge	OXO	3	1	0.5	0.5	>9
	OXXO	4	2	0.7	0.5	>12
Cantilever bridge	OOX	3	Free-end pontic	0.7	0.5	>12
	OOOX	4				
	OXOX	4	1 + free-end pontic	0.7	0.5	>12

- O** Abutment crown  
**X** Bridge unit

A material file with these parameters is available for the Dental Wings CAD programs.

## Procedure when fabricating dental framework

For selection of the blank, the parameters and the ID code are read into the Lava™ CAM software as usual. The individual scaling factor for Ceramill Sintron is printed on the blank label and must be entered manually into the respective field of the blank dialogue.

The screenshot shows a 'New blank' dialog box with the following fields and values:

- Enter Par Code:** 7800000000055292
- Par Code:** 7800000000055292
- Scaling X/Y/Z:** 1,1205
- Uniform scaling**
- Blank type data:**
  - Material:** Sintron
  - Geometry:** Sintron Multi 18 mm
  - Color:** 00
- Enter ID Code:** 1944904538351809
- ID Code:** 1944904538351809

Buttons: OK, Cancel

Fig. 2

▷ Fit the Lava™ milling machines with the tools for Lava™ Wax. These are:

- Type 11
- Type 12
- Type 13
- Type 18

### Preparation for the final sintering

- ▷ Remove the Ceramill Sintron frames from the blank using a cross-toothed carbide cutter.
- ▷ After removal from the blank, brush off adhering milling dust from the frame with a suitable brush.

#### NOTE:

In the blank condition (unsintered), the frame may not be steamed off or subjected to other moist environments (e.g. ultrasonic bath) under any circumstances.

- ▷ Place the Ceramill frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.
- ▷ Mount sintering bowl and sintering aids together.
- ▷ Place sintering bowl incl. the sintering aids into the furnace.

## Final sintering

For the sintering, use only the Ceramill Argotherm high-temperature furnace with protective gas connection (please also refer to Ceramill Argotherm Operating Instructions).

### NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

### NOTE:

Faulty sintering results!

- ▷ To ensure proper sintering results, pay attention to the correct arrangement of the sintering aids. Please also refer to the Ceramill Argotherm Operating Instructions.
- ▷ Perform the final sintering of the frames in the Ceramill Argotherm.
    - Program duration: approx. 5 hours
  - ▷ Upon completion of the program, remove the sintering aids using pliers.

**NOTE:**

Faulty sintering results!

- ▷ Normally, the sintered bridges have a silver colour. When the frame or certain areas show signs of oxidation, select the sintering result using the document for evaluation of the sintering results (see Ceramill M-Center homepage, protected download area).

**Post-processing****Post-processing fully anatomical dentures**

After the final sintering has taken place:

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 µm aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework/adapt the frame using a cross-toothed carbide cutter.
- ▷ For fully anatomical restorations, polish using commercially available polishing aids (rubber polishers, brushes, NEM polishing paste).

**Post-processing/preparation for the veneering**

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 µm aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework the frame using a cross-toothed carbide cutter.
- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250 µm aluminium oxide at a pressure of 3 - 4 bar.

- ▷ Grasp the frame with artery forceps; do not touch with fingers.

- ▷ Steam off frame with steam cleaner (e. g. Steamer X3 (116910)).
- ▷ Oxidation firing is not necessary. Should an optical check nevertheless be desired, perform oxidation firing at 980 °C with 1 minute dwell time. Afterwards, thoroughly blast off one more time.

## Finish

- ▷ Gloss up metal edges using silicone polish. Prevent heat generation.

## Laser welding

- ▷ Prepare connecting location in the form of a modified X seam.
- ▷ Blast off and clean the spot to be laser welded (using 110 - 250 µm pure aluminium oxide).
- ▷ Use CoCr laser welding wire (721130 ff).

After laser welding the Ceramill Sintron, the workpiece can be ceramic veneered.

## Soldering

- ▷ Solder using CoCrMo filler (781630 Giosolder) and the corresponding flux.

After soldering the Ceramill Sintron, the workpiece can be ceramic veneered.

## Veneer ceramics

As a general rule, all conventional metal veneer ceramics can be used.

- ▷ Observe and adhere the respective manufacturer's information!

- ▷ Match the thermal expansion coefficient (CTE) of the ceramic to the CTE of the alloy.

When using Creation (788000) ceramic material, do not apply bonder:

- ▷ Use Pulveropaker opaque powder.

## Fixation

### Conventional cementation

Due to the high strength and stability of Ceramill Sintronframes, conventional fixation with commercially available cement is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

## Download Infos

Further instructions/manuals and video tutorials can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com) and/or [www.amanngirrbach.com](http://www.amanngirrbach.com).

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**Amann Girrbach AG**
**Amann Girrbach GmbH**