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www.amangirrbach.com/instruction-manuals.

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Table of Contents

Explanation of Symbols	19
General Safety Instructions	21
Suitable Personnel	21
Properties	21
Warranty/Exclusion of Liability	24
Testing the Blanks	24
Mounting to the Worktable	25
Application	25
Environmental Protection	30
Accessories	31
Download information	31

Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Date of manufacture
	Used by YYYY-MM-DD or YYYY-MM
Rx only	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.

General Safety Instructions

CAUTION:

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

Suitable Personnel

NOTE:

The product may only be used by trained dental technicians.

Properties

Intended use

Zirconium-oxide blanks for permanent and removable dental prosthetics.

Product description

Ceramill ZI are blanks made of zirconium oxide (Y-TZP ZrO_2) for type II, class 6 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g. Ceramill Motion) or manual copy milling machines (e.g. Ceramill Base and Multi-x). After the specified final sintering, Ceramill ZI meets the requirements of the standard DIN EN ISO 6872:2009, and after appropriate further processing (see page 26 ff.) the requirements under ISO 13356:2008 are also complied with.

Safety data sheet/Declaration of conformity

The data sheet on safety and the declaration of conformity are available on request under www.amanngirrbach.com.

Technical data

	Unit	Value
Bending strength (4-point)	MPa	> 1200
E-module	GPa	> 200
Grain size	µm	≤ 0.6
Density	g/cm ³	≥ 6.07
Open porosity	%	0
Thermal expansion coefficient (CTE) (25 - 500°C)	1/K	10.4 ± 0.5 × 10 ⁻⁶
Chemical solubility	µg/cm ²	< 5
Radioactivity	Bq/g	< 0.2

Chemical composition

Oxide	Mass percentage
ZrO ₂ + HfO ₂ + Y ₂ O ₃	> 99.0
Y ₂ O ₃	4.5 - 5.6
HfO ₂	< 5
Al ₂ O ₃	< 0.5
Other oxides	< 0.5

Shelf life of the blanks

When appropriately stored, Ceramill ZI zirconium-oxide blanks can be used five years after the manufacturing date.

Storage

Store Ceramill ZI in the original packaging in a dry location.

Packaging sizes

Art. No.:	Designation	Description	Package
760101	Ceramill ZI Pre-forms	Zirconium-oxide blank, 16.5 x 16 mm	Pack of 12
760102	Ceramill ZI Preforms L	Zirconium-oxide blank, 16.5 x 20 mm	Pack of 12
760110	Ceramill ZI TC L	Zirconium-oxide blank, 16.5 x 20 mm	Pack of 12
760116	Ceramill ZI 31	Zirconium-oxide blank, 40 x 20 x 16 mm	Pack of 6
760143	Ceramill ZI 51L	Zirconium-oxide blank, 65 x 30 x 20 mm	Pack of 3
760172	Ceramill ZI 71 XS	Zirconium-oxide blank, dental arch form h = 12 mm	Pack of 1

Art. No.:	Designation	Description	Package
760173	Ceramill ZI 71 S	Zirconium-oxide blank, dental arch form h = 14mm	Pack of 1
760174	Ceramill ZI 71	Zirconium-oxide blank, dental arch form h = 16mm	Pack of 1
760176	Ceramill ZI 71 M	Zirconium-oxide blank, dental arch form h = 18mm	Pack of 1
760184	Ceramill ZI 71 L	Zirconium-oxide blank, dental arch form h = 20mm	Pack of 1
760175	Ceramill ZI 71 XL	Zirconium-oxide blank, dental arch form h = 25mm	Pack of 1
760181	Ceramill ZI 77L	Zirconium-oxide blank, 98 x 20 mm, for Ceramill Base	Pack of 1
760183	Ceramill ZI 77XL	Zirconium-oxide blank, 98 x 25 mm, for Ceramill Base	Pack of 1
760186	Ceramill ZI 77L	Zirconium-oxide blank, 98 x 20 mm, for Ceramill Multi-x	Pack of 1
760187	Ceramill ZI 77XL	Zirconium-oxide blank, 98 x 25 mm, for Ceramill Multi-x	Pack of 1

Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the Ceramill ZI blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

Application

Indication ranges

- _ Anatomically reduced crown and bridge frames in the anterior and posterior tooth range.
- _ Bridge frames with a maximum of three connected intermediate units in the anterior and two connected intermediate links in the posterior region, and a maximum anatomical length of 50 mm.
- _ Cantilever bridges with a maximum of one bridge pontic (maximum one free-end pontic and no further than the second premolar).

Indication limitation for Canada:

- _ Single crowns
- _ Front-tooth bridges
- _ Posterior bridges with up to four units

Contraindications

- _ Insufficient tooth-structure availability
- _ Insufficient preparation results
- _ Insufficient oral hygiene
- _ More than two connected bridge units in the posterior region
- _ Known incompatibilities with respect to the components

Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating Ceramill ZI zirconium oxide frames in high-density sintered condition:

	Minimum frame thickness in mm	Connector cross-section in mm ²		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill ZI	0.5	≥ 7	≥ 9	3	2

Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm ²
				incisal/occlusal	Circular	
Primary components/ Double crowns	–	1	–	0.7	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	> 7
	OXXO	4	2	0.7	0.5	> 9
	OXXXO	5	3			
Posterior bridge	OXO	3	1	0.7	0.5	> 9
	OXXO	4	2	1.0	0.7	> 12
Cantilever bridge	OOX	3	Free-end pontic	1.0	0.7	> 12
	OOOX	4				
	OXOX	4	1+ Free-end pontic	1.0	0.7	> 12

O Abutment crown

X Bridge unit

Procedure when fabricating dental framework

Observe zirconium-oxide-efficient design forms and do not create pointed edges.

The enlargement factor determined by AmannGirrbach for the Ceramill ZI blanks is provided on the blanks with a batch-specific coding.

On the Ceramill Motion:

- ▷ Enter the enlargement factor in the Ceramill Match CAM software.
- ▷ Mill the blanks with the Ceramill Motion applying the milling strategies provided for the material.

Further details on designing and manufacturing framework with the Ceramill Mind CAD software are available on the Ceramill M-Center homepage in the protected download area.

On Ceramill copy milling machines:

- ▷ Set the enlargement factor on the copy milling machines.
- ▷ Follow the procedures and milling strategies for zirconium oxide given in the training documents.

For framework design information on the Ceramill copy milling machines, please refer to the training documents and the processing instructions of the modelling plastics Ceramill Gel/Pontic (760514/760522).

- ▷ Depending on requirement, the frames can be partially or completely coloured with Ceramill Liquid dye solutions (see processing instructions for Ceramill Liquid dye solutions).

Final sintering

The Ceramill Therm high-temperature furnace (178300/178350) or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

- ▷ Place the Ceramill ZI frames for sintering into the sintering bowl (178360) filled with sintering pearls (178311).
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
- Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min
 - Dwell time at final temperature: 2 hrs
 - Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended (recommended grain approx. 40 μm).

Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 μm , pressure \leq 2 bar) and subsequent steaming off.

Veneer ceramics

Ceramill ZI is suitable for veneering with conventional zirconium-oxide veneer ceramics (e.g., Creation ZI-CT). For information on the thermal-expansion coefficient of long-term cooled Ceramill ZI blanks, please refer to the Technical data (see page 22).

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Fixation

Conventional cementation:

Due to the high strength and stability of Ceramill ZI frames, conventional fixation with zinc-oxide phosphate cement or glass ionomer cement (e.g., GC Fuji Plus) is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Adhesive fixation:

For adhesive fixation, we recommend the bonding composite Panavia® 21 or Panavia® F 2.0, owing to the excellent bond strength with the zirconium-oxide frame material.

Provisional cementation is not recommended, as the frames can be subject to damage upon removal.

Environmental Protection

Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

Accessories

The dye solutions Ceramill Liquid CL 1 - CL 4 are used for staining Ceramill ZI blanks! The instruction manual is can be found under www.amanngirrbach.com/instruction-manuals.

For more information on product-specific accessories, please refer to www.amanngirrbach.com.

Download information

Further instructions/manuals and video tutorials can be downloaded under www.ceramill-m-center.com and/or www.amanngirrbach.com.

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Amann Girrbach AG
Amann Girrbach GmbH