# **ceramill<sup>®</sup> zi**

<b>DE</b> Gebrauchsanweisung	3 - 15
<b>EN</b> Instruction Manual	16 - 27
<b>FR</b> Instructions d'utilisation	28 - 40
<b>FR</b> Istruzioni per l'uso	41 - 53
<b>ES</b> Instrucciones de uso	54 - 66

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

## **Table of Contents**

<b>Explanation of Symbols .....</b>	<b>17</b>
<b>General Safety Instructions .....</b>	<b>19</b>
<b>Suitable Personnel .....</b>	<b>19</b>
<b>Properties .....</b>	<b>19</b>
<b>Warranty/Exclusion of Liability .....</b>	<b>21</b>
<b>Testing the Blanks .....</b>	<b>21</b>
<b>Mounting to the Worktable .....</b>	<b>22</b>
<b>Application .....</b>	<b>22</b>
<b>Environmental Protection .....</b>	<b>27</b>
<b>Accessories .....</b>	<b>27</b>
<b>Download information .....</b>	<b>27</b>

## Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

## Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

## Other symbols on the product

Symbol	Meaning
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Used by YYYY-MM-DD or YYYY-MM
<b>Rx only</b>	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.

## General Safety Instructions

### CAUTION:

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

## Suitable Personnel

### NOTE:

The product may only be used by trained dental technicians.

## Properties

### Intended use

Zirconium-oxide blanks for permanent and removable dental prosthetics.

### Product description

The blanks are made of zirconium oxide (Y-TZP  $ZrO_2$ ) for type II, class 5 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing carrier structures for permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/ telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g., Ceramill Motion).

After the specified final sintering, the material meets the requirements of the standard DIN EN ISO 6872, and after appropriate further processing (see page 23 ff.) the requirements under ISO 13356 are also complied with.

### Safety data sheet/Declaration of conformity

The data sheet on safety and the declaration of conformity are available on request under [www.amanngirrbach.com](http://www.amanngirrbach.com).

### Technical data

	Unit	Value
Bending strength (3-point)	MPa	1200 ± 150
Bending strength (4-point)	MPa	1000 ± 150
E-module	GPa	> 200
Grain size	µm	≤ 0.6
Density	g/cm <sup>3</sup>	≥ 6.07
Open porosity	%	0
Thermal expansion coefficient (CTE) (25 - 500°C)	1/K	10.4 ± 0.5 × 10 <sup>-6</sup>
Chemical solubility	µg/cm <sup>2</sup>	< 5
Radioactivity	Bq/g	< 0.2

**Chemical composition**

Oxide	Mass percentage
$ZrO_2 + HfO_2 + Y_2O_3$	> 99.0
$Y_2O_3$	4.5 - 5.6
$HfO_2$	< 5
$Al_2O_3$	< 0.5
Other oxides	< 0.5

**Shelf life of the blanks**

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

**Storage**

Store the blanks in the original packaging in a dry location.

**Warranty/Exclusion of Liability**

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

**Testing the Blanks**

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

## Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

### Application

#### Indication ranges

- \_ Anatomically reduced crown and bridge frames in the anterior and posterior tooth range.
- \_ Bridge frames with a maximum of three connected intermediate units in the anterior and two connected intermediate links in the posterior region, and a maximum anatomical length of 50 mm.
- \_ Cantilever bridges with a maximum of one bridge pontic (maximum one free-end pontic and no further than the second premolar).

Indication limitation for Canada:

- \_ Single crowns
- \_ Front-tooth bridges
- \_ Posterior bridges with up to four units

#### Contraindications

- \_ Insufficient tooth-structure availability
- \_ Insufficient preparation results
- \_ Insufficient oral hygiene
- \_ More than two connected bridge units in the posterior region
- \_ Known incompatibilities with respect to the components

### Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating zirconium oxide frames in high-density sintered condition:

	Minimum frame thickness in mm	Connector cross-section in mm <sup>2</sup>		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill ZI	0.5	≥ 7	≥ 9	3	2

### Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm <sup>2</sup>
				incisal/occlusal	Circular	
Primary components/ Double crowns	–	1	–	0.7	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	> 7
	OXXO	4	2	0.7	0.5	> 9
	OXXXO	5	3			
Posterior bridge	OXO	3	1	0.7	0.5	> 9
	OXXO	4	2	1.0	0.7	> 12
Cantilever bridge	OOX	3	Free-end pontic	1.0	0.7	> 12
	OOOX	4				
	OXXO	4	1+	1.0	0.7	> 12
			Free-end pontic			

O Abutment crown

X Bridge unit

### Procedure when fabricating dental framework

Observe zirconium-oxide-efficient design forms and do not create pointed edges.

The determined enlargement factor for the blanks is provided on the blanks.

- ▷ Enter the enlargement factor in the CAM software.
- ▷ Mill the blanks applying the milling strategies provided for the material.
- ▷ Depending on requirement, the frames can be partially or completely coloured with Ceramill Liquid dye solutions (see processing instructions for Ceramill Liquid dye solutions).

### Final sintering

The Ceramill Therm high-temperature furnace or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

**NOTE:**

▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
- Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min
  - Dwell time at final temperature: 2 hrs
  - Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

**Post-processing**

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended (recommended grain approx. 40 μm).

**Frame preparation for veneering**

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110μm, pressure ≤ 2 bar) and subsequent steaming off.

### Veneer ceramics

Ceramill ZI is suitable for veneering with conventional zirconium-oxide veneer ceramics (e.g., Creation ZI-CT). For information on the thermal-expansion coefficient of long-term cooled Ceramill ZI blanks, please refer to the Technical data (see page 20).

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

### Fixation

Conventional cementation:

Due to the high strength and stability of Ceramill ZI frames, conventional fixation with zinc-oxide phosphate cement or glass ionomer cement (e.g., GC Fuji Plus) is possible in most cases.

▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Adhesive fixation:

For adhesive fixation, we recommend the bonding composite Panavia® 21 or Panavia® F 2.0, owing to the excellent bond strength with the zirconium-oxide frame material.

Provisional cementation is not recommended, as the frames can be subject to damage upon removal.

## Environmental Protection

### Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

## Accessories

For more information on product-specific accessories, please refer to [www.amanngirrbach.com](http://www.amanngirrbach.com).

## Download information

Further instructions/manuals and video tutorials can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com) and/or [www.amanngirrbach.com](http://www.amanngirrbach.com).

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**Amann Girrbach AG**
**Amann Girrbach GmbH**