ceramill® zirconia

 **ceramill® zi**

 **ceramill® zolid**

 **ceramill® zolid**
preshades

 **ceramill® zolid ht+**
white

 **ceramill® zolid ht+**
preshades

 **zolid gen-x**
multilayer

DE Gebrauchsanweisung	3 - 24
EN Instruction Manual	25 - 46
FR Instructions d'utilisation	47 - 68
IT Istruzioni per l'uso	69 - 90
ES Instrucciones de uso	91 - 110

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

Table of Contents

Explanation of Symbols	26
General Safety Instructions	28
Suitable Personnel	28
Properties	29
Warranty/Exclusion of Liability	32
Testing the Blanks	32
Mounting to the Worktable	32
Application	33
Fixation	44
Post-processing at the dentist's	45
Follow-up examination at the dentist's	45
Environmental Protection	46
Accessories	46
Download information	46

EN

Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
REF	Catalog number
LOT	Batch code
	Manufacturer
	Consult instructions for use
	Use by YYYY-MM-DD or YYYY-MM
Rx only	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.
SN	Serial number

General Safety Instructions

**CAUTION:**

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

**NOTE:**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Suitable Personnel

**NOTE:**

The product may only be used by trained dental technicians.

Properties

Intended use

Zirconium-oxide blanks for permanent and removable dental prosthetics.

Product Description

The blanks are made of zirconium oxide (Y-TZP ZrO_2) for type II, class 5 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing permanent and removable prosthetic restorations (e.g. crowns and bridges, conical/telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g. Ceramill milling machines) or manual copy milling machines (e.g. Ceramill Base and Multi-x). After the specified final sintering, the material meets the requirements of the standard DIN EN ISO 6872.

Safety data sheet/Declaration of conformity

The data sheet on safety and the declaration of conformity are available on request.

Technical data

	Unit	Ceramil		Ceramil		Ceramil		Ceramil		Ceramil		Zolid	
		Zi	Zolid / Zolid	Zolid	HT+ White	Zolid	HT+	Zolid	HT+	Zolid	HT+	Zolid	Gen-X Multilayer
Preshades													
Bending strength													
_3-point	MPa	1200±150	1100±150	1100±150	1100±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150
_4-point	MPa	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150	1000±150
E-module	GPa	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200	≥ 200
Thermal expansion coefficient (CTE) (25 - 500 °C)	10 ⁻⁶ /K	10.4±0.5	10.8±0.5	10.8±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.4±0.5	10.5±0.5
Chemical solubility	µg/cm ²	< 100	< 100	< 100	< 100	< 100	< 100	< 100	< 100	< 100	< 100	< 100	< 100
Vickers hardness	HV10	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200	1300±200

The delivery scope can vary country-specific.

Chemical composition

Oxide	Mass percentage							
	Ceramill Zi	Ceramill Zolid	Ceramill Zolid Pre-shades	Ceramill Zolid HT+	Ceramill White	Ceramill Zolid HT+ Preshades	Ceramill Zolid Gen-X	Zolid Multilayer
ZrO ₂ + HfO ₂ + Y ₂ O ₃	≥ 99,0	≥ 99,0	≥ 99,0	≥ 99,0	≥ 99,0	≥ 99,0	≥ 99,0	≥ 99,0
Y ₂ O ₃	4,5 - 5,6	4,5 - 5,6	4,5 - 5,6	6,7 - 7,2	6,0 - 7,0	6,0 - 7,0	6,0 - 7,0	6,0 - 7,0
HfO ₂	≤ 5	≤ 5	≤ 5	≤ 5	≤ 5	≤ 5	≤ 5	≤ 5
Al ₂ O ₃	≤ 0,5	≤ 0,5	≤ 0,5	≤ 0,5	≤ 0,5	≤ 0,5	≤ 0,5	≤ 0,5
andere Oxide	≤ 1	≤ 1	≤ 1	≤ 1	≤ 1	≤ 1	≤ 1	≤ 1

Shelf life of the blanks

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

Storage

Store the blanks in the original packaging in a dry location.

Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

Application

Indication ranges

- _ Anatomically reduced crown and bridge frames in the anterior and posterior tooth range, and monolithic (fully anatomical) crowns and bridges.
- _ Anatomically reduced four to multi-unit bridge frames with a maximum of three connected intermediate units in the anterior region and a maximum of two connected intermediate links in the posterior region
- _ Monolithic four to multi-unit bridges with a maximum of three connected intermediate units in the anterior region and a maximum of two connected intermediate links in the posterior region
- _ Cantilever frames and bridges with a maximum of one bridge unit (maximum one free-end pontic and no further than the second premolar).

Indication limitation for Canada:

- _ Single crowns
- _ Front-tooth bridges
- _ Posterior bridges with up to four units

Contraindications

- _ Insufficient tooth-structure availability
- _ Insufficient preparation results
- _ Insufficient oral hygiene
- _ More than two connected bridge units in the posterior region, more than three connected intermediate units in the anterior region
- _ Known incompatibilities with respect to the components
- _ Heavily discoloured hard tooth structure
- _ provisional insertion

Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating zirconium oxide frames in high-density sintered condition:

Minimum frame thickness in mm	Connector cross-section in mm ²		Maximum number of connected bridge units	
	anterior region	posterior region	anterior region	posterior region
0.5	≥ 7	≥ 9	3	2

Detailed minimum wall thicknesses and connector cross-sections

Indication	Maximum number of bridge units	Wall thickness in mm		Connector cross-section in mm ²
		incisal/occlusal	circular	
Primary components/ Double crowns	–	0.7	0.5	–
Single crown	–	0.5	0.5	–
Front-tooth bridge	1	0.5	0.5	>7
	3	0.7	0.5	>9
Posterior bridge	1	0.7	0.5	>9
	2	1.0	0.7	>12
Cantilever bridge	1	1.0	0.7	>12


Crown border parameters in the CAD software

▷ Pay attention to sufficient crown border thicknesses Recommended values:

- 1. Crown border thickness (Horizontal): 0.2 mm
- 2. Angled: 0.3 mm
- 3. Angle: 60°
- 4. Height (Vertical): 0 mm

For tangential preparations, increase the “crown border thickness” and “Angled” parameters, e.g.:

- 1. Crown border thickness (Horizontal): 0.5 - 0. mm
- 2. Angled: 0.35 - 0.4 mm


Observe zirconium-oxide-efficient design forms and do not create pointed edges.


Positioning in the CAM software

- ▷ Apply enough sufficiently thick holding ridges (at least 3 ridges for a single crown; thickness > 1.7 mm).


- ▷ Position holding ridges approx. 1 - 2 mm away from the crown margin, especially for long crown margins.

This provides support to the crown margin.

- ▷ Position the crown in the blank in such a manner that the clearance from the crown margin to the blank surface is at least 1 mm.

This ensures that the crown has no contact to the sinter skin of the blank.

Procedure when fabricating dental framework

For maximum fitting accuracy, the specific expansion/shrinkage value is given on the blanks.

When machining the blanks in milling machines, it is recommended to only use cutters with sufficient cutting edge quality.

- ▷ Depending on requirement in the respective CAM software, enter the appropriate value.

In case of uncertainties, contact your CAM supplier.

- ▷ Mill the blanks with the milling machine applying the milling strategies provided for the material.

Removal from the blank

Utmost care must be taken when removing the frames from the blank:

- ▷ Remove the frames only using suitable tools approved for zirconium oxide (avoid cutting discs).

For trimming and dressing the connectors and smoothing out the frames, we recommend use suitable tools approved for zirconium oxide.

For detailed information on the subject “Working zirconium oxide” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Colouration of the white zirconium oxide

▷ Depending on requirement, partially or completely stain the white zirconium oxide frames prior to the sintering with dye solution.

For detailed information on the subject “Colouration and finalising” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Colouration of the pre-shaded zirconium oxide

These blanks are not stained with the dye solutions, but monolithically stained upon successful final sintering with stain or glazing material or used as a frame for veneering.

Zolid Gen-X Multilayer

These blanks are available in 16 A-D Vita shades + 2 bleach colours. The integrated colour gradient gives the multilayer blanks a natural appearance immediately after sintering. The finalising takes place using glazing material, stains and/or veneer ceramics.

Ceramill Zolid HT+ Preshade

These blanks are available in 16 A-D Vita shades. The monochrome-stained blanks can be finalised using glazing materials, stains and/or veneer ceramics.

Ceramill Zolid Preshade

Blank Shade	Tooth shade						
	A1	A2	A3	A3,5	B1	B2	B3
Ceramill Zolid Preshades 0 (Bleach)	Bleach						
Ceramill Zolid Preshades 2 (A light)	●	●	○	-	-	-	-
Ceramill Zolid Preshades 3 (A medium)	-	●	●	○	-	-	○
Ceramill Zolid Preshades 1 (B light)	○	○	-	-	●	●	-

● ideal for this shade

- not suitable for this shade

○ shade can be achieved

For detailed information on the subject “Colouration and finalising” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Final sintering

A high-temperature sintering furnace, e.g. Ceramill Therm or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

Highly translucent zirconium oxide and opaque zirconium oxide can be sintered together in one sintering process.

NOTE:

When highly translucent zirconium oxide and opaque zirconium oxide are sintered together in the same sintering bowl and with the same sintering pearls, then this can lead to changes in the translucency. Therefore, it is advisable for highly translucent zirconium oxide and opaque zirconium oxide, to use individual sintering pearls in separate sintering bowls.

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
 - Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min
 - Dwell time at final temperature: 2 hrs
 - Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

For detailed information on the subject "Sintering" please refer to the brochure "Processing Techniques Zolid". These can be found under www.amanngirrbach.com in the download area.

Frames for veneering

Post-processing

If required, the zirconium oxide frames can be reworked after the final sintering without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).

Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent steaming off.

Veneer ceramics

For information on the thermal-expansion coefficient, please refer to the Technical data (see page 30).

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Monolithic dentures

Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).
- ▷ Polish the contact surfaces to the antagonist and to the neighbouring teeth to a high gloss when in high-density sintered condition.

Frame preparation for staining and glazing

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent ultrasonic treatment and/or steaming off.

The exterior crown surfaces are not corundum blasted.

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Staining and glazing

- ▷ Stain and glaze the units afterwards.

For detailed information on the subject “Sintering” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Fixation

Conventional cementation: Due to the high strength and stability of the frames, conventional fixation with zinc-oxide phosphate cement or glass ionomer cement (e.g., GC Fuji Plus) is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Adhesive fixation: For adhesive fixation, we recommend the bonding composite Panavia[®] 21 or Panavia[®] F 2,0m owing to the excellent bond strength with the zirconium-oxide frame material.

Provisional cementation is not recommended, as the frames can be subject to damage upon removal.

Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40 μm). For polishing, we recommend using diamond polishing bits.

NOTE:

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

Follow-up examination at the dentist's

It is urgently recommended to check the monolithic dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

Accessories

For more information on product-specific accessories, please refer to www.amanngirrbach.com.

Download information

Further instructions/manuals can be downloaded under www.amanngirrbach.com/instruction-manuals.

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Amann Girrbach AG
Amann Girrbach GmbH