 **ceramill® zolid**  
fx preshade

<b>DE</b> Gebrauchsanweisung	3 - 19
<b>EN</b> Instruction Manual	20 - 36
<b>FR</b> Instructions d'utilisation	37 - 55
<b>FR</b> Istruzioni per l'uso	56 - 74
<b>ES</b> Instrucciones de uso	75 - 91

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

## **Table of Contents**

<b>Explanation of Symbols</b> .....	<b>21</b>
<b>General Safety Instructions</b> .....	<b>23</b>
<b>Suitable Personnel</b> .....	<b>23</b>
<b>Properties</b> .....	<b>23</b>
<b>Warranty/Exclusion of Liability</b> .....	<b>26</b>
<b>Testing the Blanks</b> .....	<b>26</b>
<b>Mounting to the Worktable</b> .....	<b>26</b>
<b>Application</b> .....	<b>27</b>
<b>Fixation</b> .....	<b>34</b>
<b>Post-processing at the dentist's</b> .....	<b>35</b>
<b>Follow-up examination at the dentist's</b> .....	<b>35</b>
<b>Environmental Protection</b> .....	<b>36</b>
<b>Accessories</b> .....	<b>36</b>
<b>Download information</b> .....	<b>36</b>

## Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

## Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

## Other symbols on the product

Symbol	Meaning
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Used by YYYY-MM-DD or YYYY-MM
<b>Rx only</b>	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.

## General Safety Instructions

### CAUTION:

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

## Suitable Personnel

### NOTE:

The product may only be used by trained dental technicians.

## Properties

### Intended use

Zirconium-oxide blanks for permanent and removable dental prosthetics.

### Product Description

The blanks are made of zirconium oxide ( $ZrO_2$ ) for type II, class 4 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing carrier structures for permanent and removable prosthetic restorations (e.g., crowns and three-unit bridges into the molar region) using CNC milling machines (e.g., Ceramill Motion).

These blanks are not stained with the Ceramill Liquids, but monolithically stained upon successful final sintering with Ceramill Stain&Glaze or used as frame for veneering. The following Table lists the recommendations for achieving the respective tooth shade.

Blank Shade	Tooth shade							
	Bleach	A1	A2	A3	A3,5	B1	B2	B3
Ceramill Zolid FX								
Preshade Bleach	●	-	-	-	-	-	-	-
Ceramill Zolid FX								
Preshade A light	-	●	●	○	-	-	-	-
Ceramill Zolid FX								
Preshade A medium	-	-	●	●	○	-	-	-
Ceramill Zolid FX								
Preshade B light	-	-	-	-	-	●	●	○

- ideal for this shade
- shade can be achieved
- not suitable for this shade

After completion of the specified end-sintering, the material meets the requirements of DIN EN ISO 6872.

**Safety data sheet/Declaration of conformity**

The data sheet on safety and the declaration of conformity are available on request under [www.amanngirrbach.com](http://www.amanngirrbach.com).

**Technical data**

	<b>Unit</b>	<b>Value</b>
Bending strength (3-point)	MPa	700±150
Bending strength (4-point)	MPa	600±150
E-module	GPa	> 200
Grain size	µm	≤ 0.8
Density	g/cm <sup>3</sup>	≥ 6.0
Open porosity	%	0
Vickers hardness	HV10	1400 ± 200
Thermal expansion coefficient (CTE) (25 - 500°C)	10 <sup>-6</sup> /K	10.1 ± 0.5
Chemical solubility	µg/cm <sup>2</sup>	< 100
Radioactivity	Bq/g	< 1.0

**Chemical composition**

<b>Oxide</b>	<b>Mass percentage</b>
ZrO <sub>2</sub> + HfO <sub>2</sub> + Y <sub>2</sub> O <sub>3</sub>	≥ 99.0
Y <sub>2</sub> O <sub>3</sub>	9.15 - 9.55
HfO <sub>2</sub>	< 5
Al <sub>2</sub> O <sub>3</sub>	< 0.5
Other oxides	< 1

**Shelf life of the blanks**

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

**Storage**

Store the blanks in the original packaging in a dry location.

**Warranty/Exclusion of Liability**

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

**Testing the Blanks**

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

**Mounting to the Worktable**

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

## Application

### Indication ranges

- \_ Anatomically reduced crown framework and three-unit bridge frames in the anterior and posterior tooth range
- \_ Monolithic (fully anatomical) crowns and up to three-unit bridges in the anterior and posterior tooth range
- \_ Anatomically reduced bridge frames with a maximum of one intermediate link in the anterior and posterior tooth range
- \_ Monolithic bridges with a maximum of one intermediate link in the anterior and posterior tooth range

### Contraindications

- \_ Insufficient tooth-structure availability
- \_ Insufficient preparation results
- \_ Insufficient oral hygiene
- \_ More than one bridge unit
- \_ Known incompatibilities with respect to the components
- \_ Heavily discoloured hard tooth structure

### Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating zirconium oxide frames in high-density sintered condition:

	Minimum frame thickness in mm	Connector cross-section in mm <sup>2</sup>		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill Zolid FX PS	0.5	≥ 12	≥ 12	1	1

### Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm <sup>2</sup>
				incisal/occlusal	circular	
Primary components/ Double crowns	–	1	–	0.7	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	≥ 12
Posterior bridge	OXO	3	1	0.7	0.5	≥ 12


- O** Abutment crown  
**X** Bridge unit

## Border parameters in the Ceramill Mind CAD software

Pay attention to sufficient crown border thicknesses Recommended values:

- \_ 1. Crown border thickness (Horizontal): 0.2 mm
- \_ 2. Angled: 0.3 mm
- \_ 3. Angle: 60°
- \_ 4. Height (Vertical): 0 mm

For tangential preparations, increase the “crown border thickness” and “Angled” parameters, e.g.:


- \_ 1. Crown border thickness (Horizontal): 0.5 - 0.3 mm
- \_ 2. Angled: 0.35 - 0.4 mm

Observe zirconium-oxide-efficient design forms and do not create pointed edges.

## Positioning in the CAM software

- ▷ Apply enough sufficiently thick holding ridges (at least 3 ridges for a single crown; thickness > 1.7 mm).


- ▷ Position holding ridges approx. 1 - 2 mm away from the crown margin, especially for long crown margins.

This provides support to the crown margin.

- ▷ Position the crown in the blank in such a manner that the clearance from the crown margin to the blank surface is at least 1 mm.

This ensures that the crown has no contact to the sinter skin of the blank.

### Procedure when fabricating dental framework

The determined enlargement factor for the blanks is provided on the blanks.

When processing Ceramill Zolid FX Preshades in the milling machine, it is recommended to use only mills, whose service life is still in the green range of the wear indication.

- ▷ Enter the enlargement factor in the CAM software.
- ▷ Mill the blanks applying the milling strategies provided for the material.

Further details on designing and manufacturing framework with the Ceramill Mind CAD software are available on the Ceramill M-Center homepage in the protected download area.

- ▷ Clean the frames from residual zirconia dust using a brush or compressed air.

### Removal from the blank

Utmost care must be taken when removing the frames from the blank:

- ▷ Separate the frames by grinding them out using a turbine grinder (without water-cooling) and a pointed diamond bit (do not use cutting discs).

For trimming and dressing the connectors and smoothing out the frames, we recommend using the Amann Girrbach polishing set.

### Final sintering

The Ceramill Therm high-temperature furnace or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

Highly translucent zirconium oxide and opaque zirconium oxide can be sintered together in one sintering process.

#### NOTE:

When highly translucent zirconium oxide and opaque zirconium oxide are sintered together in the same sintering bowl and with the same sintering pearls, then this can lead to changes in the translucency. Therefore, it is advisable for highly translucent zirconium oxide and opaque zirconium oxide, to use individual sintering pearls in separate sintering bowls.

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.

- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

**NOTE:**

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
- Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min
  - Dwell time at final temperature: 2 hrs
  - Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

## Frames for veneering

### Post-processing

If required, the zirconium oxide frames can be reworked after the final sintering without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).

### Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent steaming off.

## Veneer ceramics

Ceramill Zolid FX Preshades is suitable for veneering with conventional zirconium-oxide veneer ceramics (e.g., Creation ZI-CT). For information on the thermal-expansion coefficient of the blanks, please refer to the Technical data (see page 25).

## Monolithic dentures

### Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked using a water-cooled lab turbine and suitable diamond grinding points (grain approx. 40 µm).
- ▷ Polish the contact surfaces to the antagonist and to the neighbouring teeth to a high gloss when in high-density sintered condition. We recommend using the Amann Girrbach polishing set for this (please also refer to the processing instructions of the polishing set).

### Preparation for staining and glazing

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent ultrasonic treatment and/or steaming off.

The exterior crown surfaces are not corundum blasted.

For bridges with massive intermediate units, we recommend step-by-step long-term cooling to 500 °C.

### Staining and glazing

- ▷ Stain and glaze the units afterwards. We recommend using the Amann Girrbach stain and glaze kit for this (please also refer to the processing instructions of the stain and glaze kit).

## Fixation

Due to the high translucency of this zirconium oxide, the selection of cement shade and possible discolourations of the tooth stump are to be observed, as these may affect the final colour result.

Conventional cementation:

Due to the high strength and stability of the restorations, conventional fixation with zinc-oxide phosphate cement or glass ionomer cement (e.g., GC Fuji Plus) is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Adhesive fixation:

For adhesive fixation, we recommend the bonding composite Panavia<sup>®</sup> 21 or Panavia<sup>®</sup> F 2.0, owing to the excellent bond strength with the zirconium-oxide frame material.

Provisional cementation is not recommended, as the frames can be subject to damage upon removal.

## Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40  $\mu\text{m}$ ). For polishing, we recommend using the Amann Girrbach diamond polishing bits (please also refer to the processing instructions of the polishing kit).

### NOTE:

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

## Follow-up examination at the dentist's

It is urgently recommended to check the monolithic dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

## Environmental Protection

### Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

## Accessories

For more information on product-specific accessories, please refer to [www.amanngirrbach.com](http://www.amanngirrbach.com).

## Download information

Further instructions/manuals can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com) and/or [www.amanngirrbach.com](http://www.amanngirrbach.com).

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**Amann Girrbach AG**
**Amann Girrbach GmbH**