ceramill® zolid

Gebrauchsanweisung
Instruction Manual
Instructions d'utilisation
Istruzioni per l'uso
Instrucciones de uso

Deutsch	3 - 21
English	22 - 41
Français	42 - 61
Italiano	62 - 81
Español	82 - 101

AMANN GIRRIBACH

- Original Gebrauchsanweisung -

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www.amanngirrbach.com/instruction-manuals.

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- Translation of the original Operating Instructions -

Table of Contents

Explanation of Symbols	23
General Safety Instructions	25
Suitable Personnel	25
Properties	25
Warranty/Exclusion of Liability	28
Testing the Blanks	28
Mounting to the Worktable	28
Application	29
Fixation	39
Post-processing at the dentist's	39
Follow-up examination at the dentist's	40
Environmental Protection	40
Accessories	41
Download information	41

Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Date of manufacture
	Used by YYYY-MM-DD or YYYY-MM
Rx only	Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist.

General Safety Instructions

CAUTION:

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

Suitable Personnel

NOTE:

The product may only be used by trained dental technicians.

Properties

Intended use

Zirconium-oxide blanks for permanent and removable dental prosthetics.

Product description

Ceramill Zolid are blanks made of zirconium oxide (Y-TZP ZrO₂) for type II, class 6 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g. Ceramill Motion) or manual copy milling machines (e.g. Ceramill Base and Multi-x). After the specified final sintering, Ceramill Zolid meets the requirements of the standard DIN EN ISO 6872, and after appropriate further processing (see page 31 ff.) the requirements under ISO 13356 are also complied with.

Safety data sheet/Declaration of conformity

The data sheet on safety and the declaration of conformity are available on request under www.amanngirrbach.com.

Technical data

	Unit	Value
Bending strength (4-point)	MPa	1200±200
E-module	GPa	> 200
Grain size	µm	≤ 0.6
Density	g/cm ³	≥ 6.05
Open porosity	%	0
Thermal expansion coefficient (CTE) (25 - 500°C)	1/K	11.0 ± 0.5 × 10 ⁻⁶
Chemical solubility	µg/cm ²	< 5
Radioactivity	Bq/g	< 0.2

Chemical composition

Oxide	Mass percentage
$ZrO_2 + HfO_2 + Y_2O_3$	≥ 99.0
Y_2O_3	4.5 - 5.6
HfO_2	< 5
Al_2O_3	< 0.5
Other oxides	< 0.5

Shelf life of the blanks

When appropriately stored, Ceramill Zolid zirconium-oxide blanks can be used five years after the manufacturing date.

Storage

Store Ceramill Zolid in the original packaging in a dry location.

Packaging sizes

Art. No.	Designation	Description	Package
760162	Ceramill Zolid 71 XS	Zirconium-oxide blank, dental arch form, height = 12 mm	Pack of 1
760163	Ceramill Zolid 71 S	Zirconium-oxide blank, dental arch form, height = 14 mm	Pack of 1
760164	Ceramill Zolid 71	Zirconium-oxide blank, dental arch form, height = 16 mm	Pack of 1
760165	Ceramill Zolid 71 M	Zirconium-oxide blank, dental arch form, height = 18 mm	Pack of 1
760166	Ceramill Zolid 71 L	Zirconium-oxide blank, dental arch form, height = 20 mm	Pack of 1
760167	Ceramill Zolid 71 XL	Zirconium-oxide blank, dental arch form, height = 25 mm	Pack of 1

Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the Ceramill Zolid blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

Application

Indication ranges

- _ Anatomically reduced crown and bridge frames in the anterior and posterior tooth range, and monolithic (fully anatomical) crowns and bridges.
- _ Anatomically reduced bridge frames with a maximum of three connected intermediate units in the anterior and a maximum of two connected intermediate links in the posterior region, and a maximum anatomical length of 50 mm.
- _ Monolithic bridges with a maximum of two connected intermediate units in the posterior region, and a maximum anatomical length of 50 mm.
- _ Cantilever frames and bridges with a maximum of one bridge pontic (maximum one free-end pontic and no further than the second premolar)

Indication limitation for Canada:

- _ Single crowns
- _ Front-tooth bridges
- _ Posterior bridges with up to four units

Contraindications

- _ Insufficient tooth-structure availability
- _ Insufficient preparation results
- _ Insufficient oral hygiene
- _ More than two connected bridge units in the posterior region
- _ Known incompatibilities with respect to the components
- _ Heavily discoloured hard tooth structure

Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating Ceramill Zolid zirconium oxide frames in high-density sintered condition:

	Minimum frame thick- ness in mm	Connector cross-sec- tion in mm ²		Maximum number of connected bridge units	
		Anterior	Posterior	Anterior	Posterior
Ceramill Zolid	0.5	≥ 7	≥ 9	3	2

Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm ²
				incisal/occlusal	Circular	
Primary components/Double crowns	–	1	–	0.7	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	>7
	OXXO	4	2	0.7	0.5	>9
	OXXXO	5	3			
Posterior bridge	OXO	3	1	0.7	0.5	>9
	OXXO	4	2	1.0	0.7	>12
Cantilever bridge	OXX	3	Free-end pontic	1.0	0.7	>12
	OOOX	4				
	OXOX	4	1 + free-end pontic	1.0	0.7	>12

O Abutment crown

X Bridge unit


Border parameters in the Ceramill Mind CAD software

For Ceramill Zolid, it should be observed that the crown border thicknesses are sufficient. Recommended values:

- _ 1. Crown border thickness (Horizontal): 0.2 mm
- _ 2. Angled: 0.3 mm
- _ 3. Angle: 60°
- _ 4. Height (Vertical): 0 mm

For tangential preparations, increase the “crown border thickness” and “Angled” parameters, e.g.:


- _ 1. Crown border thickness (Horizontal): 0.25 - 0.3 mm
- _ 2. Angled: 0.35 - 0.4 mm


Observe zirconium-oxide-efficient design forms and do not create pointed edges.

Positioning in the CAM software

- ▷ Apply enough sufficiently thick holding ridges (at least 3 ridges for a single crown; thickness > 1.7 mm).


- ▷ Position holding ridges approx. 1-2 mm away from the crown margin, especially for long crown margins.

This provides support to the crown margin.

- ▷ Position the crown in the blank in such a manner that the clearance from the crown margin to the blank surface is at least 1 mm.

This ensures that the crown has no contact to the sinter skin of the blank.

Procedure when fabricating dental framework

The enlargement factor determined by Amann Girrbach for the Ceramill Zolid blanks is provided on the blanks with a batch-specific coding.

On the Ceramill Motion:

When processing Ceramill Zolid preshade in Ceramill Motion, it is recommended to use only mills, whose service life is still in the green range of the wear indication.

- ▷ Enter the enlargement factor in the Ceramill Match CAM software.
- ▷ Mill the blanks with the Ceramill Motion applying the (Zolid) milling strategies provided for the material.

Further details on designing and manufacturing framework with the Ceramill Mind CAD software are available on the Ceramill M-Center homepage in the protected download area.

On Ceramill copy milling systems:

- ▷ Set the enlargement factor on the Ceramill copy milling system.
- ▷ Follow the procedures and milling strategies for zirconium oxide given in the training documents.

For framework design information on the Ceramill copy milling systems, please refer to the training documents and the processing instructions of the modelling plastics Ceramill Gel/Pontic (760514/760522).

Removal from the blank

Utmost care must be taken when removing the frames from the blank:

- ▷ Separate the frames by grinding them out using a turbine grinder (without water-cooling) and a pointed diamond bit (do not use cutting discs).

For trimming and dressing the connectors and smoothing out the frames, we recommend using the Amann Girrbaach polishing set.

Colouring prior to sintering

NOTE:

Ceramill Zolid absorbs the dye liquid more intensely than Ceramill ZI. The frames get darker.

For more details on this and for further processing information, see the processing instructions of the Ceramill Liquid dye solutions.

- ▷ Clean the frames from residual zirconia dust using a brush or compressed air.
- ▷ Depending on requirement, the frames or the fully anatomical denture can be partially or completely coloured with Ceramill Liquid dye solutions (see processing instructions for Ceramill Liquid dye solutions).

Final sintering

The Ceramill Therm high-temperature furnace (178300/178350) or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

For sintering, a separate sintering bowl (178370) and sintering pearls (178311) are used.

**NOTE:**

Loss of translucence!

- ▷ The sintering bowl with the pearls are to be used only for Ceramill ZOLID and **not** for Ceramill ZI applications!
- ▷ Do not sinter Ceramill ZOLID and Ceramill ZI frames together in one sintering bowl!

Ceramill ZOLID frames and ZI frames may be sintered together in one sintering process in separate, stacked sintering bowls. Sintering bowls can only be stacked in the Ceramill Therm II (178350).

- ▷ Place the Ceramill Zolid frames for sintering into the sintering bowl (178370) filled with sintering pearls (178311).
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

**NOTE:**

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
 - Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min

- Dwell time at final temperature: 2 hrs
- Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

Frames for veneering

Post-processing

If required, the zirconium oxide frames can be reworked after the final sintering without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).

Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent steaming off.

Veneer ceramics

Ceramill Zolid is suitable for veneering with conventional zirconium-oxide veneer ceramics (e.g., Creation ZI-CT). For information on the thermal-expansion coefficient of Ceramill Zolid blanks, please refer to the Technical data (see page 26).

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Monolithic dentures

Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended (recommended grain approx. 40 µm).
- ▷ Polish the contact surfaces to the antagonist and to the neighbouring teeth to a high gloss when in high-density sintered condition. We recommend using the Amann Girrbach polishing set for this (please also refer to the processing instructions of the polishing set).

Frame preparation for staining and glazing

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 110 µm, pressure ≤ 2 bar) and subsequent ultrasonic treatment and/or steaming off.

The exterior crown surfaces are not corundum blasted.

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Staining and glazing

- ▷ Stain and glaze the units afterwards. We recommend using the Amann Girrbach stain and glaze kit for this (please also refer to the processing instructions of the stain and glaze kit).

Fixation

Conventional cementation:

Due to the high strength and stability of Ceramill Zolid frames, conventional fixation with zinc-oxide phosphate cement or glass ionomer cement (e.g., GC Fuji Plus) is possible in most cases.

▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

Adhesive fixation:

For adhesive fixation, we recommend the bonding composite Panavia® 21 or Panavia® F 2.0, owing to the excellent bond strength with the zirconium-oxide frame material.

Provisional cementation is not recommended, as the frames can be subject to damage upon removal.

Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40 µm). For polishing, we recommend using the Amann Girrbach diamond polishing bits (please also refer to the processing instructions of the polishing kit).

NOTE:

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

Follow-up examination at the dentist's

It is urgently recommended to check the monolithic dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

Environmental Protection

Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

Accessories

The dye solutions Ceramill Liquids new formula were developed especially for Ceramill Zolid. They have a changed formula and are available in the 16 Vita classical shades (A-D).

**NOTE:**

Do not mix Ceramill Liquids CL1-CL4 (Article No. 760470) with the new Ceramill Liquids A-D (Article No. 760469).

The instruction manual is can be found under
www.amangirrbach.com/instruction-manuals.

For more information on product-specific accessories, please refer to *www.amangirrbach.com*.

Download information

Further instructions/manuals can be downloaded under
www.ceramill-m-center.com and/or *www.amangirrbach.com*.

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Amann Girrbach AG
Amann Girrbach GmbH