Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones de
empleo

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Original-Gebrauchsanweisung

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INSTRUCTIONS FOR USE

Translation of the original instructions for use

Girosolder is a dental solder for cobalt-based alloys. It complies with DIN ISO 9333.

Contents: 5 soldering rods
Weight: 3.6 g

Girosolder

Melting range 1112 - 1138 °C
Working temperature 1160 °C

Composition in % by mass

Cobalt (Co)	70.3
Chromium (Cr)	18
Molybdenum (Mo)	5
Silicon (Si)	2.5
Manganese (Mn)	2.5
Boron (B)	1.7

No nickel (Ni), beryllium (Be) or cadmium (Cd) content in accordance with DIN ISO 9333.

Recommended fluxes

Fluxes suitable for soldering non-precious cobalt-based alloys.

Soldering instructions

PREPARING THE SURFACES FOR SOLDERING

Remove any obstructions which may hinder the flow of the solder, e.g. grooves caused by grinding. This involves creating very smooth metal surfaces at the solder joint to ensure it is optimally wetted with solder. Smooth the surfaces to be soldered with crosscut tungsten carbide cutters and, if necessary, rubber polishers and polish to a high lustre. The surface area of the solder joints should be as large as possible and parallel. The solder gap should not be greater than 0.2 mm to ensure a good capillary effect for the solder. When soldering cavities or margin extensions, place platinum foil under the defect being soldered (adapt to the model).

FABRICATING A SOLDERING MODEL

Ensure that the solder joint is kept as free as possible so that only minimum heat is lost. The soldering investment model should be as small as possible. Bench dry the soldering model (do not dry in the furnace).

CHEMICAL CONDITIONING OF THE SOLDER JOINT

Formation of oxides

Oxidation, which begins at low temperatures, can hinder efficient diffusion.

Prevention

Flux is required to protect the surface against oxidation and prevent the formation of oxide as far as possible. Apply flux to the solder joint and solder: both the solder and solder joint should be generously coated with flux to ensure successful metal bonding.

Adequate supply of heat


A hot, large soldering flame is also essential to ensure that the whole soldering area is uniformly heated to the correct flow temperature of the solder.

Heating the solder model

After bench drying the solder model, brush heat uniformly and thoroughly with a sufficiently large flame, then concentrate the heat on the solder joint without removing the flame. Apply the solder already precoated with flux to the solder joint. Do not remove the flame as this may cause oxygen to become entrapped. It is essential that the solder flows immediately. If it forms a ball, the metal is not hot enough. Add more flux and repeat the procedure.

Shelf life

Shelf life of up to 5 years from the date of manufacture.

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