Manufacturer | Hersteller

**Merz Dental GmbH**  
Eetzweg 20  
24321 Lütjenburg | Germany  
Fon +49 4381 403-0  
Fax +49 4381 403-403  
[www.merz-dental.de](http://www.merz-dental.de)

ISO 13485

Distribution | Vertrieb

**Amann Girrnbach AG**  
Herrschaftswiesen 1  
6842 Koblach | Austria  
Fon +43 5523 62333-105  
Fax +43 5523 62333-5119  
[austria@amanngirrnbach.com](mailto:austria@amanngirrnbach.com)

[www.amanngirrnbach.com](http://www.amanngirrnbach.com)

0482

Rx only

Made in the European Union

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## INSTRUCTIONS FOR USE Ceramill TEMP, TOOTH SHADES

Please read carefully! Caution: Federal law restricts this device to sale by or on the order of a licensed dental practitioner

### INTENDED USE

Ceramill TEMP blanks are milling blanks from which resin crowns, bridges and frameworks can be fabricated for dental restorations in temporary, pre-prosthetic treatment in the anterior and posterior regions.

### TECHNICAL DATA

PHYSICAL VALUES FOR CERAMILL TEMP

#### Test according to

Density	1,19 g/cm <sup>3</sup>	EN ISO 10477
Solubility	0,2 µg/mm <sup>3</sup>	EN ISO 10477
Flexural strength	91,5 MPa	EN ISO 10477
Water uptake	26,5 µg/mm <sup>3</sup>	EN ISO 10477
Modulus of elasticity	2773 MPa	EN ISO 20795-1:2009
Flexural strength	96,6 MPa	EN ISO 20795-1:2009

### COMPOSITION

- \_ Polymethyl methacrylate (PMMA) and methacrylic acid ester-based cross-linked polymers
- \_ Dye
- \_ Residual peroxide
- \_ Methyl methacrylate (MMA) may be contained up to a maximum of 1% as residual monomer

### SHADES

Available in 3 different dentine shades (light, middle, dark).

Blank	Tooth shade
light	A1; B1; D2; C1
middle	A2; A3; B2; B3; B4; C2; C3; D3; D4
dark	A3,5; A4; C4

### STORAGE

No special measures required.

### HAZARD WARNINGS

Dust is created during milling of Ceramill TEMP blanks that can cause irritation of the eyes, skin and airways. Therefore always ensure that the extractor on your milling machine is functioning perfectly. Further safety information can be obtained from the safety data sheet.

### USER INFORMATION

#### Design:

- \_ For crowns, do not fall below the minimum wall thickness of 1mm.
- \_ When designing bridges with one connected bridge unit, the connection cross-section between crown and connected bridge unit may not be less than 9mm<sup>2</sup>.
- \_ The user must decide in each case whether the PMMA - OMP-N® material guarantees adequate masticatory stability for the restoration according to the indication by selection of the appropriate luting system.

#### Veneering Ceramill TEMP frameworks:

- \_ A milled Ceramill TEMP framework can be faced using conventional veneering materials. MMA-based materials are the preferred choice for veneering.
- \_ When veneering with composite, Ceramill TEMP must always be bonded as described in the instructions of the respective manufacturer.

### FINISHING

- \_ Cross-cut tungsten carbide cutters or separating discs, which are suitable for acrylic, should be used for trimming and separating frameworks from the blank.
- \_ Excessive heat build-up should be avoided during finishing and polishing to prevent fitting inaccuracies after fabrication of the framework.

### RESTORING

Frameworks, full crowns and bridges, which have been fabricated using Ceramill TEMP blanks, can be restored and adjusted at any time using any conventional MMA-based cold-curing resin such as artConcept®.

### OPTIONS FOR LUTING INTRAORALLY

- \_ for temporary restorations: temporary cements
- \_ for long-term temporary restorations: glass ionomer cements without bonder

### SPECIAL INFORMATION

Enter the LOT No. with each procedure that requires identification of the material.

### DISPOSAL

Ceramill TEMP is insoluble in water, inactive, does not pose a risk for the ground water and can therefore be disposed of with normal household waste.

### SAFETY DATA SHEETS / DECLARATION OF CONFORMITY

Safety data sheets and declaration of conformity are available on request and under [www.amanngirrnbach.com](http://www.amanngirrnbach.com)

Information as of 2013-10

**The product properties are based on complying with and observing these instructions for use.**