ceramill® motion 2

ceramill® mikro

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

CHECKING THE BLANK ZERO POINT


The milling test with the Ceramill TEST blank is used to check the zero point and consequently the milling accuracy of the Ceramill milling machine. There are programmes stored in the software for milling a total of 8 test blocks. The test blocks should be arranged as follows on the test blank.

EN


Preparation

- _ The Ceramill Roto tools should be placed with the tool shank facing towards the outside in the Ceramill 71 blank holder.
- _ Place the Ceramill Roto 2.5 mm on tool position 1
- _ Place the Ceramill Roto 1.0 mm on tool position 2
- _ Insert the test blank Z1 71 in the blank holder and tighten using the screws.
- _ Log onto the PC as user with administrator rights.

Open the Ceramill Motion software in the start menu or from the desktop.


LOADING THE TEST BLOCK PROGRAMME AND MILLING

_Load the programme - by clicking on the folder symbol in the menu bar

_Select the calibration

There is a separate NC file for each test block.

5X: 8 test blocks can be milled from each blank.

5AXIS_TESTBLOC-1.NC to
5AXIS_TESTBLOC-8.NC

The two test blocks of a series that lie furthest apart must always be milled to check the zero point of the milling machine.

Test block No. 1 and test block No. 3 or
Test block No. 4 and test block No. 6

If the Z zero point is to be checked, all the other test blocks can also be used.

4X: 7 test blocks can be milled from each block.


Motion 2 (4X):

4AXIS_TESTBLOC-1.NC to
4AXIS_TESTBLOC-7.NC


Mikro:

CERAMILLMIKRO_TESTBLOC-1.NC to
CERAMILLMIKRO_TESTBLOC-7.NC

All test blocks can be used to check the zero point.

Transfer the programme to the milling machine.


Click the start button to start the milling procedure.

The machine begins with milling the test blocks.

Figure Ceramill Motion 2 (5X)

After the milling procedure is complete, twist the respective test block manually from the test blank or break it out.

Figure Ceramill Motion 2 (5X)

MEASURING THE TEST BLOCKS

Measure all test blocks using the caliper gauge supplied in the pack contents of the milling machine.

_Close the caliper gauge and set to 0

The following values must be measured and checked for each test block:

5X:

- _Length of the "straight side" - without number
- _Length of the "bevel side" - with number

4X:

- _Length of the "long side" - with number

Note the determined values. A table is included at the end of this document for entering the measurement values.


The measured values must lie in a range of 14.00 mm +/- 0.05 mm!

If the measured values do not correspond with the guidelines, contact the Ceramill Helpdesk.


DETERMINING THE TEST BLOCK OFFSET SITUATION

In rare cases positional errors of the axes may occur. This manifests itself by a visible and measurable offset in the region of the attachment points of the test blocks. Two different offset situations can occur.


Offset situation identification number above:

The front section of the test block with milled identification number is ABOVE the back section of the test block.

Offset situation identification number below:

The front section of the test block with milled identification number is BELOW the back section of the test block.

Measurement of the offset situation - for all milled test blocks respectively.


Figure example for 5X machine.
With 4X machines there is no bevel.

_ Measure the height of the test block in the middle section - Value = A

_ Measure the height of the test block in the back section - Value = B

_ Measure the height of the test block in the front section - Value = C

_ Note the determined measurement values

_ Give the determined measurement values to the Ceramill Helpdesk

Middle section - A = In the middle section of the test block the front section overlaps with the back section

Back section - B = milled section of the test block without milled identification number

Front section - C = middle section of the test block with milled identification number

Example:

Test block height A = 9.98 mm

Test block height B = 9.94 mm

Test block height C = 9.96 mm

Attachment: Table for determined measurement values

| Test block | Measurement 1 "without number" only (5X) | Measurement 2 "with number" | Measurement 3 Value = A | Measurement 4 Value = B | Measurement 5 Value = C |
|------------|--|--------------------------------|----------------------------|----------------------------|----------------------------|
| 1 | | | | | |
| 2 | | | | | |
| 3 | | | | | |
| 4 | | | | | |
| 5 | | | | | |
| 6 | | | | | |
| 7 | | | | | |
| 8 | | | | | |

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119

Distribution | Vertrieb D/A

Amann Girrbach GmbH
Dürrenweg 40
75177 Pforzheim | Germany
Fon +49 7231 957-100
Fax +49 7231 957-159

austria@amanngirrbach.com
germany@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

AMANNGIRRBACH