# **ceramil sintron®**

|                                      |         |
|--------------------------------------|---------|
| <b>DE</b> Gebrauchsanweisung         | 3 - 20  |
| <b>EN</b> Instruction Manual         | 21 - 37 |
| <b>FR</b> Instructions d'utilisation | 38 - 55 |
| <b>IT</b> Istruzioni per l'uso       | 56 - 73 |
| <b>ES</b> Instrucciones de uso       | 74 - 90 |

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

## Table of Contents

|  |    |
|--|----|
| Explanation of Symbols .....                 | 22 |
| General Safety Instructions .....            | 25 |
| Suitable Personnel .....                     | 25 |
| Properties .....                             | 25 |
| Warranty/Exclusion of Liability .....        | 28 |
| Testing the Blanks .....                     | 29 |
| Mounting to the Worktable .....              | 29 |
| Application .....                            | 29 |
| Post-processing at the dentist's .....       | 36 |
| Follow-up examination at the dentist's ..... | 36 |
| Environmental Protection .....               | 37 |
| Accessories .....                            | 37 |
| Download information .....                   | 37 |

EN

## Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

**Other symbols in the Manual**

| <b>Symbol</b> | <b>Meaning</b>   |
|---------------|--|
| ▷             | Item of an operation description   |
| –             | Item of a list   |
| ▪             | Subitem of an operation description or a list                            |
| [3]           | Numbers in square brackets refer to position numbers in graphics/figures |

## Other symbols on the product

| Symbol  | Meaning  |
|---|--|
|  | Medical device   |
|  | Unique Device Identifier   |
|  | Catalog number   |
|  | Batch code   |
|  | Manufacturer   |
|  | Consult instructions for use   |
|  | Use by YYYY-MM-DD or YYYY-MM   |
| <b>Rx only</b>  | Symbol that may be used in place of statement “CAUTION”. US Federal law restricts this device to sale by or on the order of a dentist. |
|  | Serial number  |

## General Safety Instructions

**CAUTION:**

Possible health impairment from CoCrMo dust!

▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).

**NOTE:**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Suitable Personnel

**NOTE:**

The product may only be used by trained dental technicians.

## Properties

### Intended use

CoCrMo blanks for manufacturing permanent and removable dental prosthetics.

### Product description

Ceramill Sintron are blanks made of a CoCrMo alloy.

They are used for manufacturing carrier structures for permanent and removable prosthetic restorations (e.g., crowns and bridges, conical/ telescopic crowns, supra-constructions, etc.) using CNC milling machines (e.g., Ceramill Motion).

The material dries in the blank condition, meaning, it is machined without the use of coolant/lubricant (C/L); afterwards, it is final-sintered exclusively in the high-temperature sintering furnace (Ceramill Argotherm), inclusive of a special sintering structure (Ceramill Argovent) at a pre-programmed temperature that is adapted to this material, while under admission of argon.

**NOTE:**

A different sintering furnace than the Ceramill Argotherm may not be used. Should a different sintering furnace be used, it cannot be ensured that the technical properties of the sintered frames are achieved.

Ceramill Sintron is a Class IIa medical product; upon completion of the specified final-sintering, it meets the requirements of a Type 4 dental alloy according to DIN EN ISO 22674 for permanent and removable dentures.

**Safety data sheet/Declaration of conformity**

The data sheet on safety and the declaration of conformity are available on request under *www.amangirrbach.com*.

**Technical data**

|   | <b>Unit</b> | <b>Value</b>          |
|---|-------------|-----------------------|
| Tensile strength ( $R_m$ )                          | MPa         | 900                   |
| 0.2 % proof stress ( $R_{p0,2}$ )                   | MPa         | 450                   |
| E-module (E)  | GPa         | 200                   |
| Ductile yield                                       | %           | 30                    |
| Vickers hardness                                    | HV 10       | 270                   |
| Thermal expansion coefficient (CTE)<br>(25 - 500°C) | 1/K         | $14.5 \times 10^{-6}$ |
| Density   | $g/cm^3$    | 7.9                   |
| Open porosity                                       | %           | 0                     |
| Colour  | –           | Silver                |
| Oxidation colour                                    | –           | Grey-green            |
| Body-compatibility:                                 |             |                       |
| Corrosion resistance                                | –           | DIN EN ISO 10271      |
| Bio-compatibility                                   | –           | ISO 10993-1           |

## Chemical composition

| Sintron  | Mass percentage |
|--|-----------------|
| Cobalt (Co)                                    | 66.0            |
| Chromium (Cr)                                  | 28.0            |
| Molybdenum (Mo)                                | 5.0             |
| Further elements (Mn, Si, Fe)                  | < 1             |
| Further elements (C)                           | < 0.1           |
| Organic binder (for blanks in blank condition) | 1 - 2           |

The alloy is free of nickel, beryllium, gallium and cadmium according to DIN EN ISO 22674.

## Shelf life of the blanks

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

## Storage

Store the blanks in the original packaging in a dry location.

Store already milled blanks in a plastic bag together with a desiccant bag.

## Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

## Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

## Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

## Application

Course information and an instruction video on manufacturing large-span bridges are available under [www.amanngirrbach.com](http://www.amanngirrbach.com).

## Indication ranges

- \_ Anatomically reduced and fully anatomical crown and bridge frames in the anterior and posterior tooth range
- \_ Bridge frames with a maximum of two connected intermediate units in the anterior and posterior region
- \_ Cantilever bridges with a maximum of one bridge unit (maximum one free-end pontic and no further than the second premolar)

## Contraindications

- \_ Known incompatibilities with respect to the components

## Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating Ceramill Sintron frames in high-density sintered condition:

|          | Minimum frame thickness<br>in mm | Connector cross-section in mm <sup>2</sup> |           | Maximum number of connected bridge units |           |
|----------|----------------------------------|--|-----------|--|-----------|
|          |                                  | Anterior                                   | Posterior | Anterior                                 | Posterior |
| Ceramill | 0.4                              | ≥ 7  | ≥ 9       | 2  | 2         |
| Sintron  |                                  |  |           |  |           |

## Detailed minimum wall thicknesses and connector cross-sections

| Indication                       | Number of bridge units | Wall thickness in mm |          | Connector crosssection in mm <sup>2</sup> |
|----------------------------------|------------------------|----------------------|----------|---|
|                                  |                        | incisal/occlusal     | Circular |   |
| Primary components/Double crowns | –                      | 0.5                  | 0.5      | –   |
| Single crown                     | –                      | 0.5                  | 0.5      | –   |
| Front-tooth bridge               | 1                      | 0.5                  | 0.5      | > 7                                       |
|                                  | 2                      | 0.7                  | 0.5      | > 9                                       |
| Posterior bridge                 | 1                      | 0.5                  | 0.5      | > 9                                       |
|                                  | 2                      | 0.7                  | 0.5      | > 12                                      |
| Cantilever bridge                | 1                      | 0.7                  | 0.5      | > 12                                      |

## Procedure when fabricating dental framework

The determined enlargement factor for the blanks is provided on the blanks.

- ▷ Enter the enlargement factor in the CAM software.

### NOTE:

Insufficient milling results!

- ▷ Machine the Ceramill Sintron blanks exclusively with the Ceramill Motion 2 using the dry-milling procedure without the use of coolant/lubricant (C/L).

- ▷ Mill the blanks applying the milling strategies provided for the material.

Further details on designing and manufacturing framework with the Ceramill Mind CAD software are available on the Ceramill M-Center homepage in the protected download area.

## Preparation for the final sintering

- ▷ Remove the frames from the blank using a cross-toothed carbide cutter.
- ▷ After removal from the blank, brush off adhering milling dust from the frame with a suitable brush.

### NOTE:

In the blank condition (unsintered), the frame may not be steamed off or subjected to other moist environments (e.g. ultrasonic bath) under any circumstances.

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.
- ▷ Mount sintering bowl and sintering aids together.
- ▷ Place sintering bowl incl. the sintering aids into the furnace.

## Final sintering

For the sintering, use only the Ceramill Argotherm high-temperature furnace with protective gas connection (please also refer to Ceramill Argotherm Operating Instructions).

### NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

**NOTE:**

Faulty sintering results!

- ▷ To ensure proper sintering results, pay attention to the correct arrangement of the sintering aids. Please also refer to the Ceramill Argotherm Operating Instructions.

- ▷ Perform the final sintering of the frames in the Ceramill Argotherm.
  - Program duration: approx. 5 hours
- ▷ Upon completion of the program, remove the sintering aids using pliers.

**NOTE:**

Faulty sintering results!

- ▷ Normally, the sintered bridges have a silver colour. When the frame or certain areas show signs of oxidation, select the sintering result using the document for evaluation of the sintering results (see Ceramill M-Center homepage, protected download area).

## Post-processing

### Post-processing fully anatomical dentures

After the final sintering has taken place:

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250  $\mu\text{m}$  aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework/adapt the frame using a cross-toothed carbide cutter.

- ▷ For fully anatomical restorations, polish using commercially available polishing aids (rubber polishers, brushes, NEM polishing paste).

### Post-processing/preparation for the veneering

- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250  $\mu\text{m}$  aluminium oxide at a pressure of 3 - 4 bar.
- ▷ If required, rework the frame using a cross-toothed carbide cutter.
- ▷ Thoroughly blast off the frame with a fine-blasting unit using pure 110 - 250  $\mu\text{m}$  aluminium oxide at a pressure of 3 - 4 bar.

- ▷ Grasp the frame with artery forceps; do not touch with fingers.

- ▷ Steam off frame.
- ▷ Oxidation firing is not necessary. Should an optical check nevertheless be desired, perform oxidation firing at 980 °C with 1 minute dwell time. Afterwards, thoroughly blast off one more time.

### Finish

- ▷ Gloss up metal edges using silicone polish. Prevent heat generation.

### Laser welding

- ▷ Prepare connecting location in the form of a modified X seam.
- ▷ Blast off and clean the spot to be laser welded (using 110 - 250  $\mu\text{m}$  pure aluminium oxide).
- ▷ Use CoCr laser welding wire (721130 ff).

After laser welding the Ceramill Sintron, the workpiece can be ceramic veneered.

## Soldering

- ▷ Solder using CoCrMo filler (781630 Giosolder) and the corresponding flux.

After soldering the Ceramill Sintron, the workpiece can be ceramic veneered.

## Veneer ceramics

As a general rule, all conventional metal veneer ceramics can be used.

- ▷ Observe and adhere the respective manufacturer's information!
- ▷ Match the thermal expansion coefficient (CTE) of the ceramic to the CTE of the alloy.

When using Creation ceramic material, do not apply bonder:

- ▷ Use Pulveropaker opaque powder.

## Fixation

### Conventional cementation

Due to the high strength and stability of Ceramill Sintronframes, conventional fixation with commercially available cement is possible in most cases.

- ▷ When applying conventional cementation, observe sufficient retention and a minimum abutment height of 3 mm!

## Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40  $\mu\text{m}$ ). For polishing, we recommend using diamond polishing bits.

**NOTE:**

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

## Follow-up examination at the dentist's

It is urgently recommended to check the monolithic dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

## Environmental Protection

### Packaging

In terms of packaging, Amann Girrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Disposal of the blanks

May not be disposed of together with household waste.

Do not dispose of into sewer systems.

Always dispose of according to official regulations.

## Accessories

For more information on product-specific accessories, please refer to [www.amanngirrbach.com](http://www.amanngirrbach.com).

## Download information

Further instructions/manuals and video tutorials can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com) and/or [www.amanngirrbach.com](http://www.amanngirrbach.com).

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**Amann Girrbach AG**
**Amann Girrbach GmbH**