Ceramill Sintron (ZZ)

for the Zirkonzahn® CAD/CAM system

Processing tips

Ceramill Sintron blanks comprise a CoCrMo alloy in the form of bounded powder metal. The material is dry milled in the green body state without the use of coolant lubricant (CL) and then densely sintered only in a special high-temperature sintering furnace (Ceramill Argotherm + Argovent sinter chamber) using the pre-programmed sintering programme with argon supply. We recommend observing the following points when processing Ceramill Sintron:

CAD DESIGN SOFTWARE

The fit of frameworks can be influenced by the adjustable parameters of the cement gap and the XY extension in the design software. Material-specific standard values are generally stored in your system, which can be varied according to the milling machine. Firstly, check using test crowns which values are required for your machine to produce the desired fit.

CAM SOFTWARE

Ceramill Sintron ZZ can be set up as a sinter metal blank. The shrinkage factor S , which can be found on the packaging label and product label, should be used for the Zirkonzahn milling machine.

MILLING MACHINE

Standard zirconia cutters can be used for milling Ceramill Sintron. It is recommended only to use these for Sintron. More detailed instructions regarding further processing can be found on our website: www.amanngirrbach.com/instruction-manuals

Note:

A powerful extractor is required when milling Ceramill Sintron. Centralised extraction is not recommended due to the often inadequate suction power!

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Amann Girrbach AG
Amann Girrbach GmbH