- DE** Gebrauchsanweisung
- EN** Instructions for use
- FR** Mode d'emploi
- IT** Istruzioni per l'uso
- ES** Instrucciones de uso
- PT** Instruções de utilização
- PL** Instrukcja użycia
- NL** Gebruiksaanwijzing
- SV** Bruksanvisning
- DA** Brugsanvisningen
- LT** Informacinis Lapelis
- LV** Lietošanas informācija
- HU** Használati tájékoztató
- CS** Informace k použití
- EL** Οδηγίες χρήσης
- RO** Instrucțiuni de utilizare
- USA** Instructions for use

Manufacturer | Hersteller

Distribution | Vertrieb

**Merz Dental GmbH**  
Eetzweg 20  
24321 Lütjenburg | Germany  
Fon +49 4381 403-0  
Fax +49 4381 403-403

[www.merz-dental.de](http://www.merz-dental.de)

907168/2018-10/04

**Amann Girrbach AG**  
Herrschaftswiesen 1  
6842 Koblach | Austria  
Fon +43 5523 62333-105  
Fax +43 5523 62333-5119  
[austria@amanngirrbach.com](mailto:austria@amanngirrbach.com)  
[www.amanngirrbach.com](http://www.amanngirrbach.com)

Made in the European Union

**CE** 0482

Rx Only

EN ISO 13485

32528+FB 2018-10-01

**AMANNGIRRBACH**

# Ceramill TEMP

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Hersteller

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## Ceramill TEMP

PMMA milling blank without fillers and fibres

### Instructions for use, please read carefully!

#### INTENDED USE

Ceramill TEMP, tooth colour, for functional restoration by replacement of lost tooth substance in permanent dentition in the form of temporary crowns and bridges in the anterior and posterior regions as well as for planning aesthetic corrections with the aid of temporary non-prep test veneers in the anterior region.

#### INDICATION

\_ Loss of hard tooth substance and resulting aesthetic impairment

#### CONTRAINDICATION

\_ Hypersensitivity to acrylate compounds

#### TECHNICAL DATA

PHYSICAL VALUES FOR CERAMILL SPLINTEC

##### Test according to

Solubility	0,2 µg / mm <sup>3</sup>	EN ISO 10477
Flexural strength	91,5 MPa	EN ISO 10477
Water uptake	26,5 µg / mm <sup>3</sup>	EN ISO 10477
Modulus of elasticity	2773 MPa	EN ISO 20795-1:2013
Flexural strength	96,6 MPa	EN ISO 20795-1:2013

#### COMPOSITION

- \_ Polymethylmethacrylate (PMMA) and cross-linkes polymers based on methacrylic acid esters
- \_ Colourants
- \_ Residual peroxide as dibenzoyl peroxide
- \_ Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%

#### SHADES

Available in 4 different dentine shades (bleach, light, middle, dark).

Blank	Tooth shade
bleach	BL1
light	A1; B1; D2; C1
middle	A2; A3; B2; B3; B4; C2; C3; D3; D4
dark	A3,5; A4; C4

#### WEARING PERIOD

Can be worn up to 3 years.

#### STORAGE

No special measures required.

#### WARNINGS

Processing Ceramill TEMP frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment. Further safety information can be obtained from the safety data sheet.

#### GENERAL INFORMATION ON CERMILL TEMP CONSTRUCTIONS, TOOTH COLOURED

\_ Minimum wall thickness after polishing ≥ 1 mm

#### SPECIAL INFORMATION ON CERMILL TEMP CONSTRUCTIONS, TOOTH COLOURED

- \_ In the construction of bridges with one pontic, the connection cross-section between the crown and pontic must not be less than 9 mm<sup>2</sup>.
- \_ In the construction of bridges spanning up to two pontics, the connection cross-section between crowns and pontics must not be less than 12 mm<sup>2</sup>.
- \_ By selecting the respective fixation system according to the indication, the user has to decide case-by-case whether the OMP-N® (PMMA) material ensures sufficient chewing stability for the restoration.

## TRIMMING

- \_ Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework and separate frameworks, crowns, bridges, denture bases, bite splints and drilling templates from the blank.
- \_ To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

## RESTORATION

Frameworks, crowns, bridges, denture bases, bite splints and drilling templates, which have been fabricated from Ceramill Splintec blanks can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Combipress N/LM, Weropress® or artConcept®. Prior conditioning of the surface is necessary, preferably with artConnect or a similar MMA-based bonding agent.

## FIXATION IN THE MOUTH

- \_ For temporary restorations: temporary cements
- \_ For long-term temporary restorations: glass ionomer cements without bonder

## NOTES

Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery. Enter the LOT No. with each procedure that requires identification of the material. To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

## DISPOSAL

Residues from processing the Ceramill Splintec have to be disposed of as residual waste.

## SAFETY DATA SHEETS

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from [www.merz-dental.de](http://www.merz-dental.de).

Date of information 2018-07

**The product properties are based on complying with and observing these instructions for use.**

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#### CAUTION

Federal (USA) law restricts this device by or on the order of a licensed dental practitioner.

#### NOTE

Ceramill Splntec is designed for milling equipment from Amann Girrbach or Equivalent.

#### INTENDED USE

Ceramill TEMP ,tooth coloured, indicated for use as a milling blank in the fabrication of provisional crowns and bridges.

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- \_ In the construction of bridges with one pontic, the connection cross-section between the crown and pontic must not be less than 9 mm<sup>2</sup>.
- \_ In the construction of bridges spanning up to two pontics, the connection cross-section between crowns and pontics must not be less than 12 mm<sup>2</sup>.
- \_ By selecting the respective fixation system according to the indication, the user has to decide case-by-case whether the OMP-N<sup>®</sup> (PMMA) material ensures sufficient chewing stability for the restoration.

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#### RESTORATION

Crowns, bridges, denture bases, bite splints and drilling templates, which have been fabricated from Ceramill Splntec can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Weropress<sup>®</sup>. Prior conditioning of the surface is necessary, preferably with a standard MMA-based bonding agent.

#### FIXATION IN THE MOUTH

- \_ For the temporary restoration: temporary cements
- \_ For long-term temporary restorations: glass ionomer cements without bonding

## Wear period

Can be worn up to 3 years

## STORAGE

No special measures required. See label for expiration date.

## NOTES

Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery. Record the LOT number with every process that requires identification of the material. To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

## RISK PHRASES

Processing Ceramill Splintec may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual postprocessing and that you are using personal protective equipment.

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 LOT	Batch code
 REF	Catalog number
	Manufacturer
	Follow the instructions for use
	Usable up to