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Dos and Don'ts of material-appropriate processing of glass-ceramics

Glasses and glass-ceramic are brittle materials. Special attention must be paid to material-appropriate design characteristics, positioning in the blank and the type of processing (wet grinding) when handling such brittle materials. The following information is an aid for processing this type of material.

CERAMILL MIND

_ Correct material and indication selection in the database

It is important to grind the defined indication with exactly the same material, which was selected in the database. The selection of the indication and material may not be changed later. This is the only way the correct grinding strategy will be used with the correct cutter radius adjustment.

_ In the design of anterior crowns with incisal edges that taper to point a 2-3° block-out angle for undercuts is recommended when creating the crown base.


_ Forms of design that are currently still critical

- x Inlays with recesses in the cavity can result in increased cutter wear or even cutter fracture.
- x Steep cusp slopes should be avoided. The cusps should be designed with a gentle slope.

_ Design of the fissure

- ✓ A rounded design (v-form) of the fissures is recommend. If needed use smoothing function.
- x Do not design a parallel and deep fissure. Deepening the fissure with the knife-point-function is not recommended.

CERAMILL MATCH 2

_ Selecting the correct positioning in the blank

- ✓ If possible place the front end of the design at the edge of the blank. This reduces the amount of grinding and therefore the cutter wear. Attach the pin orally or buccally.
- x Do not attach the pin mesially or distally; do not position the crown near to the adapter.

- ✓ Align veneers upright! The path of insertion selected in the CAD also defines the alignment of the design in the CAM.
- x Do not align veneers flat!

Dos and Don'ts of material-appropriate processing of glass-ceramics

- ✓ The occlusal aspect of the design must face upwards (the screw side at the holder).
- ✗ The cavity aspect of the design should not face upwards.


- ✓ Position the pin in the centre! (The midpoint of the pin should be level with the anatomical equator).
- ✗ Do not place the pin near the crown margin.

_ Inclining designs

Always incline the path of insertion of the design away from the holder - basically, from both the tool holder and the slot carrier.


_ Selecting the correct block

Two sizes of block are available:

1. Description I12 - height 10 mm
2. Description C14/I14 - height 12 mm

It is essential to ensure that the selected block size is clamped in the Ceramill Motion 2.


_ The tick for the Ceramill Coolstream must be activated during calculation of the milling path

Calculation and grinding should only be completed with cooling.

_ Placement of the Connector

The connector needs to be positioned exactly in the middle of the block, so that it is scrubbed from both sides equally to the highest point of the connector.

12-COMPARTMENT MULTIBLOCK HOLDER


_ Inclining designs to the different slot positions

Please note that it is not possible to incline a design on the first slot of each row (11, 21, 31). This is for safety reasons to prevent the spindle colliding with the tool holder. If a design is inclined on these slots, the inclination is automatically reset when the tool path is to be created or the connector position refreshed.

Dos and Don'ts of material-appropriate processing of glass-ceramics

If a design must be inclined, this is possible in slot 2 - 4 in all three rows.

CERAMILL COOLSTREAM


_ Coolant water supply

It is essential to ensure that sufficient coolant water is delivered to the instrument from the coolant nozzles. To check this, simply fit the wet operation cup and, with the safety guard (door) closed, keep the right key (S9) pressed until the coolant comes through (Important: no programmes should be active during the check). If insufficient coolant spray is delivered from the nozzles to the tip of the instrument, please check the following points.

1. Remove cover cap, clean (e.g. steam clean) and

re-attach.

2. Remove the water tank, check the filter and clean, if necessary. Also check whether grinding sludge has deposited at the bottom of the tank. If it has, drain the coolant and clean the entire tank. Then renew the coolant (mixing ratio 1:20) and replace the water tank in the Ceramill Coolstream.
3. Clean the coolant circulation. The cycle is restarted via the service and control function of the Motion software (the wet processing cup must also be fitted for this).
4. If the instrument has still not been rinsed sufficiently despite all the measures taken, please contact the helpdesk.

GENERAL

It is essential to follow the material manufacturer's instructions for use!