- DE** Gebrauchsanweisung
- EN** Instruction Manual
- FR** Instructions d'utilisation
- IT** Istruzioni per l'uso
- ES** Instrucciones de uso

Dieses Produkt enthält Nickel.
This product contains nickel.
Ce produit contient du nickel.
Questo prodotto contiene nichel.
Este producto contiene níquel.

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
www.amanngirrbach.com

9781935-IN-2020-06-20

AMANNGIRRBACH

DE - Original Gebrauchsanweisung -

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EN

- Translation of the original Operating Instructions -

Dental melt-on alloy based on NiCrMo, type 3 according to ISO 22674 for permanent dentures

For crowns and bridges

EXPLANATION OF SYMBOLS

Warning indications

- Warning indications in the text are marked with a triangle and boxed.
- In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- NOTE** means that property damage can occur.
- CAUTION** means that light to fairly serious personal injury can occur.
- WARNING** means that serious personal injury can occur.
- DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
>	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
MD	Medical device
UDI	Unique Device Identifier
REF	Catalog number
LOT	Batch code
	Manufacturer
	Consult instructions for use
	Date of manufacture

GENERAL SAFETY INSTRUCTIONS

CAUTION:
When processing constructions made of Girobond CBS, dusts can develop that may lead to mechanical irritation of the eyes and respiratory tract.
> Always make sure that the extraction system for the milling machine and the workplace for individual reworking operates properly.
> When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
> Further safety-relevant information can be found in the safety data sheet.

NOTE:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUITABLE PERSONNEL

NOTE:

The product may only be processed by trained dental technicians/dentists.

PROPERTIES

Contraindication

> In the case of known intolerances with the ingredients.

Side effects

Allergies against the alloy components or electrochemical-based paresthesia is possible in rare cases.
> Do not use Girobond CBS for patients with known nickel allergy.

Safety data sheet/Declaration of conformity

The safety data sheet is available for download on the Amann Girrbach website under *Services > Downloads > Additional Documents*. The product's declaration of conformity is available upon request from the manufacturer.

Technical Data

(ISO 22674 and ISO 9693)

	Unit
Proof stress (Rp 0.2)	335 MPa
Tensile strength (Rm)	520 MPa
E-module (E)	180 GPa
Ductile yield (A)	45 %
Vickers hardness	175 HV10
Melting range	1350 °C (solidus) - 1370 °C (liquidus)
Castling temperature	1410 °C
Density	8,4 g/cm ³
Thermal expansion coefficient	
– 25 - 500 °C	13,8 x 10 ⁻⁶ K ⁻¹
– 25 - 600 °C	14,0 x 10 ⁻⁶ K ⁻¹
Colour	Silver
Oxidation colour	Grey
Body-compatibility	Corrosion resistance according to ISO 10271 Biocompatibility according to ISO 10993-1 Tarnish resistance according to ISO 22674

Chemical composition in percent by mass

Nickel	Ni	63,5
Chromium	Cr	24,0
Molybdenum	Mo	10,0
Silicon	Si	1,5
further elements	Nb, Mn	< 1,0

Shelf life

Shelf life 5 years from date of manufacture.

Storage

> Store at room temperature in a dry location.

WARRANTY/EXCLUSION OF LIABILITY

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are based on our own experience and testing, and can therefore only apply as a guideline. Our products are subject to continuous further development. Subject to changes in design and composition.

PRODUCTION OF A CROWN/BRIDGE FOR CERAMIC VENEERING

Modelling

> Cap dipping; wall thickness 0.3 - 0.5 mm.
> Aim for an even ceramic layer of max. 2 mm above the modelling.

Castling-channel (sprue) arrangement

> Use cross beam or ring channel with 5 mm thickness.
> Fasten sprue channel with Ø 3 mm x 2 mm length at the highest point. Form soft transitions.

Investment

> Line out investment ring. For large investment rings (size 6 and 9), use two layers.
> Position modelling (wax-up) centred in the investment ring: Place cross beam in the heat centre with even distance to rim of investment ring.
> Process investment material according to manufacturer's specifications. Recommended investment material: Giroinvest Speed and Giroinvest Super.

Preheating

> Expel wax.
> Select heating rate according to specifications for processing of the investment materials.
> Select the dwell time at final temperature depending on the size of the investment ring. Select final temperature according to manufacturer's specifications (950 - 1050 °C).

Castling

(Castling temperature approx. 1420 °C, melting interval 1270 - 1356 °C)

When the melting pots are preheated in the furnace, the melting process can be shortened.
> Do not use melting pots for other alloys.
> Do not overheat alloy.
> Only cast new metal to ensure alloy properties and for clear batch traceability.

> Place Girobond CBS casting cylinder in a ceramic melting pot (1 cube - 6 g).

For open flame casting:

- Adjust fuel gas/oxygen ratio:
 - Fuel gas/oxygen ratio: Propane-oxygen ~0.7:1.4 bar
 - acetylene-oxygen ~0.7:0.7 bar
- Do not move flame away from the alloy during the melting process.
- When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: draw off immediately.

For high-frequency casting machines:

• When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: initiate casting process.

For vacuum pressure-casting machines:

> Melt-up alloy without vacuum.
Set the model casting program and switch off the vacuum or reduce it as far as possible.
> When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: initiate casting process.
> After the casting, allow the investment ring to slowly cool in the air.

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