9178000-NIN 2020-06-20

CE 0123
ISO 13485
ISO 9001

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girschbach AG
Herrschaftswiesen 1
6842 Koblach | Austria
Tel +43 5523 62333-105
Fax +43 5523 62333-519
austria@amanngirschbach.com
www.amanngirschbach.com

AMANNGIRSBACH

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Dental melt-on alloy based on CoCrM, type 4 according to ISO 22674 for permanent and removable dentures

For crowns, bridges and model castings

EXPLANATION OF SYMBOLS

Warning indications
Warning indications in the text are marked with a triangle and boxed.
In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.
- **NOTE** means that property damage can occur.
- **CAUTION** means that light to fairly serious personal injury can occur.
- **WARNING** means that serious personal injury can occur.
- **DANGER** means that life-threatening personal injury can occur.

Important information
Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
	Item of an operation description
	Item of a list
	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
	Medical device
	Unique Device Identifier
	Catalog number
	Batch code
	Manufacturer
	Consult instructions for use
	Date of manufacture

GENERAL SAFETY INSTRUCTIONS

CAUTION:
When processing constructions made of Girobond NB, dusts can develop that may lead to mechanical irritation of the eyes and respiratory tract.
- Always make sure that the extraction system for the milling machine and the workplace for individual reworking operates properly.
- When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- Further safety-relevant information can be found in the safety data sheet.

NOTE:
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUITABLE PERSONNEL

NOTE:
The product may only be processed by trained dental technicians/dentists.

PROPERTIES

Contraindication
In the case of known intolerances with the ingredients.

Side effects
Allergies against the alloy components or electrochemical-based paresthesia is possible in rare cases.

Safety data sheet/Declaration of conformity
The safety data sheet is available for download on the Amann Girschbach website under *Services > Downloads > Additional documents*. The product's declaration of conformity is available upon request from the manufacturer.

Technical Data
(ISO 22674 and ISO 9693)

	Unit
Proof stress (Rp 0,2)	530 MPa
Tensile strength (Rm)	760 MPa
E-module (E)	210 GPa
Ductile yield (A)	10 %
Vickers hardness	310 HV10
Melting range	1370 °C (solidus) - 1400 °C (liquidus)
Casting temperature	1510 °C
Density	8,5 g/cm ³
Thermal expansion coefficient (CTE) (25 - 500°C)	14,6 × 10 ⁻⁶ K ⁻¹
Colour	Silver
Oxidation colour	Grey-green
Body-compatibility	Corrosion resistance according to ISO 10271 Biocompatibility according to ISO 10993-1 Tarnish resistance according to ISO 22674

Chemical composition in percent by mass

Cobalt	Co	62,0
Chromium	Cr	25,0
Molybdenum	Mo	5,0
Tungsten	W	5,0
Silicon	Si	1,2
further elements	Ce, Fe, Nb, N	<1,0

This product is "nickel-free" which means that it contains less than or equal to 0.1 % (mass percentage) of Ni as defined in ISO 22674.

Shelf life
Shelf life 5 years from date of manufacture.

Storage
- Store at room temperature in a dry location.

WARRANTY/EXCLUSION OF LIABILITY

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are based on our own experience and testing, and can therefore only apply as a guideline. Our products are subject to continuous further development. Subject to changes in design and composition.

PRODUCTION OF A CROWN/BRIDGE FOR CERAMIC VENEERING

Modelling
- Cap dipping; wall thickness 0.3 - 0.5 mm.
- Aim for an even ceramic layer of max. 2 mm above the modelling.

Casting-channel (sprue) arrangement
- Use cross beam or ring channel with 5 mm thickness.
- Fasten sprue channel with Ø 3 mm x 2 mm length at the highest point. Form soft transitions!

Investment
- Line out investment ring. For large investment rings (size 6 and 9), use two layers.
- Position modelling (wax-up) centred in the investment ring; Place cross beam in the heat centre with even distance to rim of investment ring.
- Process investment material according to manufacturer's specifications. Recommended investment material - Giroinvest Super.

Preheating
- Expel wax.
- Select heating rate according to specifications for processing of the investment materials.
- Select the dwell time at final temperature depending on the size of the investment ring. Select final temperature according to manufacturer's specifications (950 - 1050 °C).

Casting
(Casting temperature approx. 1510 °C, melting interval 1370 - 1400 °C)

When the melting pots are preheated in the furnace, the melting process can be shortened.
- Do not use melting pots for other alloys.
- Do not overheat alloy.
- Only cast new metal to ensure alloy properties and for clear batch traceability.

- Place Girobond NB casting cubes in a ceramic melting pot (1 cube ~ 6 g).

For open flame casting:
- Adjust fuel gas/oxygen ratio:
Propane-oxygen ~0,7 : 1,4 bar
Acetylene-oxygen ~0,7 : 0,7 bar
- Do not move flame away from the alloy during the melting process.
- When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: draw off immediately.

For high-frequency casting machines:
- When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: initiate casting process.

For vacuum pressure-casting machines:
- Melt-up alloy without vacuum.
- Set the model casting program and switch off the vacuum or reduce it as far as possible.
- When the last casting cube has joined the molten material and the molten material tears apart and the surface is glossy: initiate casting process.
- After the casting, allow the investment ring to slowly cool in the air.

Divesting / Finishing
- Blast off with 110 - 250 µm aluminium oxide.
- Cut off the casting channels with a cutting disc.
- Rework the frame using a cross-toothed carbide cutter.

Metal treatment before ceramic veneering
- Blast off the frame with a fine-blasting unit using pure 110 - 250 µm aluminium oxide at a pressure of 3 - 4 bar.
- Grasp the frame with artery forceps; do not touch with fingers.
- Steam off frame with steam cleaner.
- If necessary, carry out oxide firing for surface blast at 1040 °C with 1 - 2 min holding time. Afterwards, thoroughly inspect off again (see above).

NOTE:
Possible frame damage through improper handling!
- Do **not** blow frame dry with compressed air (oil, debris).
- Do **not** place frame into hydrofluoric acid.

Ceramic veneering
- Observe manufacturer's specifications for the ceramic being used.

Finish
- Gloss up metal edges using silicone polish. Prevent heat generation.

PRODUCTION OF MODEL CASTING WORK

Spruing technique
- Provide for casting channels to the most massive modelling areas.

Investment
- Process phosphate-bonded model-cast investment material according to manufacturer's specifications.

Preheating
- Expel wax.
- Select heating rate according to specifications for processing of the investment materials.
- Select the dwell time at final temperature depending on the size of the investment ring. The final temperature is 980 - 1050 °C.

Casting
When the melting pots are preheated in the furnace, the melting process can be shortened.
- Do not use melting pots for other alloys.
- Do not overheat alloy.
- Only cast new metal to ensure alloy properties and for clear batch traceability.

- Place casting cylinder in a ceramic melting pot (1 cube ~ 6 g).
- When the last casting cylinder has joined the molten material and the molten material tears apart, draw off immediately.
- After the casting, allow the investment ring to slowly cool in the air.

Divesting / Finishing
- Blast off frame with 110 - 250 µm aluminium oxide.
- Blast off critical, delicate areas with 50 µm aluminium oxide.
- Cut off the casting channels with a cutting disc.
- Work out frame using carbide bur (cross-toothed) or ceramic-bonded stone.
- Polish frame.

GENERAL

Laser welding
- Blast off and clean the spot to be welded (using 110 - 250 µm aluminium oxide).
- Prepare connecting location in the form of a modified X seam.
- Use CoCr laser welding wire (72H30 Hf) as core-up material.

After welding, the workpiece can be ceramic veneered.

Soldering
- Before veneering: solder with CoCrMo solder and corresponding flux.
- After veneering: in oven (at 860 °C) solder with white gold solder and corresponding flux.

Fixation
Due to the high strength and stability of the frames, conventional fixation with commercially available cement is possible in most cases.

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