giroinvest speed

Gebrauchsanweisung
User Manual
Mode d'emploi
Istruzioni d'uso
Modo de empleo

AMANN GIRRBACH

Deutsch	04 - 06
English	07 - 09
Française	10 - 12
Italiano	13 - 15
Español	16 - 18

GIROINVEST SPEED

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GIROINVEST SPEED

Giroinvest Speed is a phosphate-bonded precision investment for model-casting restorations made with the speed technique. It also suits the linear preheating as well as the conventional preheating with delays.

Giroinvest Speed can either be used with casting ring or ringless.

Giroinvest Speed:

Order no.:

724070 Pack contents 5x 4 kg

724072 Pack contents 100x 200 g bag
(= 20 kg)

Giroinvest Liquid:

Order no.:

724081 Pack contents 1l

Technical Data:

Setting expansion 1,7 %

Thermal expansion 1,3 %

Overall expansion 3,0 %

Figures relate to 100 % Giroinvest Speed Liquid

Mixing ratio

100 g Giroinvest Speed: 18 ml liquid

Working time 4 min

Mixing:

15 s manual premix

45 s evacuate without mixing

Setting time 20 min

Times relate to 20 °C working temperature

MIXING

Mixing ratio:

100 g Giroinvest Speed 18 ml liquid

200 g Giroinvest Speed/1 bag 36 ml liquid

_ Always use distilled water for diluting Giroinvest Liquid mixing concentrate.

Liquid concentration:	Ratio of Giroinvest Liquid	Ratio of distilled water
Model-casting restoration with cast clasps CoCrMo alloy	70 - 75 %	30 - 25 %
combination technique CoCrMo alloy	80 - 85 %	20 - 15 %

All ratios relate to 100g Giroinvest Speed: 18ml liquid at 20 °C working temperature. The ratios are approximate. Different techniques, accessories and equipment can alter the results.

Working time:

4 min (working temperature 20 °C)

- _ Rinse a clean mixing bowl with water and wipe it before filling with investment. Dirty or dry mixing bowls reduce the liquid content of the investment.
- _ Pour in the liquid and then sprinkle in the powder.
- _ Manually mix the powder and liquid thoroughly with a spatula for 15 s
- _ Use a vacuum mixer for further mixing. Mix under vacuum for 45 s

Setting time:

20 min (working temperature 20 °C)

SPEED CAST TECHNIQUE

- _ Adhere strictly to the setting time and temperature for placing the casting ring in the furnace with the speed cast technique.

	Preheating speed cast
Insertion temperature	850 °C
Final temperature max.	1050 °C
Hold time at final temperature	app. 60 min depending on size and number of casting rings

- _ Adding casting rings later reduces the furnace temperature and increases the preheat time!
- _ Speed cast technique: Place all casting rings quickly into the furnace and do not open the furnace for 15 min!

CONVENTIONAL PREHEATING WITH DELAYS

	1. delay	2. delay	3. delay
Heat rate	2 °C/min	4 °C/min	5 °C/min
Temperature	270 °C	580 °C	depending on the alloy (however max. 1050 °C)
Hold time	app. 30 min depending on size and number of casting rings	app. 30 min depending on size and number of casting rings	app. 30 min depending on size and number of casting rings

Linear preheating:

_ linear temperature increase from room temperature to final temperature (without delays) at a heat rate of 5-7 °C/min.

After casting:

- _ Allow the casting ring to bench cool to hand heat after casting.
- _ Do not quench in water!

!	<p>WARNINGS AND SAFETY INSTRUCTIONS</p> <ul style="list-style-type: none"> _ The investment contains quartz and cristobalite. Do not inhale dust, as there is a risk of damage to the lungs (silicosis, lung cancer)! <p>Recommendations:</p> <ul style="list-style-type: none"> _ Use a protective face mask, type FFP 2-EN 149:2001. _ Open the foil bag with scissors and avoid creating dust when filling the mixing bowl. _ Rinse the empty foil bag with water before disposal. <ul style="list-style-type: none"> _ Always use a damp cloth to remove investment dust or dried mixing liquid from the working area. _ Soak the cooled casting ring in water before devesting the casting. _ Do not open the preheat furnace during the heating stage, as the wax vapours given off can ignite with air! Furnace gases released should be extracted and directed outside. _ Speed cast technique: Place all casting rings quickly into the furnace and do not open the furnace for 15 min!
----------	--

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Amann Girrbach AG
Amann Girrbach GmbH