# starterkit frästechnik

*Milling | Fraisage*

*Fresatura | Técnica de fresado*

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## STARTERKIT FOR MILLING TECHNIQUE

- 1 Transfer support
- 2 Undercut gauge
- 3 Milling oil
- 4 Training cylinder
- 5 Training cylinder holder
- 6 Plaster plate
- 7 Red and blue carbon markers
- 8 Carbon marker holder
- 9 Milling set 12-part
- 10 Conical cutter, cross-cut, round tip
- 11 Parallel cutter, cross-cut, round tip
- 12 Conical cutter, non-cross-cut, round tip
- 13 Parallel cutter, non-cross-cut, round tip
- 14 Twist drill, thin
- 15 Twist drill, thick
- 16 Shoulder cutter
- 17 Centring drill
- 18 Wax cutter, conical
- 19 Wax cutter, parallel
- 20 Channel cutter
- 21 Wax carver

### 1 TRANSFER SUPPORT

The device is for transferring milling work from the master model to the plaster table or the plaster milling support. It is available with a 2,35 mm or 3,0 mm shaft.

### 2 UNDERCUT GAUGE

This set contains 3 undercut gauges of different sizes for indicating undercuts on teeth. All gauges have a 3,00 mm diameter shank.

### 3 MILLING OIL

The red milling oil in this set comprises a combination of high grade lubricant and coolant and was developed specifically for use with modern milling techniques.

### 4 TRAINING CYLINDER

The two brass training cylinders are for learning how to operate the new milling machine. They are also perfect for teaching beginners how to operate milling machines and rotary instruments. The type 1 training cylinder includes a symmetric occlusal surface for simulating various types of milling procedure. The type 2 training cylinder can also be used for teaching drilling techniques and milling brace supports etc.

### 5 TRAINING CYLINDER HOLDER

The training cylinder holder is for gripping the training cylinder (described above). A locking screw grips the training cylinder in the holder firmly. The training cylinder holder is then clamped on the model platform and the training session can begin.

### 6 PLASTER PLATE

The plaster plate is used for producing custom milling bases which provide easy access to the restoration being milled. This allows the restoration to be milled perfectly without damaging the alveolar ridge or model.

### 7 RED AND BLUE CARBON MARKERS

Carbon markers are used for marking the line of maximum bulbosity on the tooth. The set includes red and blue markers to improve the contrast when marking the design of the restoration.

### 8 CARBON MARKER HOLDER

The carbon marker holder grips the carbon marker in the head of the surveyor or milling machine. The marker is simply inserted into the hollow holder and gripped in the surveying head with the chuck.

## 9 MILLING SET 12-PART

### 10 CONICAL CUTTER, CROSS-CUT, ROUND TIP

Conical cutters are used for fabricating precious metal conical crowns. The cross-cut blades remove large amounts of material from the restoration. A speed of approximately 10.000 r.p.m. is recommended for fine-finish milling. Superfinishing, precision grinding and polishing should be carried out with the help of cotton wool and milling oil at a speed of approximately 3.000 r.p.m.

### 11 PARALLEL CUTTER, CROSS-CUT, ROUND TIP

Parallel cutters are used when fabricating precious metal channelshoulder and channel-shoulder-pin precision attachments. The crosscut blades reduce large amounts of material. The round tip of the cutter produces a radius shoulder on the restoration. For coarse and fine-finish milling it is advisable to use milling oil and run the cutter at a speed of approximately 10.000 r.p.m. Superfinishing, precision grinding and polishing should be carried out with the help of cotton wool and milling oil at a speed of approximately 3.000 r.p.m.

### 12 CONICAL CUTTER, NON-CROSS-CUT, ROUND TIP

(Refer to section 10) Non-cross-cut conical cutters are mostly used for dressing, i.e. superfinishing milled surfaces.

### 13 PARALLEL CUTTER, NON-CROSS-CUT, ROUND TIP

(Refer to section 12) Non-cross-cut parallel cutters are mostly used for dressing, i.e. superfinishing milled surfaces.

## 14 TWIST DRILL, THIN

Twist drills are used for drilling holes when fabricating channelshoulder and channel-shoulder-pin precision attachments as well as Interlock or "T" attachments. The hole should be drilled at approximately 10.000 r.p.m., using milling oil and exerting only minimal pressure - the chips must be removed frequently.

**CAUTION!** The drill may break if excessive pressure is exerted.

## 15 TWIST DRILL, THICK

The hole should be drilled at approximately 10.000 r.p.m., using milling oil and exerting only minimal pressure - the chips must be removed frequently.

**CAUTION!** The drill may break if excessive pressure is exerted.

## 16 SHOULDER CUTTER

Shoulder cutters are used for counterboring the concentric shoulder of an Interlock attachment. Forward feed is achieved with the milling spindle. When fabricating channel-shoulder and channel-shoulder-pin precision attachments, the shoulder is milled on the attachment. In both cases, the cutter should be run at a speed of approximately 3.000 r.p.m. and used with milling oil.

## 17 CENTRING DRILL

Centring drills are used to mark the sites for drilling holes when fabricating channel-shoulder and channel-shoulder-pin precision attachments as well as Interlock or "T" attachments. The purchase point centres the twist drill. The drill should be run at approximately 5.000 r.p.m. and forward feed is achieved with the milling spindle.

## 18 WAX CUTTER, CONICAL

Wax cutters are used for trimming the precision attachment pattern to produce conical surfaces. The wax is cut-down milled from the pattern at a speed of approximately 3.000 r.p.m. Cut-down milling prevents wax chips sticking to the milled surface.

## 19 WAX CUTTER, PARALLEL

(Refer to section 18) Wax cutters are used for trimming the precision attachment pattern to produce parallel surfaces. The wax is cut-down milled from the pattern at a speed of approximately 3.000 r.p.m. Cut-down milling prevents wax chips sticking to the milled surface.

## 20 CHANNEL CUTTER

Channel cutters are used for cutting the channels in patterns for channel-shoulder and channel-shoulder-pin precision attachments. The channels should be milled at a speed of approximately 3.000 r.p.m. and separating agent applied to the wax to prevent it smearing. The milling spindle provides for forward feed.

## 21 WAX CARVER

Wax carvers are used for trimming patterns for milling. The wax should be coated with separating agent to prevent it smearing and then trimmed off the pattern.

### Accessories/Components:

Order No.:

873001	Starterkit for milling technique, 12 parts
177810	1 Transfer support 2,35 mm
177820	1 Transfer support 3,00 mm
177830	2 Undercut gauge 0,25 mm
177840	2 Undercut gauge 0,50 mm
177850	2 Undercut gauge 0,75 mm
177656	3 Milling oil
177653	4 Training cylinder (Typ 1)
177654	4 Training cylinder (Typ 2)
177655	5 Training cylinder holder
171930	6 Plaster plate
177880	7 Mine (red 12-part)
177870	7 Mine (blue 12-part)
177860	8 Carbon marker holder
873002	9 Milling set 12-part
873022	10 Conical cutter, cross-cut, round tip
873020	11 Parallel cutter, cross-cut, round tip
873021	12 Conical cutter, non-cross-cut, round tip
873019	13 Parallel cutter, non-cross-cut, round tip
873017	14 Twist drill, thin
873016	15 Twist drill, thick
873018	16 Shoulder cutter
873015	17 Centring drill
873013	18 Wax cutter, conical
873012	19 Wax cutter, parallel
873014	20 Channel cutter
873100	21 Wax carver

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**Amann Girrbach AG**
**Amann Girrbach GmbH**