



Long-Life Tools

3mm	Art. Nr. Ring- farbe	Beschreibung	Maschinen										Material																										
			Mikro 4x	Mikro 5x	Mikro IC	Motion 2	Motion 3/Motion 3 Air	Matik	Motion DRS	Matron/Matron Core	Motion Air	Zirkonoxid	Sinterlegierung	PMMA	TEMP	TEMP ML	Splintec	PEEK	WAX	Test PMMA	MOGU Kunststoff	Modell PU	Basis Wachs	Basis Kunststoff	Prothesen Zähne	Hybride Mill	Hybride Grind	e.max	nice	Glaskeramik	Premills	CoCr hart	Kalibrierefräsen	Zahnkranz	PMMA-Monoblock	Composite	Titan hart	Testmaterial	
	751021	Roto RFID D3 grey 2,5	•	•	•	•	•	•	•	•	•	•	•												•	•								•	•				
	751022	Roto RFID D3 grey 1,0	•	•	•	•	•	•	•	•	•	•	•												•	•								•	•				•
	751023	Roto RFID D3 grey 0,6	•	•	•	•	•	•	•	•	•	•	•												•	•								•	•				
	751024	Roto RFID D3 grey 0,3	•	•	•	•	•	•	•	•	•	•	•												•														
	751025	Roto RFID D3 green 2,5	•	•		•	•	•		•	•			•	•	•	•	•	•	•	•	•	•	•											•				
	751026	Roto RFID D3 green 1,5	•	•		•	•	•		•	•											•	•	•											•				
	751027	Roto RFID D3 green 1,0	•	•		•	•	•		•	•			•	•	•	•	•	•	•	•	•	•																
	751028	Roto RFID D3 green 0,6	•	•		•	•	•		•	•			•	•	•	•	•	•	•	•	•	•																
	751029	Roto RFID D3 green 0,3	•	•		•	•	•		•	•							•				•																	

6mm

	759600	Roto RFID Grey 2,5								•		•	•												•	•									•	•			•	
	759601	Roto RFID Grey 1,0								•		•	•												•	•									•	•			•	
	759602	Roto RFID Grey 0,6								•		•	•												•	•									•	•			•	
	759603	Roto RFID Grey 0,3								•		•	•												•	•									•	•			•	
	759604	Roto RFID Green 2,5								•				•	•	•	•	•	•	•	•	•	•	•											•					
	759605	Roto RFID Green 1,0								•				•	•	•	•	•	•	•	•	•	•	•																
	759606	Roto RFID Green 0,6								•				•	•	•	•	•	•	•	•	•	•																	
	759607	Roto RFID Green 0,3								•									•			•																		
	759627	Roto RFID Green 1,5								•												•	•	•											•					
	759615	Roto RFID Brown R 1,2 EM								•																							•				•			
	759616	Roto RFID Brown F 1,2 EM								•																								•			•			
	759619	Roto RFID Brown 3,0								•																								•			•			
	759611	Roto RFID Brown 2,5								•																														
	759612	Roto RFID Brown 2,0								•																														
	759617	Roto RFID Brown 1,5 TS								•																														
	759618	Roto RFID Brown 1,5 SM								•																														
	759613	Roto RFID Brown 1,5								•																														
	759614	Roto RFID Brown 1,0								•																														
	751033	Roto RFID brown 1,2 EM								•		•	•	•	•	•		•	•	•													•			•				
	751034	Roto RFID brown 0,5 EM								•		•	•	•	•	•		•	•	•														•			•			
	751035	Roto RFID brown 0,5 SM								•		•	•	•	•	•		•	•	•															•			•		
	759622	Roto RFID Purple 2,0 TS								•		•	•																											
	759624	Roto RFID Purple 1,5 TS								•		•	•		•			•	•																					
	759625	Roto RFID Purple 1,2 EM								•		•	•		•	•		•	•	•					•											•				
	759626	Roto RFID Purple 0,6 EM								•		•	•		•	•		•	•	•																				
	759628	Roto RFID purple 1,5 EM								•		•	•		•	•		•	•	•					•															
	759623	Roto RFID Purple 1,8 TS								•					•	•		•																						
	759621	Roto RFID Calibration EM								•																													•	
	759620	Roto RFID Calib. 6 mm								•																														
	759608	Roto RFID Orange G 2,0								•																	•	•	•	•										
	759609	Roto RFID Orange G 1,0								•																	•	•	•	•							•			
	759610	Roto RFID Orange G 0,4								•																	•	•	•	•							•			

Classic Tools

3mm

	760660	Roto RFID 2,5 ZI	•	•		•	•	•		•	•	•																									
	760661	Roto RFID 1,0 ZI	•	•		•	•	•		•	•	•																									•
	760662	Roto RFID 0,6 ZI	•	•		•	•	•		•	•	•																									
	760663	Roto RFID 0,3 ZI	•	•		•	•	•		•	•	•																									
	751010	Roto RFID 2,5 DC BLM							•	•	•		•													•											
	751011	Roto RFID 1,0 DC BLM							•	•	•		•													•											
	751012	Roto RFID 0,6 DC BLM							•	•	•		•													•											

3mm	Art. Nr. Ring- farbe		Beschreibung	Maschinen								Material																																			
				Mikro 4x	Mikro 5x	Mikro IC	Motion 2	Motion 3/Motion 3 Air	Matik	Motion DRS	Matron/Matron Core	Motion Air	Zirkonoxid	Sinterlegierung	PMMA	TEMP	TEMP ML	Splintec	PEEK	WAX	Test PMMA	MOGU Kunststoff	Modell PU	Basis Wachs	Basis Kunststoff	Prothesen Zähne	Hybride Mill	Hybride Grind	e.max	nice	Glaskeramik	Premills	CoCr hart	Kalibrierfräsen	Zahnkranz	PMMA-Monoblock	Composite	Titan hart	Testmaterial								
	760652	Roto WMB DC 2,5			•							•			•	•										•	•																				
	760653	Roto WMB DC 1,0			•							•			•	•										•	•																				
	760654	Roto WMB DC 0,6			•							•			•	•										•																					
	760655	Roto WMB DC 0,3			•							•			•	•										•																					
	760684	Roto RFID 2,5 DMB DC	•	•	•	•	•	•	•	•	•				•											•	•																				
	760685	Roto RFID 1,0 DMB DC	•	•	•	•	•	•	•	•	•				•											•	•																				
	760686	Roto RFID 0,6 DMB DC	•	•	•	•	•	•	•	•	•				•											•																					
	760687	Roto RFID 0,3 DMB DC	•	•	•	•	•	•	•	•	•															•																					
	760668	Roto RFID 2,5 PMMA/Wax	•	•		•	•	•		•	•				•	•	•	•	•	•	•	•																									
	760669	Roto RFID 1,0 PMMA/Wax	•	•		•	•	•		•	•				•	•	•	•	•	•	•	•		•																							
	760670	Roto RFID 0,6 PMMA/Wax	•	•		•	•	•		•	•				•	•	•	•	•	•	•	•																									
	760671	Roto RFID 0,3 PMMA/Wax	•	•		•	•	•		•	•								•					•																							
	760664	Roto RFID 2,5 Sintron	•	•		•	•	•		•	•		•																																		
	760665	Roto RFID 1,0 Sintron	•	•		•	•	•		•	•		•																																		
	760666	Roto RFID 0,6 Sintron	•	•		•	•	•		•	•		•																																		
	760667	Roto RFID 0,3 Sintron	•	•		•	•	•		•	•		•																																		
	760680	Roto RFID 2,5 SC				•	•		•	•					•	•	•		•		•			•																							
	760681	Roto RFID 1,0 SC				•	•		•	•					•	•	•		•		•			•																							
	760682	Roto RFID 0,6 SC				•	•		•	•					•	•	•		•		•			•																							
	760683	Roto RFID 0,3 SC				•	•		•	•									•																												
	751030	Roto RFID D3 brown 1,2 EM					•		•			•	•	•	•	•		•	•	•																						•	•				
	751031	Roto RFID D3 brown 0,5 EM		•		•	•	•		•	•		•	•	•	•		•	•	•						•																	•	•			
	751032	Roto RFID D3 brown 0,5 SM					•		•			•	•	•	•	•		•	•	•					•																						
	760688	Roto RFID 2,5 COCR							•	•																																					
	760691	Roto RFID 1,5 COCR Telescope					•		•																																						
	760689	Roto RFID 1,5 COCR					•		•																																						
	760690	Roto RFID 1,0 COCR					•		•																																						
	751014	Roto RFID 2,0 T-Shape		•		•	•	•		•	•		•	•																																	
	751003	Roto RFID 1,5 T-Shape		•	•		•	•	•		•	•		•	•			•	•		•	•																									
	751015	Roto RFID 1,8 T-Shape			•		•	•	•		•	•			•	•		•																													
	751001	Roto RFID 1,2 ST		•	•		•	•	•		•	•		•	•		•	•	•		•	•	•				•																				
	751036	Roto RFID D3 purple 1.5 ST		•	•		•	•	•		•	•		•	•		•	•	•		•	•	•				•																				
	751002	Roto RFID 0,6 ST		•	•		•	•	•		•	•		•	•		•	•	•		•	•	•																								
	751020	Roto RFID 2.5 Calibration								•	•																																			•	
	181205	Roto RFID Calib. 3 mm		•		•	•	•	•	•	•																																				
	760702	Roto RFID 3,0 Titan					•		•																																					•	
	760703	Roto RFID 2,0 Titan					•		•																																					•	
	760708	Roto RFID 1,5 T-Shape Titan					•		•																																					•	
	760704	Roto RFID 1,5 Titan					•		•																																					•	
	760707	Roto RFID 1,2 R Titan					•		•																																					•	
	760705	Roto RFID 1,2 F Titan					•		•																																					•	
	760706	Roto RFID 1,0 Titan					•		•																																					•	
	760672	Roto RFID 2,0 Ti			•	•	•	•	•	•																																					
	760673	Roto RFID 1,0 Ti			•	•	•	•	•	•																																					
	760676	Diamond RFID 1,8			•	•	•	•	•	•																		•	•	•	•	•															
	760677	Diamond RFID 1,4			•	•	•	•	•	•																		•	•	•	•	•															
	760678	Diamond RFID 1,0			•	•	•	•	•	•																																					