

Ceramill Matron: Serious Fun

Supercharge your workflows.

Ceramill Matron is the fifth generation milling machine from Amann Girrbach, with superior properties compared to machines of older generations. Its extraordinary product features make it stronger, faster and, smarter than conventional milling machines – offering maximum ease of use and simply being more fun to work with. Whether you are a CAD/CAM expert or newcomer, you will benefit from the outstanding precision and quality milling results across a wide range of materials, particularly in the field of carbide processing.



Learn more about true high performance metal milling with Ceramill Matron:
bit.ly/3S87u5x



The advantages at a glance

- ✓ High performance metal milling with 3mm or 6mm collet
- ✓ RFID enabled performance 8x blank changer
- ✓ 45° A/B axes and 5-axis+
- ✓ SnapMag: Magnetic Tool Tray
- ✓ Guided workflows on the 21.5" display



Stress Less

At Amann Girrbach, we believe successful dental workflows depend on more than product performance. Seamless digital integration between dental laboratories and practices is equally essential—especially where both work closely together.

We focus on providing safe, validated workflows with perfectly matched system components to ensure outstanding results. At the same time, our systems remain open to selected third-party integrations, giving users maximum flexibility.

This philosophy shaped the development of Ceramill Matron. Intelligent automation, simple operation, and an interactive user experience enable smooth, stress-free workflows—even under high daily workloads.



Ordering information

	ARTICLE NO.	PRODUCT DESIGNATION
MACHINE	181371	Ceramill Matron
ACCESSORIES	181376	Ceramill CoolantFlow
	181375	Ceramill Tank Matron
	181251	Collet wrench torque Syc.3/6
	181378	Ceramill Matron SnapMag
	181389	Auto-calibration kit
CART	181377	Ceramill Cart Matron
EXTRACTION	178800	Renfert Silent powerCAM EC
HOLDER	181380	Blank holder C-Clamp Matron
	181381	Blank holder 98 Matron
	181383	Blank holder blocks 4x UN Matron
	181384	Blank holder 4x Ti-Preforms Matron
	181385	Blank holder metal Matron



“As a laboratory that has worked extensively with cobalt chrome, we are very impressed by the smoothness of the Ceramill Matron while milling this extremely hard material. This results in long milling cutter service lives and hence a plus in the efficient use of resources.”

Andreas Zeller CDT & Moritz Zeller CDT,
Zeller Zahntechnik, Ulm, Germany

Stress less with Ceramill Matron.