



Ref. Certif. No.

DE 3 - 30808

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME

SYSTEME CEI D'ACCEPTATION MUTUELLE DE CERTIFICATS D'ESSAIS DES EQUIPEMENTS ELECTRIQUES (IECEE) METHODE OC

CB TEST CERTIFICATE CERTIFICAT D'ESSAI OC

Product

Laboratory Equipment
(dental scanner)

Name and address of the applicant

Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach
AUSTRIA

Name and address of the manufacturer

Amann Girrbach AG
Herrschaftswiesen 1, 6842 Koblach, AUSTRIA

Name and address of the factory

smart optics Sensortechnik GmbH
Lise-Meitner-Allee 10, 44801 Bochum, GERMANY

Ratings and principal characteristics

Rated voltage: 100-240V
Rated frequency: 50-60Hz
Rated power consumption: 60W
Protection class: I
Laser class: 1

Trade mark (if any)

Amann Girrbach

Model/type Ref.

Ceramill MAP600, IDC S3

Additional information (if necessary)

Certificate DE 3 - 30754 issued 2018-07-18 is replaced by this version due to technical changes

A sample of the product was tested and found to be in conformity with

IEC 61010-1:2010

as shown in the Test Report Ref. No. which forms part of this certificate

028-713134717-100

This CB Test Certificate is issued by the National Certification Body
Ce Certificat d'essai OC est établi par l'Organisme **National de Certification**

CB 019612 0061 Rev. 00

Date, 2018-10-10

(Abdul Sabbagh)



Product Service

Page 1 of 1

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany