

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT
(IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Laboratory Equipment

Name and address of the applicant

Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach
AUSTRIA



AMANNGIRRBACH

Name and address of the manufacturer

Amann Girrbach AG
Herrschaftswiesen 1, 6842 Koblach, AUSTRIA

Name and address of the factory

Amann Girrbach AG
Schweizerstrasse 96, 6830 Rankweil, AUSTRIA

Ratings and principal characteristics

| | | |
|---------------------------------|----------------|----------------------|
| Model reference: | Ceramill Matik | Ceramill Cleanstream |
| Rated voltage: | 100-230V | 100-230V |
| Rated frequency: | 50-60Hz | 50-60Hz |
| Rated power consumption: | 400W | 200W |
| Protection class: | I | I |

Trade mark (if any)

AMANNGIRRBACH

Customer's Testing Facility (CTF) Stage used

CTF STAGE 1

Model/type Ref.

Ceramill Matik, Ceramill Cleanstream

A sample of the product was tested and found to be in conformity with

IEC 61010-1:2010
IEC 61010-1:2010/AMD1:2016

as shown in the Test Report Ref. No. which forms part of this certificate

028-713165776-000

This CB Test Certificate is issued by the National Certification Body

CB 019612 0066 Rev. 00

Date, 2019-09-23

(Abdul Sabbagh)



Product Service