



CERAMILL NP M UNITS

PRODUCT DESCRIPTION

Dental metal ceramic alloy based on CoCr, Type 4 according to DIN ISO 22674 for fixed and removable prosthodontics.


MANUFACTURER

Ceramill NP M corresponds to Colado® CAD CoCr4 from Ivoclar Vivadent AG. Further information can be found in the manufacturer's original manual under the following link

<http://www.ivoclarvivadent.com/en/download-center/>

DELIVERY

The framework is delivered in the rough, unprepared condition.





CERAMILL NP M UNITS

PRODUKTBESCHREIBUNG

Ceramill NP M ist eine dentale Aufbrennlegierung auf CoCr-Basis, Typ 4 gemäß DIN ISO 22674 für festsitzenden und herausnehmbaren Zahnersatz.


HERSTELLER

Ceramill NP M entspricht Colado® CAD CoCr4 von Ivoclar Vivadent AG.
Weitere Informationen sind in der Originalanleitung des Hersteller unter folgendem Link zu finden:

<http://www.ivoclarvivadent.com/de/downloadcenter/>

ANLIEFERUNG

Das gefräste Gerüst wird im rohen, unbearbeiteten Zustand geliefert.



2017-06-01