



CERAMILL TI ALLOY UNITS

PRODUCT DESCRIPTION

Dental titanium milling alloy, Ti Al6 V4, alloy type Grade 5 Type 4 according to DIN ISO 22674 for fixed porcelain-bonded restorations.


MANUFACTURER

Ceramill TI alloy corresponds to Colado® CAD Ti5 from Ivoclar Vivadent AG. Further information can be found in the manufacturer's original manual under the following link

<http://www.ivoclarvivadent.com/en/download-center/>

DELIVERY

The milled substructure is delivered in a raw, unfinished state.





CERAMILL TI ALLOY UNITS

PRODUKTBESCHREIBUNG

Dentale Titanfräslegierung, Ti Al6 V4, Legierungstyp Titan Grade 5 Type 4 gemäß DIN ISO 22674 für festsitzenden metallkeramischen Zahnersatz.


HERSTELLER

Ceramill Ti alloy entspricht Colado® CAD Ti5 von Ivoclar Vivadent AG
Weitere Informationen sind in der Originalanleitung des Hersteller unter folgendem Link zu finden:

<http://www.ivoclarvivadent.com/de/downloadcenter/>

ANLIEFERUNG

Das gefräste Gerüst wird im rohen, unbearbeiteten Zustand geliefert.



2017-06-01