



CERAMILL TI UNITS

PRODUCT DESCRIPTION

Pure dental titanium, titanium Grade 2 Type 3 according to DIN ISO 22674 for fixed porcelain-bonded restorations.


MANUFACTURER

Ceramill Ti corresponds to Colado® CAD Ti2 from Ivoclar Vivadent AG
Further information can be found in the manufacturer's original manual under the following link

<http://www.ivoclarvivadent.com/en/download-center/>

DELIVERY

The milled substructure is delivered in a raw, unfinished state.





CERAMILL TI UNITS

PRODUKTBESCHREIBUNG

Dentales Reintitan, Titan Grade 2 Type 3 gemäß DIN ISO 22674 für festsitzenden metallkeramischen Zahnersatz.


HERSTELLER

Ceramill TI entspricht Colado® CAD Ti2 von Ivoclar Vivadent AG
Weitere Informationen sind in der Originalanleitung des Hersteller unter folgendem Link zu finden:

<http://www.ivoclarvivadent.com/de/downloadcenter/>

ANLIEFERUNG

Das gefräste Gerüst wird im rohen, unbearbeiteten Zustand geliefert.



2017-06-01